

Product datasheet for **SC326192**

EXOC7 (NM_001145299) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EXOC7 (NM_001145299) Human Untagged Clone
Tag:	Tag Free
Symbol:	EXOC7
Synonyms:	2-5-3p; BLOM4; EX070; EXO70; Exo70p; EXOC1; NEDSEBA; YJL085W
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326192 representing NM_001145299.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATCCCCACAGGAGCATCCGCTCGACGGCGGGAGATTGAGGACAAGCTGAAGCAGGAGGAGGAG
ACTCTGTCTTCATCCGAGACAGCCTGGAGAAGAGCGACCAGCTCACTAAGAACATGGTGTCTATCTTA
TCATCCTTTGAGAGCCGCTTATGAAGCTGGAGAAGCTCCATCATCCCTGTGCACAAGCAGACGGAGAAT
CTGCAGCGGCTGCAGGAGAATGTTGAGAAGACGCTGTCTGCCTGGACCATGTATCAGCTACTACCAT
GTGGCCAGTGACACTGAGAAGATCATCAGAGAGGGCCCCACAGGTAGGCTGGAAGAGTACCTGGGAAGC
ATGGCCAAGATTCAGAAGGCAGTGGAGTATTTCCAGGACAACAGCCAGACAGCCCGGAACCTCAACAAA
GTGAAACTGCTCTTTGAGCGCGGAAGGAGGCCCTGGAGTCCGAATTTCCAGCCTGATGACGCGGCAC
AGTAAGGTCGTCTCGCCGTGCTCATCTTGGATCTGATCAGTGGTACGATGATCTGGAGGCCAGGAG
GACGTGACCTGGAGCACCTGCCGAGAGCGTGTCCAGGATGTCATTGCATCTCCCGTGGCTGGTG
GAATATGGCCGAACCAAGATTTTCAAGCTCTACTACCAGATACGCTCCAGCCAGCTGGACCGCTCC
ATCAAAGGACTGAAGGAGCATTTCCATAAGAGCAGTTCTTCTCTGGGGTTCCCTACTCCCTGCTATC
CCCAACAAGAGGAAGACACACCTACCAAGAAGCCAGTCAAGCGGCCAGGACGATCCGTAAAGGCTCAG
AACCTTCTGAAACAGTATTTCCAGCATGGTCTAGATGGGAAAAAGGGGGCTCTAACCTCATTCTCTG
GAAGGTCACGAGCATGATTTCCGAGTTAAGCACCTGTCCGAGGCCCTGAACGACAAGCAGGGCCGCTG
GCCGGGAGAGATGACATGCTGGACGTGGAGACCGATGCCTACATCCACTGCGTCAGTGCCTTCGTCAAG
CTGGCGCAGAGCGAGTACCAGCTGCTGGCCGACATCATCCCCGAGCACCACAGAAGAAGACCTTCGAC
TCCCTGATACAGGATGCCCTGGATGGGCTGATGCTTGAAGGGGAGAACATCGTGTCTGCTGCCGGGAAG
GCCATTGTGCGACACGACTTCTCCAGGTGCTCACCGTCTTCCCCATCCTGCGACACCTCAAGCAGACC
AAGCCTGAGTTTGACCAGGTGCTCCAGGGCACGGCTGCCAGCACAAAGAACAAGCTGCCTGGCCTCATC
ACATCCATGGAGACCATCGGTGCCAAAGCGCTGGAGGACTTCGACAGACAACATCAAGAATGACCCGGAC
AAGGAGTACAACATGCCGAAGGACGGCACCGTACACGAGCTCACCAGCAATGCCATCCTCTTCTGCAG
CAGCTTTTGGACTTCCAGGAGACGGCAGGCGCCATGCTGGCCTCCAAGAGACCAGCTCTTCGGCCACC
AGCTACAGCTCTGAGTTCAGCAAGCGGCTGCTAAGCACCTATATCTGTAAGTGTGGGCAACCTGCAG
TTGAACCTGCTGAGCAAGTCCAAGGTGTACGAGGACCCAGCTCTGAGCGCCATCTTCTGCACAACAAC
TACAATTACATCCTCAAGTCCCTGGAGAAGTCTGAAGTATCCAGCTGGTGGCAGTGACACAGAAGACT
GCTGAGCGCTCCTACCGGAGCACATTGAGCAGCAGATCCAGACCTACCAGCGCAGCTGGTTAAAGGTG
ACTGATTACATCGCAGAGAAGAATCTACCTGTGTTCCAGCCGGGAGTCAAGCTCCGGGACAAGGAGCGG
CAGATTATCAAGGAGCGTTTTAAGGGCTTCAATGATGGCCTCGAAGAAGTGTGCAAAATCCAGAAGGCC
TGGGCTATTCCAGACACAGAGCAGAGGGACAGGATTTCGCCAGGCCAGAAGACCATTGTCAAAGGAGACC
TACGGGGCCTTTTACAGAAGTTTGGCAGCGTGCCTTCCACCAAGAACCAGGAGAAGTACATCAAGTAC
GGGGTGGAGCAGGTGGGCGACATGATCGATCGCTTTTTCGACACCTCTGCCGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001145299

Insert Size: 2124 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145299.3](#)

RefSeq Size: 4874 bp

RefSeq ORF: 2124 bp

Locus ID: 23265

UniProt ID: [Q9UPT5](#)

Cytogenetics: 17q25.1

Protein Families:	Druggable Genome
Protein Pathways:	Insulin signaling pathway
MW:	80.6 kDa
Gene Summary:	<p>The protein encoded by this gene is a component of the exocyst complex. The exocyst complex plays a critical role in vesicular trafficking and the secretory pathway by targeting post-Golgi vesicles to the plasma membrane. The encoded protein is required for assembly of the exocyst complex and docking of the complex to the plasma membrane. The encoded protein may also play a role in pre-mRNA splicing through interactions with pre-mRNA-processing factor 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Nov 2011]</p> <p>Transcript Variant: This variant (6) includes an alternate in-frame coding exon, compared to variant 1, resulting in a longer protein isoform (6).</p>