

## Product datasheet for **SC326202**

### **GIT2 (NM\_001135214) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	GIT2 (NM_001135214) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIT2
Synonyms:	CAT-2; CAT2; PKL
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001135214, the custom clone sequence may differ by one or more nucleotides

```

ATGTCGAAACGGCTCCGGAGCAGCGAGGTGTGCGCTGACTGCAGCGGGCCGGATCCTTCC
TGGGCATCAGTAAATAGGGGAACGTTTTTATGTGATGAGTGCTGCAGTGTCCATCGGAGT
CTAGGGCGCCATATCTCCAAGTGAAGCATCTGAAACACACACCGTGGCCTCCAACACTG
CTTCAGATGGTTGAGACCTTGTATAATAACGGTGCTAACTCTATATGGGAGCATTCTTTG
CTGGACCCTGCGTCTATTATGAGTGAAGACGTAAGCTAATCCACAGGATAAAGTACAT
CCCAATAAAGCGGAATTCATCAGAGCCAAGTATCAGATGTTAGCGTTCGTCCATCGTTG
CCCTGCCGGGATGACGATAGTGTGACTGCCAAAGATCTTAGCAAGCAACTCCATTCGAGC
GTGAGAACAGGAATCTTGAAACCTGTTTGAAGTCTTTAGGAGCACAAGCCAAC
TTCTTTCATCCTGAAAAAGGAAACACCCCACTCCATGTTGCCCTCAAAGCAGGGCAGATT
TTACAGGCTGAATTATTGGCAGTATATGGAGCAGACCCAGGCACACAGGATTCTAGTGGG
AAAACCTCCGTTGATTATGCAAGCAAGGAGGGCACCATGAGCTGGCAGAGCGCCTCGT
GAAATACAGTATGAGCTAACGGACAGACTAGCCTTCTATCTCTGTGGCAGGAAACCAGAT
CACAAAAATGGACAGCACTTTATAATACCTCAAATGGCAGACAGCAGCCTGGATTGTCT
GAATTGGCAAAAGCTGCTAAGAAGAAACTCAATCTCTAAGTAATCATTGTGTTGAAGAA
CTTGCCATGGATGTGTACGATGAAGTTGACAGGCGAGAGACGGATGCAGTCTGGCTTGCC
ACGCAAAACCACAGCGCCTGGTAACCGAGACAACGGTCGTCCCTTTCTCCGGTCAAT
CCTGAGTACTCATCAACACGAAATCAGGGCAGACAGAAGTTAGCTCGGTTCAACGCCAT
GAGTTTGGCCACGCTGGTCATTGACATTCTCAGTGACGCCAAGAGGAGACAGCAGGGCAGT
TCTCTCTCGGTTCAAAGACAATGTGGAGTCTACTGAAAACCATCAATAACCAGCAC
ACAGATTTGGAAACCACTGCAAGCAAAACAAACCGGCAGAAGAGCCTAGATTGATTGA
TCAGATGGACCAGTCACTGTACAGGAATTTATGGAGGTCAAAAACGCTCTAGTGGCTTCT
GAGGCCAAGATACAGCAGCTAATGAAGGTGAATAACAACCTTGAGTGACGAGCTGAGAATT
ATGCAGAAAAAGCTTCAAACACTCCAGAGTGAAAATTCGAACCTCAGGAAACAGGCCACA
ACCAATGTATATCAGGTGCAAACCTGGTTCTGAGTACACAGACACTTCCAACCACTTTCC
TTAAAGAGACGTCCTGCTGCCGGGGCAGTAGGCCCATGTCCATGTACGAGACGGATCA
GGTCAGAAACCATATCTCCAATGGGAGAAGCGAGCCGCCCGAAGAGAGCAGGATGAGA
CTCCAGCCCTTCCCGCGCACGCATCCAGGCTGGAGAAGCAGAACAGCACACCTGAGAGT
GACTACGACAACACTCCCAACGACATGGAGCCAGATGGCATGGGGTCAAGCCGAAAGGGA
CGGCAAAAGAGTATGGTGTGGCCAGGGGATGGCTTGGTACCAGACACAGCAGAACCCCAT
GTGGCCCCAAGCCCCACTCTCCCTAGCACCGAAGATGTCATCAGGAAGACTGAACAGATC
ACCAAAAACATACAGGAGCTCTTAAGAGCAGCCCAAGAAAATAAACATGACAGTTATATT
CCCTGCTCAGAGAGGATACACGTAGCTGTTACAGAAATGGCAGCATTATCCCCAAAAAA
CCCAAGTCTGATATGGTGAAGACTTCCCTTCGTTTACTGACGTCAGTGCCTACCGACTG
CAGTCAGAGTGAAGAAGACCCTCCAGGGGACCCCGGCTCACCCACAGACGTTCCAGCTG
GTCACGCAGCAGGTATCCAGTGTGCGTACGACATCGCCAAGGCTGCCAAGCAGCTGGTT
ACCATCACCAAGAGAACAACAAC

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_001135214

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001135214.1, NP_001128686.1</u>
<b>RefSeq Size:</b>	5553 bp
<b>RefSeq ORF:</b>	2190 bp
<b>Locus ID:</b>	9815
<b>UniProt ID:</b>	<u>Q14161</u>
<b>Cytogenetics:</b>	12q24.11
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Endocytosis
<b>Gene Summary:</b>	<p>This gene encodes a member of the GIT protein family, which interact with G protein-coupled receptor kinases and possess ADP-ribosylation factor (ARF) GTPase-activating protein (GAP) activity. GIT proteins traffic between cytoplasmic complexes, focal adhesions, and the cell periphery, and interact with Pak interacting exchange factor beta (PIX) to form large oligomeric complexes that transiently recruit other proteins. GIT proteins regulate cytoskeletal dynamics and participate in receptor internalization and membrane trafficking. This gene has been shown to repress lamellipodial extension and focal adhesion turnover, and is thought to regulate cell motility. This gene undergoes extensive alternative splicing to generate multiple isoforms, but the full-length nature of some of these variants has not been determined. The various isoforms have functional differences, with respect to ARF GAP activity and to G protein-coupled receptor kinase 2 binding. [provided by RefSeq, Sep 2008]</p> <p>Transcript Variant: This variant (5) lacks one in-frame exon in the 3' coding region, compared to variant 1. The resulting isoform (5) is missing one fragment but has the same N- and C-termini as isoform 1.</p>