

Product datasheet for **SC326209**

ZMYND15 (NM_001136046) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZMYND15 (NM_001136046) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZMYND15
Synonyms:	SPGF14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326209 representing NM_001136046.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGTTTGTGTCTGGATACCGGATGAGTTCCTTGATTTCACTGCCCTTCTCTTCGGCTGGTTCCGA
AAGTTTGTGGCAGAGCGTGGAGCTGTAGGACTAGCCTTGAGGGCCGCTGCCGGCAGCTGGAGGCCAG
ATCAGAAGGCTACCCAGGACCTGCCCTTTGGGTGCTCCATGTCTGCCCAACCATAGTGTGGGATC
AGCCTGGGCAAGGGGCAGAACCAGTCTGGACCAGGCTGGGACTGCCTGGCTCCTGGGAGACAAC
CCTCCACTCCACCTGCGAGACCTGAGCCCTACATCAGCTTTGTCAGCCTAGAGGATGGGAGGAAGGG
GAGGAGGAAGAGGAGGAAGATGAAGAAGAAGAGAAGAGAGAGGACGGGGTGCAGGCAGCACAGAGAAG
GTGGAACAGAGGAGGACCGGGAGCTAGCCCTACCAGCAGGGAGTCCCCCAGGAAACAAACCTCCA
GGAGAGTCAGAGGAGGCTGCCGGGAGGACGAGGTGGCAAGGATGGCTGCCGAGAGGACAGGCTGGAG
AACGAAACAAGACCCAGAAGAGGAAGGACAGAGGAGTGGGCTGCCCCCTGCACGTTTCTGTCTC
TTACTTGTGACGGATGAGCATGGCACCATTGGGCATTGATCTGCTAGTGGATGGAGCCAGGGAACC
GCAAGCTGGGGCTCAGGGACCAAGGACCTGGCTCCTTGGGCTATGCTCTCTCTGTACAGCATGGCC
TGTCCTATGGGCTCTGGGATCCCCGAAAGCCCCGACAGCTTACTGTGGGAGATGCCCGGCTGCATCGA
GAGCTGGAGAGCTTGGTCCCAAGGCTAGGTGTGAAGTTAGCCAAAACCCCAATGCCGACATGGGGTCCC
CGGCCAGGCTTACCTTTGCTTCCCTTCGTGCTCGAACCTGCCATGTGTGTACAGGCACAGCTTTGAA
GCGAAGCTGACACCTGCCCCCAGTGTAGTGCTGTCTTGTATTGTGGAGAGGCTTGTCTCCGGGCTGAC
TGGCAGCGGTGCCAGATGATGTGAGTCACCGATTTTGGTGCCCAAGGCTTGACGCTTCATGGAGCGG
GCAGGAGAACTGGCAACCCTGCCTTTACCTACACCGCAGAGGTGACCACTGAAACCTTCAACAAAGAG
GCCTTCTGGCTCTCGGGCTCACTCGTGGCTATTGGACCCAGCTCAGCATGCTGATTCCAGGCCCG
GGCTTCTCCAGACACCCGAGGCAACACGCGCATCCCTCAGCCTTCTTCGCGGTGGAGACCCTACCAG
CTTCTCCAGGAGACGGGACTGCCCTGATGCTCCTGTGCCCCACATCCACCCGGGTGTTTTTGGC
TCATGGCAGGATTACTACACATGGCGGGCCTCAGCTTGGACTCCCCATAGCCGTGCTTCTCACCTAC
CCGCTGACCGTGTACTACGTATCACCCACCTGGTGCCCGAGTCTTCCCTGAGCTCAACATCCAAAC
AAACAGTCACTGAAGATCCACGTGGTGGAGGCCGGAAGGAGTTTGACCTTGTATGGTGTGTTTGGGAG
CTTTTGGTCTGCTCCCCATGTGGCCTGGAGCTGCAGTTTGTAGGTGATGGCTGCCCGCGAAAGC
GACGAGCAGATTTTACCCTGCAGAGGGACAGCCTGGAGGTGTCTGTCCGGCTGGTTCGGCATATCA
GCACGGCCAGCTCTGGCACTAAGGAGAAAGGGGCCGAGGACCTGCAGATCAAGGTGTGAGCAAGG
CCCTACCACCTGTTCCAGGGGCCAAGCCTGACCTGGTTATTGGATTTAACTCCGGGTTTGTCTCAAG
GATACGTGGCTGAGGTCTCTGCCCGGTTACAGTCCCTCCGAGTGCCAGCCTTCTTACCAGAGAGAGC
GAGTACAGCTGTGTGATGGACGGCCAGACCATGGCGGTGGCCACTGGAGGGGGACACAGCCCTCCCCAG
CCCAACCCCTTCCGCTCCCCCTTCGCCTCAGAGCGGCCGACAACTGCATGTCTGGTACTGCAATGCC
TTCATCTTCCACCTGGTTTACAAGCCTGCTCAAGGAGCGGGGCCCGCCGGCGCCGGGCCCCACCC
CCATCCCCAACTCCCTCTGCTCCTCCTGCCCCACCCGAAGGCGCGGAGGAGAAAAGAAACCTGGGCGG
GGGCCCCGCGCGGAAATGA
ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001136046

Insert Size: 2229 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001136046.2</u>
RefSeq Size:	2424 bp
RefSeq ORF:	2229 bp
Locus ID:	84225
UniProt ID:	<u>Q9H091</u>
Cytogenetics:	17p13.2
MW:	81.9 kDa
Gene Summary:	This gene encodes a MYND-containing zinc-binding protein with a nuclear localization sequence. A similar gene in mice has been shown to act as a testis-specific transcriptional repressor by recruiting histone deacetylase enzymes to regulate spatiotemporal expression of many haploid genes. This protein may play an important role in spermatogenesis. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jun 2012] Transcript Variant: This variant (1) encodes isoform (1).