

Product datasheet for **SC326222**

PWWP2A (NM_001130864) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PWWP2A (NM_001130864) Human Untagged Clone
Tag:	Tag Free
Symbol:	PWWP2A
Synonyms:	MST101
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326222 representing NM_001130864.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGGCGCCGTGGCTGCAGAGGCGGCAGCGACTGCAGCGTCCCCGGGGAGGGGGCGCCGGCGAGGCC
GAGCCGGAGATGGAGCCCATCCCCGGCAGTGAGCCGGCACTGACCCCTCCCGTACAGGCCACTGAA
CGGTCTGTGCCGGATGGCGAGACTGACGGGCAGCAATCCGCTCCTCAGGCCGACGAGCCGCGCTCCCG
CCGCCACCGCCGCGCCGGGGAGCTCGCCCGCAGCCAGAGGCGGTGGGCGCGAGCTGGAGGCTGAG
GAGAACTGTCCGTTCCGGTGGCGGAGTCGGCGGCAGCCGCGCTCAGGGAGGGCCGGAATTCACCT
TCTCTGCATCGCCGCCGAGCAGCCCCGGCTCCCGAGGAGCGCAGGAGCCGCGCTGCCTCAGCCC
GTAGCCCCGGCGCTCGTGCCGCCGGCGGGGACTCCACGGTGTGCAACTGATCCGGGGCTCGGAG
GTGCGGGTCACGCTGGACCACATCATTGAGGACGCGCTTGTCTGTCTGTTCCGCTTCGGGGAGAAGCTC
TTCTCCGGGCTCCTCATGGATCTGTCCAAAGGTTTGGGCCCCATGGTATCCCTGTGACAGTATTTCCC
AAAAGGAATATAAGGATAAACCAGAAGCCATGCCGCTCCAAAGTAATACATTCCAAGAAGGGACAGAA
GTCAAGTGTGAAGCAAATGGTGTCTTCCCGATGACCTTCTCCTGTCCCGCATCCCGAGCTGAGCTTG
GCTGAAAGCCTGTGGACTTCAAACCACCACTCTCTTCATGAAGGAGCACCTTATCTCCCCCTTTG
TTTATCAGGGACACATATAACCAATCAATACCTCAGCCACCTCCTCGGAAAATTAAGCGACCCAAACGA
AAAATGTACAGGAAGAACCCACTTCAATAATGAATGCTATTAACACGACCCAGGCAAGTTCTGTGT
GATAAATGTAAAAACAGTGTTGTTGCTGAAAAAAGGAAATTAGAAAAGGTAGTAGTGCAACTGACTCT
TCTAAATATGAAGATAAAAAACGGAGAAATGAAAGTGTAACACTGTGAACAAAAAAGTGAAGTAC
CATAAAGTGGATGGGAAAAACCAAAATGAAAGCCAGAAAAAGAAATGCTGTGGTTAAGGTTTCAAATATT
GCTCACAGCAGAGGCAGAGTAGTAAAAGTTTCTGCTCAGGCAAAATACATCAAAGCTCAGTTAAGTACT
AAAAAAGTTCTCCAGAGTAAGAACATGGATCATGCGAAAAGCTCGGGAAGTGTAAAAATTGCCAAAGAA
AAGGCACAAAAGAAAGCAAAATGAAACCTCTACTTCCAAAATGCACATTCAAAGTCCATTTACACGT
CGATATCAGAATCCTAGCTCAGGTTCCCTTCCACCCCGGGTTCGTTTAAACACAGAGGTACAGGAAT
GAAGAAAATGACTCTTCTCTGAAGACAGGACTTGAGAAAATGCGGAGTGGCAAGATGGCACCCAAGCCC
CAGTCTCGCTGCACCTCTACCCGCTCAGCAGGTGAGGCCCTTCAGAAAATCAGAGTCCCTCAAAGGC
CCTGAAGAGGCCAGCAGTGAGTTTCCAGGACACAAATGAAGTGCATGTGCTGGTGATCAGGATGAACCA
CAGACATTGGGCAAAAAGGCAGCAAAAACAATATCTCTGTTTATATGACCTAAATCAAAGAAATCT
GACTCTTCCAGTGCTTCACTGTGTAGCATTGATAGCACAGATGATTTGAAATCTTCAACTCTGAGTGT
AGTTCTTCTGAAAGCTTTGATTTTCTCCAGGCAGTATGCATGCACCTTCCACCTCTCCACTTCTCC
TCTTCAAAGGAAGAGAAAAAGCTCAGTAATTCCTTGAAAATGAAAGTCTTTTCCAAAACGCTCTCTAAA
TGCGTCACACCAGATGGCAGGACCATATGTGTAGGGGACATTGTTTGGGCCAAGATATATGGCTTCCT
TGGTGGCCAGCCCGTATTCTTACTATAACTGTGAGCCGGAAGATAACGGCCTTTTAGTCCGACAGGAG
GCCCGTATTTTATGGTTTGGGTCTCCAACAACATCTTCTTGTCTTTTCAACTCTCCCCCTTTTA
GAAAACCTCCAGTCACGCTTTAATAAGAAGAGAAAGGGCCTGTATCGAAGGCTATCAGAGGCGAGCT
AAGGCTGCCAAGCAGCTGACCCCGAAGTGCGGGCTTTGTTGACACAGTTTGAAACGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001130864

Insert Size: 2268 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001130864.1</u>
RefSeq Size:	3394 bp
RefSeq ORF:	2268 bp
Locus ID:	114825
UniProt ID:	<u>Q96N64</u>
Cytogenetics:	5q33.3
MW:	82 kDa
Gene Summary:	H2A.Z-specific chromatin binding protein which may play an important role in the neural crest stem cell migration and differentiation during early development. Also required for proper mitosis progression.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate splice site at its 3'-terminal exon, compared to variant 1, which results in a protein (isoform B) with a longer and distinct C-terminus, compared to isoform A.