

## Product datasheet for **SC326232**

### Breast cancer suppressor candidate 1 (VWA5A) (NM\_001130142) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Breast cancer suppressor candidate 1 (VWA5A) (NM_001130142) Human Untagged Clone
Tag:	Tag Free
Symbol:	VWA5A
Synonyms:	BCSC-1; BCSC1; LOH11CR2A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC326232 representing NM\_001130142.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGTGCACTTCTGTGGCCTACTCACCTCCACCGGGAGCCAGTGCCGCTGAAGAGTATCTCTGTGAGC
GTGAACATTTACGAGTTTGTGGCTGGTGTCTGCAACTTTGAACTACGAGAATGAGGAGAAAGTTCTT
TTGGAGGCCTTCTTTGTGTTCCCATGGATGAAGACTCTGCTGTTTACAGCTTTGAGGCCTTGGTGGAT
GGGAAGAAAATTGTAGCAGAATTACAAGACAAGATGAAGGCCCGACCAACTATGAGAAAAGCCATCTCC
CAGGGCCACCAGGCCTTCTTATTGGAGGGGACAGCAGCTCCAGGGATGTCTTCTTGTCAATGTGGGT
AACCTCCAACCTGGTGAAGGGCGCAGTACCCTGAAGTATGTGCAGGAGCTGCCTCTGGAAGCAGAT
GGGGCTCTGCGCTTTGTGCTCCAGCTGTCTGAATCCTAGATACCAGTTCTCTGGGTGCTAAGGAC
AGTTGCCTAATGTGAAGACTCCTATAGCCCTGTGGAGGACCTGCCCTACACACTCAGCATGGTCGCC
ACCATAGATTCCAGCATGGCATTGAGAAGGTCCAATCCAACCTGCCCTTGAGTCCTACCGAGTACCTA
GGAGAGACAAGACTTCTGCTCAGGTTTCCCTGGCTGCTGGACACAAGTTTGATCGGGACGTGGAACCTC
TGATTTACTACAATGAGGTGCATACCCCAAGCGTGGTTTTGGAGATGGGGATGCCTAACATGAAGCCA
GGTCATTTGATGGGAGATCCATCTGCAATGGTGAGTTTCTATCCAAATATCCCAGAAAGATCAACCATCA
AATACCTGTGGAGAGTTTATCTTCTCATGGACCGCTCGGGAAGTATGCAGAGCCCCATGAGTAGCCAG
GATACATCTCAGCTGCGAATACAGGCAGCCAAGGAAACACTGATTTTGTGCTGAAGAGTTTACCTATA
GGCTGTTATTTCAACATCTATGGATTTGGCTCTTCTATGAGGCATGCTTTCCGGAGAGTGTGAAGTAC
ACTCAGCAAACAATGGAGGAGGCTCTGGGAGAGTGAAGCTTATGCAGGCCGACCTAGGGGGCACTGAA
ATCTTGGCACCCTCCAGAACATTTACAGGGGACCTCCATCCCAGGCCACCCCTACAGCTTTTGTGTC
TTACAGATGGAGAAGTTACAGACAGTTTGTGTAATTAAGAAGTTAGGATCAACAGACAGAAACAC
AGGTGTTTCTCATTGGTATTGGAGAAGGCACCTCCACCAGCCTAATAAAAGGTATTGCCCGGCATCA
GGGGGCACCTCAGAATTTATCACAGGCAAAGACAGGATGCAGTCCAAGGCTCTCAGGACTCTGAAACGC
TCTCTGCAGCCTGTGGTAGAGGATGTCTCTGAGCTGGCATTTCCTCTGCTGTCTGCTAAAATG
CTTTCCCAAGACAGACTGTCATCTTTAGGGGTGAGAGATTAATCAGCTATGCCAGCTGACCGGGAGG
ATGCCAGCAGCAGAGACAACAGGAGAAGTATGCCTCAAATATACACTCCAGGGCAAGACTTTTGAGGAT
AAGGTGACATTTCTCTACAACCAAGCCTGATGTCAACCTCACCATTACCGCCTTGCTGCCAAGTCC
TTGCTCCAGACCAAGGACATGGGCCTCAGGGAGACTCCAGCAAGTGATAAAAAAGATGCATTGAACCTT
AGCCTTGAGTCTGGTGTGATAAGCTCCTCACAGCTTTCATTGCTATCAATAAGGAGCTCAACAAGCCG
GTTCAGGGGCCTCTGGCTCATAGGGACGTCCCAAGGCCAATTCTGTTGGGTGCTTCTGCCCATTAAG
ATAAATGCCAATCAGGTTTTCGAAAGGCCTTACACTCTGACCGTCTCCTTCTGCATCTCAGCCAGA
GGGGAACCTATGTGTTATAAGGCCAAGACATTCAGATGGACGATTACAGTCTCTGTGGGTTGATAAGT
CACAAGGACCAGCACAGTCCAGGCTTTGGAGAGAATCACCTTGTGCAGCTGATTTACCACCAAAATGCA
AATGGTTCTGGGATCTGAATGAAGATCTAGCCAAGATCCTAGGTATGAGTTTGAAGAAATAATGGCT
GCACAGCCTGCCGAGCTTGTGGATTCTCAGGCTGGGCCACCATCCTGGCCGTGATCTGGCTGCACAGC
AATGGTAAGGACTTGAAGTGTGAATGGGAGCTTCTGAAAAGGAAGGCCGTGGCCTGGATGCGTGCCCAT
GCAGGCTCCACCATGCCCTCGTTGTGAAAGCTGCTATTACTTCTGAAGTCATCTGTGGATCTGTCT
ATCTTTGCCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001130142

**Insert Size:** 2361 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001130142.1](#)

**RefSeq Size:** 3552 bp

**RefSeq ORF:** 2361 bp

**Locus ID:** 4013

**UniProt ID:** [O00534](#)

**Cytogenetics:** 11q24.2

**MW:** 86.5 kDa

**Gene Summary:** May play a role in tumorigenesis as a tumor suppressor. Altered expression of this protein and disruption of the molecular pathway it is involved in, may contribute directly to or modify tumorigenesis.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (3) represents the longest transcript and encodes the longer isoform (1).