

Product datasheet for **SC326235**

TRPC4 (NM_001135958) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRPC4 (NM_001135958) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRPC4
Synonyms:	HTRP-4; HTRP4; TRP4
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001135958, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTCAGTTCTATTACAAAAGAAATGTTAATGCTCCCTATAGAGACCGCATCCCTCTA
AGGATAGTAAGAGCAGAATCAGAACTCTCGCCATCAGAAAAAGCCTACTTGAATGCTGTG
GAAAAGGGAGATTATGCCAGTGTCAAGAAATCCCTAGAGGAAGCTGAAATTTATTTTAAA
ATCAATATTAATTGCATTGATCCTCTCGGAAGAAGCTCTCCTCATTGCAATTGAAAAAT
GAGAAGTTGGAGCTCATCGAACTACTCTTAAGCTTTAATGTCTATGTTGGAGATGCTCTA
TTACATGCTATCAGAAAAGAAAGTCGTCGGAGCTGTTGAGCTGTTATTGAACCACAAAAAA
CCTAGTGGAGAAAAACAGTTTGTGCCCCAGCCAAATTGTCAACAGCTGCTGGCATCTCGC
TGGTACGATGAGTTCCAGGCTGGAGGAGAAGACACTGGGCAGTGAAGATGGTGACATGT
TTCATAATAGGACTTCTTTTCTGTCTCTGTGTGCTACCTGATAGCTCCCAAAAGC
CCACTTGGACTGTTTCATCAGGAAGCCATTTATCAAGTTTATCTGCCACACAGCCTCTAT
TTGACTTTTTTGTCTGCTGCTTGCCTCTCAGCACATCGACAGGTCAGACTTGAAC
AGGCAAGTCCACCACCAACCATCGTCGAGTGGATGATATTACCGTGGGCTGGGCTTC
ATATGGGGAGAAATTAACAGATGTGGGATGGCGGACTTCAAGACTACATCCATGATTGG
TGGAACTAATGGACTTTGTAATGAACTCCTTATTTAGCAACAATCTCCTTGAATAAT
GTTGCATTTGTAAGTACAGTGCCTTAATCCACGAGAATCATGGGACATGTGGCATCCC
ACTCTGGTGGCAGAGGCTTTATTTGCTATTGCAAACATCTTCAGTTCTCTGCGTCTGATC
TCACTGTTTACTGCAAAATCTCACCTGGGACCTCTGCAAAATATCTCTGGGAAGAATGCTC
CTGGACATTTTGAAGTTTCTATTCATATACTGCCTTGTGTTGCTAGCATTTGCAAAATGGC
CTAAATCAATTGTACTTCTATTTAAGAAACGAAAGGGTTAACCTGCAAAGGCATAAGA
TGTGAAAAGCAGAATAATGCATTTTCAACGTTATTTGAGACACTGCAGTCCCTGTTTTGG
TCAATATTTGGGCTCATCAATTTATATGTGACCAATGTCAAAGCACAGCATGAATTTACT
GAGTTTTGTTGGTGCCACCATGTTTGGGACATACAATGTCATCTCTCTGGTTGTTCTACTC
AACATGTTAATAGCTATGATGAATAATTCTTACCAACTGATTGCTGACCATGCAGATATA
GAATGGAATTTGCACGAACAAAGCTTTGGATGAGTTATTTGAAGAAGGAGGACTCTG
CCTACTCCCTTCAATGTCATCCCGAGCCCCAAGTCTCTCTGGTACCTGATCAAATGGATC
TGGACACACTTGTGCAAGAAAAAGATGAGAAGAAAGCCAGAAAGTTTGGAAACAATAGGG
AGGCGAGCTGCTGATAACTTGAAGACATCACCATACCAAGAAGTTATGAGGAACCTG
GTGAAGCGATACGTTGCTGCAATGATTAGAGATGCTAAAAGTGAAGAAGGCTGACCGAA
GAGAAGTTAAGGAACTAAAGCAAGACATTTCTAGTTTCCGCTTTGAAGTCTGGGATTA
CTAAGAGGAAGCAAATTTCCACAATAAATCTCGGAATGCCTCGAAGGAGTCTTCAAAT
TCGGCAGACTCAGATGAAAAGAGTGATAGCGAAGGTAATAGCAAGGACAAGAAAAAGAA
TTCAGCCTTTTTGATTTAACCACCTGATTTCATCCGAGATCAGCAGCAATTGCCTCTGAA
AGACATAACATAAGCAATGGCTCTGCCCTGGTGGTTTCAGGAGCCGCCAGGGAGAAGCAG
AGAAAAGTGAATTTGTGACCGATATCAAAAATTTGGGTTATTTTCATAGACGATCAAAA
CAAAATGCTGCTGAGCAAAATGCAAACCAATCTCTCTGTTTCAAGAAGAAGTTGCTCGT
CAACAGGCTGCAGGACCACTTGAAGAAATATCAACTGGAATCTCGAGGATTAGCTTCA
CGGGTGACCTGAGCATTCCCGGTCTCAGTGAACAATGTGTGTTAGTAGACCATAGAGAA
AGGAATACGGACACACTGGGTTACAGGTAGGAAAGAGAGTGTGTCCATTCAAGTCAGAG
AAGGTGGTGGTGGAGGACACGGTTCCTATAATACCAAAGGAGAAACATGAAAAGAGAG
GACTCTAGTATAGACTATGATCTAAACCTCCAGACACAGTCACCCACGAAGATTACGTG
ACCACAAGATTG

```

Restriction Sites: Please inquire

ACCN: NM_001135958

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135958.1 , NP_001129430.1
RefSeq Size:	3015 bp
RefSeq ORF:	2415 bp
Locus ID:	7223
UniProt ID:	Q9UBN4
Cytogenetics:	13q13.3
Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane
Gene Summary:	<p>This gene encodes a member of the canonical subfamily of transient receptor potential cation channels. The encoded protein forms a non-selective calcium-permeable cation channel that is activated by Gq-coupled receptors and tyrosine kinases, and plays a role in multiple processes including endothelial permeability, vasodilation, neurotransmitter release and cell proliferation. Single nucleotide polymorphisms in this gene may be associated with generalized epilepsy with photosensitivity. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2011]</p> <p>Transcript Variant: This variant (6, also known as zeta) lacks an alternate in-frame segment in the 5' coding region, compared to variant 1 (alpha). The encoded isoform (zeta) is shorter, compared to isoform alpha. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>