

Product datasheet for **SC326238**

NELL2 (NM_001145109) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NELL2 (NM_001145109) Human Untagged Clone
Tag:	Tag Free
Symbol:	NELL2
Synonyms:	NRP2
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001145109, the custom clone sequence may differ by one or more nucleotides

```

ATGGAAACGGGACTGGGTGCTCCGCTATTTAAAGCTTGGCTGCTGATCTCAGTTTGGGG
CTTGGTGTGGACCCTTCCCTACAGATTGACGTCTTAACAGAGTTAGAAGTTGGGGAGTCC
ACGACCGGAGTGCGTCAGGTCCCGGGGCTGCATAATGGGACGAAAGCCTTCTCTTTCAA
GATACTCCAGAAAGCATAAAAGCATCCACTGCTACAGCTGAACAGTTTTTTCAGAAGCTG
AGAAATAAACATGAATTTACTATTTTGGTGACCCTAAAACAGACCCACTTAAATTCAGGA
GTTATTCTCTCAATTCACCACTTGGATCACAGGTACCTGGAACTGGAAAGTAGTGGCCAT
CGGAATGAAGTCAGACTGCATTACCGCTCAGGCAGTCACCGCCCTCACACAGAAGTGTTT
CCTTACATTTTGGCTGATGACAAGTGGCACAAGCTCTCCTTAGCCATCAGTGCTTCCCAT
TTGATTTTACACATTGACTGCAATAAAATTTATGAAAGGGTAGTAGAAAAGCCCTCCACA
GACTTGCCTCTAGGCACAACATTTTGGCTAGGACAGAGAAATAATGCGCATGGATATTTT
AAGGGTAAATGCAAGATGTCCAATTACTTGTGTCATGCCCCAGGGATTTATTGCTCAGTGC
CCAGATCTTAATCGCACCTGTCCAATTGCAATGACTTCCATGGACTTGTGCAGAAAATC
ATGGAGCTACAGGATATTTTAGCCAAAACATCAGCCAAGCTGTCTCGAGCTGAACAGCGA
ATGAATAGATTGGATCAGTGCTATTGTGAAAGGACTTGCACCATGAAGGGAACCACCTAC
CGAGAATTTGAGTCTGGATAGACGGCTGTAAGAAGTGCACATGCCTGAATGGAACCATC
CAGTGTGAAACTTAATCTGCCCAAATCCTGACTGCCACTTAAGTGGGCTTTGCGGTAT
GTGGATGGCAAATGCTGTAAGGAATGCAATCGATATGCCAATTTCAAGGACGAACCTAC
TTTGAAGGAGAAAGAAATACAGTCTATTCCTCTTCTGGAGTATGTGTTCTCTATGAGTGC
AAGGACCAGACCATGAACTTGTGAGAGTTTCCAGCTGTCCAGCTTTGGATTGTCCAGAG
TCTCATCAGATAACCTTGTCTCACAGCTGTTGCAAAGTTTGTAAAGGTTATGACTTTTGT
TCTGAAAGGCATAACTGCATGGAGAATTCATCTGCAGAAATCTGAATGACAGGGCTGTT
TGTAGCTGTGAGATGGTTTTAGGGCTCTTCGAGAGGATATGCCTACTGTGAAGACATC
GATGAGTGTGCTGAAGGGCGCCATTACTGTGCGTAAAATACAATGTGTGCAACCCCCG
GGTTCTTTTATGTGCATCTGCAAACTGGATACATCAGAATTGATGATTATTCATGTACA
GAACATGATGAGTGTATCACAAATCAGCACAACCTGTGATGAAAATGCTTTATGCTTCAAC
ACTGTTGGAGGACACAACCTGTGTTTGAAGCCGGGCTATACAGGGAAATGGAACGACATGC
AAAGCATTTTGAAGAGATGGCTGTAGGAATGGAGGAGCCTGTATTGCCGTAATGTGTGT
GCCTGCCACAAGGCTTCACTGGACCCAGCTGTGAAACGGACATTGATGAATGCTCTGAT
GGTTTTGTCAATGTGACAGTCGTGCTAATTGCATTAACCTGCCTGGATGGTACCCTGT
GAGTGCAGAGATGGCTACCATGACAATGGGATGTTTTACCAAGTGGAGAATCGTGTGAA
GATATTGATGAGTGTGGGACCGGAGGCACAGCTGTGCCAATGATACCATTTGCTTCAAT
TTGGATGGCGGATATGATTGTGCGATGTCCTCATGGAAAGAATTGCACAGGGGACTGCATC
CATGATGGAAAAGTTAAGCACAATGGTCAGATTTGGGTGTTGGAAAATGACAGGTGCTCT
GTGTGCTCATGTCAGAATGGATTCGTTATGTGTCGACGGATGGTCTGTGACTGTGAGAAT
CCCACAGTTGATCTTTTTGCTGCCCTGAATGTGACCCAAGGCTTAGTAGTCAGTGCCTC
CATCAAAATGGGAAAACCTTTGTATAACAGTGGTGACACCTGGGTCCAGAATTGTCAACAG
TGCCGCTGCTTGAAGGGGAAGTTGATTGTTGGCCCCTGCCTTGCCAGATGTGGAGTGT
GAATTCAGCATTCTCCAGAGAATGAGTGTGCTGCCGCGCTGTGTACAGACCCCTTGCCAG
GCTGACACCATCCGCAATGACATCACCAGACTTGCCCTGGACGAAATGAATGTGGTTTCG
TTCACCGGCTCCTTGGATCAAACATGGCACTGAGTGTACTCTGCCAGTGAAGAAT
GGCCACATCTGTTGCTCAGTGGATCCACAGTGCCTTCAGGAACTG
    
```

Restriction Sites: Please inquire

ACCN: NM_001145109

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145109.1</u> , <u>NP_001138581.1</u>
RefSeq Size:	3210 bp
RefSeq ORF:	2448 bp
Locus ID:	4753
UniProt ID:	<u>Q99435</u>
Cytogenetics:	12q12
Protein Families:	Secreted Protein, Transmembrane
Gene Summary:	<p>The protein encoded by this gene is a glycoprotein containing several von Willebrand factor C domains and epidermal growth factor (EGF)-like domains. The encoded protein acts as a homotrimer and is found in the cytoplasm. Several variants encoding a few different isoforms exist, and at least one isoform appears to be a secreted protein. Studies in mouse suggest that this protein plays a role in neural cell growth and differentiation as well as in oncogenesis. [provided by RefSeq, Feb 2009]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (c) has a shorter and distinct N-terminus compared to isoform a.</p>