

Product datasheet for **SC326263**

KCTD1 (NM_001142730) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCTD1 (NM_001142730) Human Untagged Clone
Tag:	Tag Free
Symbol:	KCTD1
Synonyms:	C18orf5
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001142730, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGAGAATGCCTGGCAGCGGGGACTGTAACACCAGCGCGGGCGGCAGCGCCAGCGCT
GCCGCCGCCGCCGAGAACAAATGGGGAGCGGGGCGAGGGCGAGCGCGCGGGGGG
CGCGGCCGCCACAGCCGTCGCACTACTGCAGCGGGGCGAGGAGGAGGAGGAAGAG
GAGGAGGAGGACGAGATCCAGGAGGTGCAGATAACGGGGGACGAGGAGGAGGAGGAGGAC
GGAGGTGGGGGGCTGGAGGAGGACGAGGAGGAGGAGGAAGAGGAGGAGATGGGGCTGGAC
TGGGACGAGCCCTGGAGCCGAGGACTCGGCCGGGAGGAGCTGGAGCCGAGCCGGTC
CATATGATCAATATGGACCAGAGCGCCGCTGGAGCCGAGGCCGCCCGCAGCTGCTG
GCGCCCCGGGCCCGGTGGGCCGCCGGGACGGCTCCGAGCTGGACCCGACGTGCTG
CAGCGCCCCGAGCGGGCCGGCTGAGCGAGAACCCCGGCTGGCCACCCGCTACGCCGTG
CGCATCTCCGGGAGTACCTGAGCGAGAAGGCGCAGAGCCCGACTTCGAGACCATGGAC
AAGGGGGCGCTGTGCCGCTGCTGCGCTCCTTCTATGCCGAGGCCGCTCCAAAAGCGGC
CAGCTCTACAGCAAGTCGTGCTCATCAGCATCCGCAGCTCCCTCAACCCTACCTCAAT
GAGCCCCGTAAGCGCAGCTCGACCTACCAAGGACCCGAGCTGCGCAGCGCCAAC
CTGACGCTGGCCCGGTCATCCGCAAGCTCGAGGAGCAGGGCGCCGGCCGGTGGTGCAG
AAGCAAGCCATCACGCGGCCGACCTGCGCAAGCTGTACACCTCCAGCGTCTTACGACCC
AACACGCCCTTCGGGCTGCTCAACAAGGTCTGGTTCGAGACGTGCATGTACTTCTGTACC
CGCGGCCCGGAGAACAGCGTGAGTTGGAGGAGGACTCTTTGGGCTGGCCATGGACGAG
GACGGTCGCAAGTTCGTACTTCAAGTCCCTCGGGCCCTACCACAAGTCGCGCTCGTGC
TCGTGGAGCAAGAAGCGCCGAGAGCAGCGACGAGGAGAAGTTCGCCCGCATGTATGAG
ACGGGCACCGAGTTCTGCCCTACGCCAGCTTCGTCAAGTACCTGCGAAACGCAACCCT
CTCTGCAAGGCTTCTTCCAGCGCCCGGGACCACTGCAGCGAGGGCGATGTGACCTGG
TACGAGAACAAAGCCATCGGCAAGAACTTGCTAGGCACTCGGATGCAGATGCTCTCCAAG
GCGGCCAAGCTCTCCAAGACCTACACCAACCACTGCATCGGCGCGTCTCCATCGCCACG
CTCAACAGCATCGCGGGCATTGGCACCAAGCTGGGCTCGCCCGCCCGCAGGGCTGCTAC
GCCGAGGCTCTGAACGGGGCGGCACGGCACCACTCCACCAACCCCAACCCATCCCTCC
CACCAACCCGCCCCAGCGCCCTCGCTGGGGAACACTTACATCCTCCCAAGACAGC
CAGGTCGGGCCGACGTGAAATCCGAGGCTGCGCCCAAGCGGCCCTGTACGAGTCTGTG
TTCGGGTCGGGGAAATCTGCGGCCCACTTCCCAAAAGACTTTGTATCCGCCCTCG
GAGCCTGTGGATCGGGTGGTGGTGGTTCCGTGAAACACGACCCCTGCCTCTTCTCCA
GAAGCCAATGGGCACAGAAGCACCAATTCTCCCAATAGTTTACCTGCTATTGTTCC
CCCACCCAGGACAGTCGGCCCAATATGTCAAGACCTCTGATCACTAGATCCCCTGCATCT
CCTACTGAACAACCAAGGCATCCCTACTCCAGCACAACCAAAAATCCAATGCGCCTGTC
CACATTGATGTGGGCGGCCACATGTACACCAGCAGCCTGGCCACCCCTACCAAAATACCCT
GAATCCAGAATCGGAAGACTTTTTGATGGTACAGAGCCCATTTGTTTGGACAGTCTCAA
CAGCACTATTTTATTGACAGAGATGGACAGATGTTTCAAGATATATCTTGAATTTTACGA
ACATCCAAACTCCTCATTCCGTGATGTTTCAAGGACTACACTTTGTTATATGAAGAGGCA
AAATATTTTTCAGTTTCAAGCCATGTTGTTGGAGATGGAAGATGGAAGCAGGACAGAGAA
ACTGGTCGATTTTCAAGGCCCTGTGAGTGCCTCGTGTGCGTGTGGCCCCAGACCTCGGA
GAAAGGATCACGCTAAGCGGTGACAAATCCTTGATAGAAGAAGTATTTCCAGAGATCGGC
GACGTGATGTGTAACCTGTCAATGCAGGCTGGAATCACGACTCGACGCACGTCATCAGG
TTTCCACTAAATGGCTACTGTCACTCACTCAGTCCAGGTCCTCGAGAGGTTGCAGCAA
AGAGGATTTGAAATCGTGGGCTCCTGTGGGGGAGGAGTAGACTCGTCCAGTTACAGCAA
TACGTCCTTCGGCGGAACTGAGGCGGACGCCCGTGTACCCTCCGTCATCCGGATAAAG
CAAGAGCCTCTGGAC
    
```

Restriction Sites: Please inquire
ACCN: NM_001142730
Insert Size: 3438 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142730.1, NP_001136202.1</u>
RefSeq Size:	3438 bp
RefSeq ORF:	3438 bp
Locus ID:	284252
UniProt ID:	<u>Q719H9</u>
Cytogenetics:	18q11.2
Protein Families:	Ion Channels: Other
Gene Summary:	<p>This gene encodes a protein containing a BTB (Broad-complex, tramtrack and bric a brac), also known as a POZ (POxvirus and zinc finger) protein-protein interaction domain. The encoded protein negatively regulates the AP-2 family of transcription factors and the Wnt signaling pathway. A mechanism for the modulation of Wnt signaling has been proposed in which the encoded protein enhances ubiquitination and degradation of the beta-catenin protein. Mutations in this gene have been identified in Scalp-ear-nipple (SEN) syndrome. [provided by RefSeq, May 2017]</p> <p>Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>