

## Product datasheet for **SC326264**

### **NELL2 (NM\_001145107) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NELL2 (NM_001145107) Human Untagged Clone
Tag:	Tag Free
Symbol:	NELL2
Synonyms:	NRP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326264 representing NM\_001145107.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGCTTCCCTCTTTGCTTAAAGGGCAGGCGACTACTCGTCCAGCCTTGCTCCTGCAGCTGG
GTGGTCTTTTTCTCCTGTCTTTCAAGACACGCGCCGAAATCGAGGGAGGGAGACGATGGACTGAG
CTGATCCGCACCATGGAGTCTCGGTCTTACTGAGAACATTCTGTTTGATCTTCGGTCTCGGAGCAGTT
TGGGGGCTTGGTGTGGACCCCTCCCTACAGATTGACGTCTAACAGAGTTAGAATTGGGGAGTCCACG
ACCGGAGTGCGTCAGGTCCCGGGCTGCATAATGGGACGAAAGCCTTTCTCTTCAAGATACTCCAGA
AGCATAAAAGCATCCACTGCTACAGCTGAACAGTTTTTTCAGAAGCTGAGAAATAAACATGAATTTACT
ATTTTGGTGACCTAAAACAGACCCACTTAAATTCAGGAGTTATTCTCTCAATTCACCACCTGGATCAC
AGGTACCTGGAAGTGAAGTAGTGGCCATCGGAATGAAGTCAGACTGCATTACCGCTCAGGCAGTCAC
CGCCCTCACACAGAAGTGTTCCTTACATTTTGGCTGATGACAAGTGGCACAAGCTCTCCTTAGCCATC
AGTGCTTCCCATTGATTTTACACATTGACTGCAATAAAATTTATGAAAGGGTAGTAGAAAAGCCCTCC
ACAGACTTGCCTTAGGCCACAACATTTTGGCTAGGACAGAGAAAATATGCGCATGGATATTTAAGGGT
ATAATGCAAGATGTCCAATTACTTGTATGCCCCAGGGATTTATTGCTCAGTGCCAGATCTTAATCGC
ACCTGTCCAATTGCAATGACTTCCATGGACTTGTGCAGAAAATCATGGAGCTACAGGATATTTAGCC
AAAACATCAGCCAAGCTGTCTCGAGCTGAACAGCGAATGAATAGATTGGATCAGTGCTATTGTGAAAGG
ACTTGCACCATGAAGGGAACCACTACCGAGAATTTGAGTCTGGATAGACGGCTGAAGAAGTGCACA
TGCTGAATGGAACATCCAGTGTGAACTCTAATCTGCCAAATCTGACTGCCACTTAAGTCGGCT
CTTGCGTATGTGGATGGCAAATGCTGTAAGGAATGCAAATCGATATGCCAATTTCAAGGACGAACCTAC
TTTGAAGGAGAAAAGAAATACAGTCTATTCCTCTTCTGGAGTATGTGTTCTCTATGAGTCAAGGACCAG
ACCATGAAACTTGTGAGAGTTCAGGCTGTCAGCTTTGGATTGTCCAGAGTCTCATCAGATAACCTTG
TCTCACAGCTGTTGCAAAGTTTGTAAAGTTATGACTTTTGTCTGAAAGGCATAACTGCATGGAGAAT
TCCATCTGCAGAAATCTGAATGACAGGGCTGTTGTAGCTGTCGAGATGGTTTTAGGGCTCTTCGAGAG
GATAATGCCTACTGTGAAGACATCGATGAGTGTGCTGAAGGGCGCCATTACTGTCGTGAAAATACAATG
TGTGTCAACACCCCGGTTCTTTTATGTGCATCTGCAAACTGGATACATCAGAATTGATGATTATTCA
TGTACAGAACATGATGAGTGTATCACAATCAGCACAAGTGTGATGAAAATGCTTTATGCTTCAACT
GTTGGAGGACACAAGTGTGTTGCAAGCCGGGCTATACAGGGAATGGAACGACATGCAAAGCATTTTGC
AAAGATGGCTGTAGGAATGGAGGAGCCTGTATTGCCGCTAATGTGTGTGCTGCCACAAGGCTTCACT
GGACCCAGCTGTGAAACGGACATTGATGAATGCTCTGATGGTTTTGTCAATGTGACAGTCGTGCTAAT
TGCAATTAACCTGCCTGGATGGTACCAGTGTGAGTGCAGAGATGGTACCATGACAATGGGATGTTTTCA
CCAAGTGGAGAATCGTGTGAAGATATTGATGAGTGTGGGACCGGGAGGCACAGCTGTGCCAATGATACC
ATTTGCTTCAATTTGGATGGCGGATATGATTGTCGATGTCCTCATGAAAAGAATTGCACAGGGGACTGC
ATCCATGATGAAAAGTTAAGCACAATGGTCAGATTTGGGTGTTGGAAAATGACAGGTGCTCTGTGTGC
TCATGTCAGAATGGATTCGTTATGTGTCGACGGATGGTCTGTGACTGTGAGAATCCCACAGTTGATCTT
TTTTGCTGCCCTGAATGTGACCCAAGGCTTAGTAGTCAGTGCCTCCATCAAAATGGGAAAATTTGTAT
AACAGTGGTGACACCTGGGTCCAGAATTGTCAACAGTCCCGCTGCTTGAAGGGGAAGTTGATTGTTGG
CCCTGCCTTGCCAGATGTGGAGTGTGAATTCAGCATTCTCCAGAGAATGAGTGTGCTGCCGCGCTGT
GTCACAGACCCCTTGCAGGCTGACACCATCCGCAATGACATCACCAAGACTTGCCTGGACGAAAATGAAT
GTGGTTCGCTTACCGGGTCTCTTGGATCAAACATGGCACTGAGTGTACTCTCTGCCAGTGCAAGAAT
GGCCACATCTGTTGCTCAGTGGATCCACAGTGCCTTCAGGAACTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001145107

**Insert Size:** 2601 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001145107.1](#)

**RefSeq Size:** 3385 bp

**RefSeq ORF:** 2601 bp

**Locus ID:** 4753

**UniProt ID:** [Q99435](#)

**Cytogenetics:** 12q12

<b>Protein Families:</b>	Secreted Protein, Transmembrane
<b>MW:</b>	96.9 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a glycoprotein containing several von Willebrand factor C domains and epidermal growth factor (EGF)-like domains. The encoded protein acts as a homotrimer and is found in the cytoplasm. Several variants encoding a few different isoforms exist, and at least one isoform appears to be a secreted protein. Studies in mouse suggest that this protein plays a role in neural cell growth and differentiation as well as in oncogenesis. [provided by RefSeq, Feb 2009]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a).</p>