

Product datasheet for **SC326268**

PIK3R5 (NM_001142633) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIK3R5 (NM_001142633) Human Untagged Clone
Tag:	Tag Free
Symbol:	PIK3R5
Synonyms:	F730038I15Rik; FOAP-2; p101; P101-PI3K
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326268 representing NM_001142633.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCAGCCAGGGCCACGACATGCACGGAGGACCGCATCCAGCATGCCCTGGAACGCTGCCTGCATGGA
CTCAGCCTCAGCCGCCCTCCACCTCCTGGTCAGCTGGGCTGTGTCTGAACTGCTGGAGCCTGCAGGAG
CTGGTCAGCAGGGACCCGGGCCACTTCTTATCCTCCTTGAGCAGATCCTGCAGAAGACCCGAGAGGTC
CAGGAGAAGGGCACCTACGACCTGCTCACCCCGCTGGCCCTGCTCTTCTATTCCACTGTTCTTTGTACA
CCACACTTCCCACCAGACTCGGATCTCCTTCTGAAGGCAGCCAGCACCTACCACCGGTTCTGACCTGG
CCTGTTCTTACTGCAGCATCTGCCAGGAGCTGCTCACCTTATTGATGCTGAACTCAAGCCCCAGGG
ATCTCCTACCAGAGACTGGTGAGGGCTGAGCAGGGCTGCCATCAGGAGTACCAGCAGCTCCACCGTC
ACCGTGTCTGTGAACCCAGTGAAGTGCAGGCCGAGTTCTTGTGTAGCCAATAAGCTGAGTACG
CCCGGACACTCGCTCACAGTGCCTACACCACCTGCTCCTGCACGCCTCCAGGCCACCTTTGGGGCC
CACTGTGACGTCCCGGCCCTGCACTGCAGGCTACAGGCCAAGACCCTGGCAGAGCTTGAGGACATCTTC
ACGGAGACCCGAGAGGCACAGGAGCTGGCATCTGGCATCGGGGATGCTGCAGAGGCCCGGGGTGGCTC
AGGACCAAGTGCAGGCCGTGGGAGAAAAAGCTGGCTTCCCTGGGGTGTAGACTGCAAAAACAGGG
AAGCTCCACACCATCCCCATCCCTGTGCCAGGTGCTACACCTACAGCTGGAGCCAGGACAGCTTTGAC
ATCCTGCAGGAAATCCTGCTCAAGGAACAGGAGCTACTCCAGCCAGGGATCCTGGGAGATGATGAAGAG
GAGGAAGAGGAGGAGGAGGAGGTGGAGGAGGACTTGAAACTGACGGGCACTGTGCCGAGAGAGATTCC
CTGCTCTCCACAGCTCTTTGGCGTCCCATGACTCCACCTTGTCCCTGCATCCTCCAGGCCTCGGGG
CCGGCCCTCTCGGCCATCTGCTGACTTCTTTGTCTCAGGCCTCTGTATGGCATGGACAGCGGCTAC
GTGGAGACAGCGAGGAGACTCCTCCGAGTGGCCTTGAGGCGTGGCAGCCAGGAAGCCGAGCCAC
CGCAGGCCTGGGAGAAAGTTCATCAGGATCTATAAACTCTTCAAGAGCACCAGCCAGCTGGTACTGCGG
AGGACTCTCGGAGCTGGAGGGCAGCTCGGACACGGCCCTGCCCTGAGGCGGGCAGGGAGCCTCTGC
AGCCCCCTGGACGAACAGTATCACCCCTTCCCGGGCCAGCGCTCCCGCTCCCTGCCACGCCAAA
CTCGGTACCCAGCTGCCAGCTGGCTTCTGGCCCTGCTCACGCCCCAGCGCCGCCCGCCCTTCTG
AGTGGAGATGAGGATCCCAAGGCTTCCACGCTACGTGTTGTGGTCTTTGGCTCCGATCGGATTTAGGG
AAGGTGGCTCGGGCTACAGCAACCTTCGGCGGCTGGAGAACAATCGCCACTCCTCACAGGTTCTTC
AACTTACAGTCTTCTACGTGCCTGTGAAGCGAAGTATGGGACCAGCCCTGGTGCCTGTCCACCCCT
CGGAGCCAGACGCCCTCACCCCGACAGACTCCCTAGGCAGCCAGCCCTGGAGAGCTGGGACACCAC
CCATGGGAGGAGAGCACC AATGACATCTCCACTACCTCGGCATGCTGGACCCCTGGTATGAGCGCAAT
GTACTGGGCTCATGCACCTGCCCCCTGAAGTCTGTGCCAGCAGTCCCTGAAGGCTGAAGCCCAGGCC
CTGGAGGGCTCCCAACCCAGCTGCCATCCTGGCTGACATGCTACTCTACTACTGCCGCTTTGCCGCC
AGACCGGTGCTGCTGCAAGTCTATCAGACCGAGCTGACCTTATCACTGGGGAGAAGACGACAGAGATC
TTCATCCACTCCTGGAGCTGGTCACTCCGCTGCCACACGTGCCATCAAGGCGTCAGGTCTGGCAGC
AAGCGGCTGGGCATCGATGGCGACCGGGAGGCTGTTCTTAACACTACAGATTATTTACAGCAAGGGG
GCCATCAGTGGACGAAGTCGCTGGAGCAACCTGGAGAAGGTCTGTACCTCCGTGAACCTCAACAAGGCC
TGCCGGAAGCAGGAGGAGCTGGATTCCAGCATGGAGGCCCTGACGCTAACTGACAGAAGTGGTGAAA
AGGCAGAACTCCAAATCCAAGAAGGGCTTTAACAGATTAGCACATCGCAGATCAAAGTGGACAAGGTG
CAGATCATCGGCTCCAACAGCTGCCCTTTGCTGTGTGCCTGGACCAGGATGAGAGAAAGATCCTGCAG
AGTGTAGTCAGATGTGAGGTCTCACCGTCTACAAGCCAGAGAAGAGCGACCTCTCCTACCACCCAG
ACGCCTCCTGACCTGCCGCCAGGCCACCTGATCTCTGCTCCCTTCTGCTGCCATCATGACT
TTCAGTGGAGCTCTGCCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001142633

Insert Size: 2643 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142633.2](#)

RefSeq Size: 4543 bp

RefSeq ORF: 2643 bp

Locus ID: 23533

UniProt ID: [Q8WYR1](#)

Cytogenetics: 17p13.1

Protein Families:	Druggable Genome
Protein Pathways:	Acute myeloid leukemia, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, Leukocyte transendothelial migration, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Phosphatidylinositol signaling system, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, VEGF signaling pathway
MW:	97.3 kDa
Gene Summary:	<p>Phosphatidylinositol 3-kinases (PI3Ks) phosphorylate the inositol ring of phosphatidylinositol at the 3-prime position, and play important roles in cell growth, proliferation, differentiation, motility, survival and intracellular trafficking. The PI3Ks are divided into three classes: I, II and III, and only the class I PI3Ks are involved in oncogenesis. This gene encodes the 101 kD regulatory subunit of the class I PI3K gamma complex, which is a dimeric enzyme, consisting of a 110 kD catalytic subunit gamma and a regulatory subunit of either 55, 87 or 101 kD. This protein recruits the catalytic subunit from the cytosol to the plasma membrane through high-affinity interaction with G-beta-gamma proteins. Multiple alternatively spliced transcript variants encoding two distinct isoforms have been found. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (1) and variant 2 encode the same protein (isoform 1).</p>