

Product datasheet for **SC326301**

NEDD4 2 (NEDD4L) (NM_001144968) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NEDD4 2 (NEDD4L) (NM_001144968) Human Untagged Clone
Tag:	Tag Free
Symbol:	NEDD4 2
Synonyms:	hNEDD4-2; NEDD4-2; NEDD4.2; PVNH7; RSP5
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001144968, the custom clone sequence may differ by one or more nucleotides

```

ATGAGACGTCTCGCATTGAGCAGGGAGAGTCCCGTATTCTCAGAGTAAAAGTTGTTTCT
GGAATTGATCTCGCCAAAAGGACATCTTTGGAGCCAGTGATCCGTATGTGAACTTTCA
TTGTACGTAGCGGATGAGAATAGAGAATTGCTTTGGTCCAGACAAAAACAATAAAAAG
ACACTGAACCCAAAATGGAATGAAGAATTTTATTTACGGGTAAACCCATCTAATCACAGA
CTCCTATTTGAAGTATTTGACGAAAATAGACTGACACGAGACGACTTCTGGGCCAGGTG
GACGTGCCCTTAGTCACCTTCCGACAGAAGATCCAACCATGGAGCGACCCTATACATTT
AAGGACTTTCTCCTCAGACCAAGAAGTCATAAGTCTCGAGTTAAGGGATTTTGGCATTG
AAAATGGCCTATATGCCAAAAATGGAGGTCAAGATGAAGAAAACAGTGACCAGAGGGAT
GACATGGAGCATGGATGGGAAGTTGTTGACTCAAATGACTCGGCTTCTCAGCACCAAGAG
GAACTTCTCCTCCTCTGCTCCCGGGTGGGAAGAAAAAGTGGACAATTTAGCCGA
ACTTACTATGTCAACCACAACAACCGGACCACTCAGTGGCACAGACCAAGCCTGATGGAC
GTGTCCTCGGAGTCGGACAATAACATCAGACAGATCAACCAGGAGGCAGCACACCCGGCC
TTCCGCTCCCGCAGGCACATCAGCGAAGACTTGGAGCCCAGCCCTCGGAGGGCGGGGAT
GTCCCGAGCCTTGGGAGACCATTTAGAGGAAGTGAATATCGCTGGAGACTCTCTCGGT
CTGGCTCTGCCCCACCACCGGCCTCCCCAGGATCTCGGACCAGCCCTCAGGAGCTGTCA
GAGGAATAAGCAGAAGGCTTCAGATCACTCCAGACTCCAATGGGGAACAGTTTCACTCT
TTGATTCAAAGAGAACCCTCCTCAAGGTTGAGGTGATGCAGTGTACCCGACGCAGTTGCA
GAACAGGGCCATCTACCACCGCCAGTGGCCAGCTGGGAGAGCGCGTTCATCACTGTC
ACGGGTGGTGAGGAACCAACGCCATCAGTGGCCTATGTACATACCACGGCCGGTCTGCT
TCAGGCTGGGAAGAAAAGAAAAGATGCTAAGGGGCGCACATACTATGTCAATCATAACAAT
CGAACCAACTTGGACTCGACCTATCATGCAGCTTGCAGAAGATGGTGGCTCCGGATCA
GCCACAACAGTAACAACCATCTAATCGAGCCTCAGATCCGCCGGCCTCGTAGCCTCAGC
TCGCCAACAGTAACTTTATCTGCCCGCTGGAGGGTGCCAAGGACTCACCCGTACGTCGG
GCTGTGAAAGACACCCTTCCAACCCACAGTCCCCACAGCCATCACCTTACAACCTCCCC
AAACCACAACAAAAGTACACAGAGCTTCTTGCACCCGGCTGGGAAATGAGGATAGCG
CCAAACGGCCGGCCCTTCTTATTGATCATAACACAAAGACTACAACCTGGGAAGATCCA
CGTTTGAAATTTCCAGTACATATGCGGTCAAAGACATCTTTAAACCCCAATGACCTTGGC
CCCCTTCTCCTGGCTGGGAAGAAAGAATTCATTGGATGGCCGAACGTTTTATATTGAT
CATAATAGCAAAATTAAGTGGGAAGACCCAAGACTGCAGAACCAGCTATTACTGGT
CCGGCTGTCCCTACTCCAGAGAATTTAAGCAGAAATATGACTACTTCAGGAAGAAATTA
AAGAAACCTGCTGATATCCCAATAGGTTTGAATGAACTTCACAGAAATAACATATTT
GAAGAGTCCTATCGGAGAATTATGTCCGTGAAAAGACCAGATGTCCTAAAAGCTAGACTG
TGGATTGAGTTTGAATCAGAGAAAGGCTTGTACTATGGGGTGTGGCCAGAGAATGGTTC
TTCTTACTGTCCAAAGAGATGTTCAACCCCTACTACGGCCTCTGTAATGAGGATCATTGTCC
GACAACTACACCCTCAGATCAACCCTAATTCAGGCCTCTGTAATGAGGATCATTGTCC
TACTTCACTTTTATTGGAAGAGTTGCTGGTCTGGCCGATTTTCATGGGAAGCTCTTAGAT
GGTTTTCATTAGACCATTTTACAAGATGATGTTGGGAAAGCAGATAACCCTGAATGAC
ATGGAATCTGTGGATAGTGAATATTACAACCTTTGAAATGGATCCTGGAGAATGACCCT
ACTGAGCTGGACCTCATGTTCTGCATAGACGAAGAAAACCTTGGACAGACATATCAAGTG
GATTTGAAAGCCCAATGGGTGAGAAATAATGGTCACAAATGAAAACAAAAGGGAATATATC
GACTTAGTCATCCAGTGGAGATTTGTGAACAGGGTCCAGAAGCAGATGAACGCCTTCTTG
GAGGGATTCACAGAACTACTTCTATTGATTTGATTTAAATTTTTGATGAAAATGAGCTG
GAGTTGCTCATGTGCGCCTCGGTGATGTGGATGTGAATGACTGGAGACAGCATTCTATT
TACAAGAACGGCTACTGCCAAACCACCCCGTCATTGAGTGGTTCTGGAAGGCTGTGCTA
CTCATGGACGCCAAAAGCGTATCCGGTACTGCAGTTTGTACAGGGACATCGCGAGTA
CCTATGAATGGATTTGCCGAACCTTATGTTCCAATGGTCTCAGCTGTTTACAATAGAG
CAATGGGGCAGTCTGAGAAACTGCCAGAGCTCACACATGCTTTAATCGCCTTGACTTA
CCTCCATATGAAACCTTTGAAGATTTACGAGAGAACTTCTCATGGCCGTGGAAAATGCT
CAAGGATTTGAAGGGTGGAT

```

Restriction Sites:	Please inquire
ACCN:	NM_001144968
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001144968.1</u> , <u>NP_001138440.1</u>
RefSeq Size:	8378 bp
RefSeq ORF:	2904 bp
Locus ID:	23327
UniProt ID:	<u>Q96PU5</u>
Cytogenetics:	18q21.31
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis, Ubiquitin mediated proteolysis
Gene Summary:	<p>This gene encodes a member of the Nedd4 family of HECT domain E3 ubiquitin ligases. HECT domain E3 ubiquitin ligases transfer ubiquitin from E2 ubiquitin-conjugating enzymes to protein substrates, thus targeting specific proteins for lysosomal degradation. The encoded protein mediates the ubiquitination of multiple target substrates and plays a critical role in epithelial sodium transport by regulating the cell surface expression of the epithelial sodium channel, ENaC. Single nucleotide polymorphisms in this gene may be associated with essential hypertension. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Mar 2012]</p> <p>Transcript Variant: This variant (g) uses an alternate first exon, compared to variant j, which results in the use of a different start codon and a protein (isoform 4) with a shorter and distinct N-terminus, compared to isoform 1. Sequence Note: The RefSeq transcript was derived from the reference genome assembly. The genomic coordinates were determined from alignments.</p>

