

Product datasheet for **SC326311**

CD51 (ITGAV) (NM_001144999) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD51 (ITGAV) (NM_001144999) Human Untagged Clone
Tag:	Tag Free
Symbol:	CD51
Synonyms:	CD51; MSK8; VNRA; VTNR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326311 representing NM_001144999.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCCTAGGCACCCTCCTTCTGATCCTGTACATCTTAATGTTGTGCCGGATGTTTCTTCTCGTGGGA
GCTCCCAAAGCAAACACCACCCAGCCTGGGATTTGGAAGGAGGGCAGGTCTCAATGTGACTGGTCT
TCTACCCGCCGGTGCCAGCCAATTGAATTTGATGCAACAGGCAATAGAGATTATGCCAAGGATGATCCA
TTGGAATTTAAGTCCCATCAGTGGTTTGGAGCATCTGTGAGGTCGAAACAGGATAAAAATTTGGCCTGT
GCCCCATTGTACCATTGGAGAAGTGAAGTGAACAGGAGCGAGAGCCTGTTGGAACATGCTTTCTCAA
GATGGAACAAAGACTGTTGAGTATGCTCCATGTAGATCACAAGATATTGATGCTGATGGACAGGGATTT
TGTCAAGGAGGATTCAGCATTGATTTTACTAAAGCTGACAGAGTACTTCTGGTGGTCTGGTAGCTTT
TATTGGCAAGGTCAGCTTATTTCCGGATCAAGTGGCAGAAATCGTATCTAAATACGACCCCAATGTTTAC
AGCATCAAGTATAATAACCAATTAGCAACTCGGACTGCACAAGCTATTTTGGTATGACAGCTATTTGGGT
TATTCTGTGGCTGTCGGAGATTTCAATGGTGTGATGGCATAGATGACTTTGTTTCAGGAGTTCCAAGAGCA
GCAAGGACTTTGGGAATGGTTTATATTTATGATGGGAAGAACATGTCCTCTTATACAATTTACTGGC
GAGCAGATGGCTGCATATTTCCGATTTTCTGTAGCTGCCACTGACATTAATGGAGATGATTATGCAGAT
GTGTTTATTGGAGCACCTCTCTTCATGGATCGTGGCTCTGATGGCAAACCTCCAAGAGGTGGGGCAGGTC
TCAGTGTCTCTACAGAGAGCTTCAGGAGACTTCCAGACGACAAAGCTGAATGGATTTGAGGTCTTTGCA
CGGTTTGGCAGTGCCATAGCTCCTTTGGGAGATCTGGACCAGGATGGTTTCAATGATATTGCAATTGCT
GCTCCATATGGGGTGAAGATAAAAAAGGAATTGTTTATATCTTCAATGGAAGATCAACAGGCTTGAAC
GCAGTCCCATCTCAAATCCTTGAAGGCGAGTGGCTGCTCGAAGCATGCCACCAAGCTTTGGCTATTCA
ATGAAAGGAGCCACAGATATAGACAAAAATGGATATCCAGACTAATTGTAGGAGCTTTGGTGTAGAT
CGAGCTATCTTATACAGGGCCAGACCAGTTATCACTGTAAATGCTGGTCTTGAAGTGTACCTAGCATT
TTAAATCAAGACAATAAAACCTGCTCACTGCCTGGAACAGCTCTCAAAGTTTCTGTTTTAATGTTAGG
TTCTGCTTAAAGGCAGATGGCAAAGGAGTACTTCCAGGAAACTTAATTTCCAGGTGGAACCTCTTTTG
GATAAACTCAAGCAAAAGGGAGCAATTCGACGAGCACTGTTTCTCTACAGCAGGTCCCAAGTCACTCC
AAGAACATGACTATTTCAAGGGGGGACTGATGCAGTGTGAGGAATTGATAGCGTATCTGCCGGATGAA
TCTGAATTTAGAGACAACTCACTCCAATTACTATTTTTATGGAATATCGGTTGGATTATAGAACAGCT
GCTGATACAACAGGCTTGCAACCCATTCTAACCAGTTCACGCCTGCTAACATTAGTCGACAGGCTCAC
ATTCTACTTGACTGTGGTGAAGACAATGTCTGTAAACCAAGCTGGAAGTTTCTGTAGATAGTGATCAA
AAGAAGATCTATATTGGGGATGACAACCCCTCTGACATTGATTGTTAAGGCTCAGAATCAAGGAGAAGGT
GCCTACGAAGCTGAGCTCATCGTTTCCATTCCACTGCAGGCTGATTTTCATCGGGTTGTCGAAACAAT
GAAGCCTTAGCAAGACTTTCCTGTGCATTTAAGACAGAAAACCAAACCTCGCCAGGTGGTATGTGACCTT
GGAAACCAATGAAGGCTGGAACCAACTCTTAGCTGGTCTTCGTTTCAGTGTGCCACCAGCAGTCAGAG
ATGGATACTTCTGTGAAATTTGACTTACAAATCCAAAGCTCAAATCTATTTGACAAAGTAAGCCAGTT
GTATCTCACAAGTTGATCTTGTGTTTTAGCTGCAGTTGAGATAAGAGGAGTCTCGAGTCTGTATCAT
ATCTTTCTCCGATTCCAAACCTGGGAGCACAAGGAGAACCCTGAGACTGAAGAAGATGTTGGGCCAGTT
GTTCAGCACATCTATGAGCTGAGAAAACAATGGTCCAAGTTCATTGCAAGGCAATGCTCCATCTTCAG
TGGCCTTACAAATATAATAATAAACAACCTGTTGTATATCCTTCAATGATATTGATGGACCAATGAAC
TGCACTTCAGATATGGAGATCAACCCCTTGGAGAATTAAGATCTCATCTTTGCAAACAACGAAAAGAAAT
GACACGGTTGCCGGCAAGGTGAGCGGGACCATCTCATCACTAAGCGGGATCTTGCCCTCAGTGAAGGA
GATATTACACTTTGGTTGTGGAGTTGCTCAGTCTTGAAGATTGTCTGCCAAGTTGGGAGATTAGAC
AGAGGAAAGAGTGAATCTTGTACGTAAGTCATTACTGTGGACTGAGACTTTTATGAATAAAGAAAAT
CAGAATCATTCTATTCTGAAGTCGTCTGCTTCATTTAATGTCATAGAGTTTCTTATAAGAATCTT
CCAATTGAGGATATCACAACCTCCACATGGTTACCCTAATGTCACCTGGGGCATTGAGCCAGCGCCC
ATGCCCTGTGCTGTGGGTGATCATTAGCAGTTCTAGCAGGATTGTTGCTACTGGCTGTTTTGGTA
TTTTGTAATGTACAGGATGGGCTTTTTTAAACGGGTCCGGCCACCTCAAGAAGAACAAGAAAGGGAGCAG
CTTCAACCTCATGAAAATGGTGAAGGAACTCAGAAACTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001144999
Insert Size:	3009 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001144999.1
RefSeq Size:	6769 bp
RefSeq ORF:	3009 bp
Locus ID:	3685
UniProt ID:	P06756
Cytogenetics:	2q32.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cell adhesion molecules (CAMs), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung cancer
MW:	111.1 kDa

Gene Summary:

The product of this gene belongs to the integrin alpha chain family. Integrins are heterodimeric integral membrane proteins composed of an alpha subunit and a beta subunit that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha V subunit. This subunit associates with beta 1, beta 3, beta 5, beta 6 and beta 8 subunits. The heterodimer consisting of alpha V and beta 3 subunits is also known as the vitronectin receptor. This integrin may regulate angiogenesis and cancer progression. Alternative splicing results in multiple transcript variants. Note that the integrin alpha 5 and integrin alpha V subunits are encoded by distinct genes. [provided by RefSeq, Oct 2015]

Transcript Variant: This variant (2) differs in the 5' UTR and 5' coding region and uses an alternate start codon, compared to variant 1. The encoded isoform (2) has a distinct and shorter N-terminus, compared to isoform 1. It is not known whether this isoform (2) is proteolytically processed in the same manner as isoform 1.