

Product datasheet for **SC326453**

PAFAH1B3 (NM_001145939) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PAFAH1B3 (NM_001145939) Human Untagged Clone
Tag:	Tag Free
Symbol:	PAFAH1B3
Synonyms:	PAFAHG
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001145939 edited CGGGAAGCGCCGGGCGAGCCCACTGTGCGGCCGTCGTGGGGGAAGCGAAGGTTCCCTGATT CCACCTACCTCCTAAGTGAAGGCTTCCGTCCCTGGAGAGGAGGCGGCCCTGTATCGGC CTTCGTCCCTCCGCGGGTGTGCTAGCGTTGGGACGGTCCTTTGTTGCCGCGAGGGGTAGGA GTGGGCGTGGCGGAGCCAGCTCCGTTCCGGAACACTCCCGGGCCGACCCGACTCGCTCATC CTGCAGGAGCTGCGGCGCAAGATGAGTGGAGAGGAGAACCCAGCCAGCAAGCCACGCC GGTGCAGGACGTACAGGGCGACGGGCGCTGGATGTCCCTGCACCATCGGTTCTGGCTGA CAGCAAAGATAAGGAACCCGAAGTCGTCTTCATCGGGGACTCCTTGGTCCAGCTCATGCA CCAGTGCAGATCTGGCGGAGCTCTTCTCCTCTGCATGCACTTAACCTTTGGCATTGG TGGTGACGGCACACAGCATGTAAGTGTGGCGGCTGGAGAATGGGGAGCTGGAACACATCCG GCCAAGATTGTGGTGGTCTGGGTGGGACCAACAACCACGGACACACAGCAGAGCAGGT GACTGGTGGCATCAAGGCCATTGTGCAACTGGTGAATGAGCGACAGCCCCAGGCCGGGT TGTGGTGTGGGCTGCTTCCGCGAGGCCAACATCCCAACCCACTTCGGGAGAAGAACC ACAGGTGAACGAGCTGGTACGGGCGGCACTGGCTGGCCACCCTCGGGCCCACTTCCCTAGA TGCCGACCCTGGCTTTGTGCACTCAGATGGCACCATCAGCCATCATGACATGTATGATTA CCTGCATCTGAGCCGCTGGGCTACACACTGTTTGGCGGGCTCTGCACTCCCTGCTTCT GCGTCTGCTGGCCCAAGACCAGGGCCAAGGTGCTCCCCTGCTGGAGCCCGCACCCCTAAGC ATCCTGCTGCCTCCCACAACATTAACCTCTCCTCCTCAGTAAAAAAAAAAAAAAAAAAAA A
Restriction Sites:	Please inquire
ACCN:	NM_001145939
Insert Size:	1000 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145939.1 , NP_001139411.1
RefSeq Size:	1108 bp
RefSeq ORF:	696 bp
Locus ID:	5050
UniProt ID:	Q15102
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome
Protein Pathways:	Ether lipid metabolism, Metabolic pathways
Gene Summary:	<p>This gene encodes an acetylhydrolase that catalyzes the removal of an acetyl group from the glycerol backbone of platelet-activating factor. The encoded enzyme is a subunit of the platelet-activating factor acetylhydrolase isoform 1B complex, which consists of the catalytic beta and gamma subunits and the regulatory alpha subunit. This complex functions in brain development. A translocation between this gene on chromosome 19 and the CDC-like kinase 2 gene on chromosome 1 has been observed, and was associated with cognitive disability, ataxia, and atrophy of the brain. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2 and 3 encode the same protein.</p>