

Product datasheet for **SC326473**

RBM17 (NM_001145547) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	RBM17 (NM_001145547) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM17
Synonyms:	SPF45
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326473 representing NM_001145547. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCCTGTACGATGACCTAGGAGTGGAGACCAGTGACTCAAAAACAGAAGGCTGGTCCAAAACCTTC
AAACTTCTGCAGTCTCAGCTTCAGGTGAAGAAGGCAGCTCTCACTCAGGCAAAGAGCCAAAGGACGAAA
CAAAGTACAGTCTCGCCCGAGTCATTGACCTGAAGCGAGGTGGCTCCTCAGATGACCGGCAAATTTGTG
GACTCCACCGCATGTAGCAGTGGGCTGAAGGATCCTGTTCCAGTGGGTTTTCTGCAGGGGAAGTT
CTGATTCCTTAGCTGACGAATATGACCCTATGTTTCTAATGATTATGAGAAAGTAGTGAAGCGCCAA
AGAGAGGAACGACAGAGACAGCGGGAGCTGGAAGACAAAAGGAAATAGAAGAAAGGGAAAAAGGCGT
AAAGACAGACATGAAGCAAGTGGGTTTGAAGGAGACCAGATCCAGATTCTGATGAAGATGAAGATTAT
GAGCGAGAGAGGAGGAAAAGAAGTATGGGCGGAGCTGCCATTGCCCCACCCACTTCTCTGGTAGAGAAA
GACAAAGAGTTACCCCGAGATTTTCTTATGAAGAGGACTCAAGACCTCGATCACAGTCTTCCAAAGCA
GCCATTCTCCCCAGTGTACGAGGAACAAGACAGACCGAGATCTCAACCGGACCTAGCAACTCTTTC
CTCGCTAACATGGGGGACGCGTGGCGACAAGATCATGCAGAAGTACGGCTTCCGGGAGGGCCAGGGT
CTGGGGAAGCATGAGCAGGGCTGAGCACTGCCTTGTGAGTGGAGAAGACCAGCAAGCGTGGCGGCAAG
ATCATCGTGGGCGACGCCACAGAGAAAGATGCATCCAAGAAGTCAAGTCAAATCCGCTGACTGAAATA
CTTAAGTGTCTACTAAAGTGGTCTTACTAAGGAACATGGTGGTGGGAGAGGTTGGATGAAGACTTG
GAAGTTGAAACCAAGGAAGAATGTGAAAAATATGGCAAAGTTGGAAAATGTGTGATATTTGAAATTCCT
GGTGCCCTGATGATGAAGCAGTACGGATATTTTGAATTTGAGAGAGTTGAATCAGCAATTAAGCG
GTTGTTGACTTGAATGGGAGGTATTTGGTGGACGGGTGTAAGGATGTTTCTACAATTTGGACAAA
TTCAGGGTCTTGATTTGGCAGAACAAGTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



Plasmid Map:


ACCN: NM_001145547

Insert Size: 1206 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145547.1](#)

RefSeq Size: 3760 bp

RefSeq ORF: 1206 bp

Locus ID: 84991

UniProt ID: [Q96I25](#)

Cytogenetics: 10p15.1

Protein Pathways: Spliceosome

MW: 45 kDa

Gene Summary: This gene encodes an RNA binding protein. The encoded protein is part of the spliceosome complex and functions in the second catalytic step of mRNA splicing. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 9 and 15. [provided by RefSeq, Mar 2009]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.