

## Product datasheet for **SC326481**

### **GATA2 (NM\_001145662) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	GATA2 (NM_001145662) Human Untagged Clone
Tag:	Tag Free
Symbol:	GATA2
Synonyms:	DCML; IMD21; MONOMAC; NFE1B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC326481 representing NM\_001145662.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGGTGGCGCCCGAGCAGCCGCGCTGGATGGCGCACCCGGCCGTGCTGAATGCGCAGCACCCCGAC
TCACACCACCCGGGCGCTGGCGCACAACTACATGGAACCCGCGCAGCTGCTGCCTCCAGACGAGGTGGAC
GTCTTCTCAATCACTCGACTCGCAGGGCAACCCTACTATGCCAACCCGCTCACGCGGGGCGCGC
GTCTCCTACAGCCCCGCGCACGCCCGCCTGACCGGAGGCCAGATGTGCCGCCCACTTGTTCACAGC
CCGGGTTTGCCCTGGCTGGACGGGGGCAAAGCAGCCCTCTCTGCCGCTGCGGCCACCACCACAACCCC
TGGACCGTGAGCCCTTCTCCAAGACGCCACTGCACCCCTCAGCTGCTGGAGGCCCTGGAGGCCACTC
TCTGTGTACCCAGGGGCTGGGGTGGGAGCGGGGAGGCAGCGGGAGCTCAGTGGCTCCCTCACCCCT
ACAGCAGCCACTCTGGTCCCACCTTTTCGGCTTCCCACCCACGCCACCCAAAGAAGTGTCTCCTGAC
CCTAGCACCGGGGCTGGCTTCCAGCCTCATCTCCGCGGGGGTAGTGCAGCCCGAGGAGAGGAC
AAGGACGGCGTCAAGTACCAGGTGTCACTGACGGAGAGCATGAAGATGGAAGTGGCAGTCCCCTGCGC
CCAGGCCCTAGCTACTATGGGCACCCAGCCTGCTACACACCACCCCATCCCACCTACCCCTCCTATGTG
CCGGCGGCTGCCACGACTACAGCAGCGGACTTTCACCCCGGAGGCTTCTGGGGGGACCGGCCCTCC
AGCTTCACCCCTAAGCAGCGCAGCAAGGCTCGTTCCTGTTCAGAAGCCGGGAGTGTGTCAACTGTGGG
GCCACAGCCACCCCTCTCTGGCGCGGGACGGCACCCGGCCACTACCTGTGCAATGCCTGTGGCCTTAC
CACAAGATGAATGGGCAGAACCAGCCACTCATCAAGCCCAAGCGAAGACTGACGACAACCACCACCTTA
TGGCGCGAAACGCCAACGGGGACCCTGTCTGCAACGCCTGTGGCCTCTACTACAAGTGCACAATGTT
AACAGGCCACTGACCATGAAGAAGGAAGGGATCCAGACTCGGAACCGGAAGATGTCCAACAAGTCCAAG
AAGAGCAAGAAAGGGGCGAGTGTCTCGAGGAGCTGCAAGTGCATGCAGGAGAAGTCAATCCCCCTTC
AGTGCAGCTGCCCTGGCTGGACACATGGCACCTGTGGGCCACCTCCCGCCCTTCAGCCACTCCGGACAC
ATCCTGCCACTCCGACGCCATCCACCCCTCCTCCAGCCTCTCCTTGGCCACCCACCCCGTCCAGC
ATGGTGACCGCCATGGGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



**ACCN:** NM\_001145662

**Insert Size:** 1401 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001145662.1</a></u>
<b>RefSeq Size:</b>	3263 bp
<b>RefSeq ORF:</b>	1401 bp
<b>Locus ID:</b>	2624
<b>UniProt ID:</b>	<u><a href="#">P23769</a></u>
<b>Cytogenetics:</b>	3q21.3
<b>Protein Families:</b>	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
<b>MW:</b>	49.1 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the GATA family of zinc-finger transcription factors that are named for the consensus nucleotide sequence they bind in the promoter regions of target genes. The encoded protein plays an essential role in regulating transcription of genes involved in the development and proliferation of hematopoietic and endocrine cell lineages. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Mar 2009]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR and uses an alternate splice site in the CDS but maintains the reading frame, compared to variant 1. This variant encodes isoform 2, which is shorter than isoform 1.</p>