

## Product datasheet for **SC326482**

### **NONO (NM\_001145409) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NONO (NM_001145409) Human Untagged Clone
Tag:	Tag Free
Symbol:	NONO
Synonyms:	MRXS34; NMT55; NRB54; P54; P54NRB; PPP1R114
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC326482 representing NM\_001145409.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGAGTAATAAACTTTTAACTTGGAGAAGCAAAACCATACTCCAAGAAAGCATCATCAACATCAC
CACCAGCAGCAGCACACCAGCAGCAACAGCAGCAGCCGCCACCACCAGCAATACCTGCAATGGGCAA
CAGGCCAGCAGCCAAAATGAAGGCTTGACTATTGACCTGAAGAATTTTAGAAAACCAGGAGAGAAGACC
TTCACCCAACGAAGCCGCTCTTTTGTGGGAAATCTTCTCCCGACATCACTGAGGAAGAAATGAGGAAA
CTATTTGAGAAATATGGAAAGGCAGGCGAAGTCTTCATTCATAAGGATAAAGGATTTGGCTTTATCCGC
TTGAAACCCGAACCCTAGCGGAGATTGCCAAAGTGGAGCTGGACAATATGCCACTCCGTGGAAAGCAG
CTGCGTGTGCGCTTTCCTGCCATAGTGCATCCCTTACAGTTCGAAACCTTCTCAGTATGTGTCCAAC
GAACTGCTGGAAGAAGCCTTTTCTGTGTTGGCCAGGTAGAGAGGGCTGTAGTCATTGTGGATGATCGA
GGAAGGCCCTCAGAAAAGGCATTGTTGAGTTCAGGGAAGCCAGCTGCTCGAAAGCTCTGGACAGA
TGCAGTGAAGGCTCCTTCTGTCTAACACATTTCTCGTCTGTGACTGTGGAGCCCATGGACCAGTTA
GATGATGAAGAGGGACTTCCAGAGAAGCTGGTTATAAAAAACCAGCAATTCACAAGGAACGAGAGCAG
CCACCCAGATTTGCACAGCCTGGCTCCTTTGAGTATGAATATGCCATGCCCTGGAAGGCACTCATTGAG
ATGGAGAAGCAGCAGCAGGACCAAGTGGACCGCAACATCAAGGAGGCTCGTGAGAAGCTGGAGATGGAG
ATGGAAGCTGCACGCCATGAGCACCAGGTCATGCTAATGAGACAGGATTTGATGAGGCGCCAAGAAGAA
CTTCGGAGGATGGAAGAGCTGCACAACCAAGAGGTGCAAAAACGAAAGCAACTGGAGCTCAGGCAGGAG
GAAGAGCGCAGGCGCGTGAAGAAGAGATGCGGCGGCAGCAAGAAGAAATGATGCGGCGACAGCAGGAA
GGATTCAGGGAACTTCCCTGATGCGAGAGAGCAGGAGATTCGGATGGGTCAGATGGCTATGGAGGT
GCTATGGGCATAAACAACAGAGGTGCCATGCCCTGCTCCTGTGCCAGCTGGTACCCAGCTCCTCCA
GGACCTGCCACTATGATGCCGGATGGAACCTTGGGATTGACCCCAACAACACTGAACGCTTTGGTCAG
GCTGCTACAATGGAAGGAATTGGGGCAATTGGTGAACCTCCTCTGCATTCAACCGTGCAGCTCCTGGA
GCTGAATTTGCCCAACAAACGTGCGCGATACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001145409

**Insert Size:** 1416 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001145409.1</a></u>
<b>RefSeq Size:</b>	3057 bp
<b>RefSeq ORF:</b>	1416 bp
<b>Locus ID:</b>	4841
<b>UniProt ID:</b>	<u><a href="#">Q15233</a></u>
<b>Cytogenetics:</b>	Xq13.1
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	54.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16. [provided by RefSeq, Feb 2009]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein (isoform 1).</p>