

## Product datasheet for **SC326492**

### ADAD2 (NM\_001145400) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ADAD2 (NM_001145400) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAD2
Synonyms:	TENRL
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001145400, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTTCGGCTTCTCAGGGCGCTGACGACGACGGCAGTCGTAGGAAGCCCCGCCTGGCT
GCATCGTTGCAGATCAGCCCCAGCCCCGCCCTGGCGACCGCTACCCGCCAGGCCAA
AGTGCCTGGGGGCCCGCGCCCGCCCGCAGCTATCGCGCGGAGGGCGGGTGGCCCCAG
GTCTCGGTGTTGAGGGACAGTGGGCCTGGGGCAGGGGCCGGAGTCGGGAACTGGGGCA
GCCCGGGCTGGGAAAACCTTGGGGGAACAGATGGGGAAGGCCCGAGGGTCCCTGTGCC
CCAGCAGGGCTCAGCCTGCCGCTCAAAGACCCACCTGCCAGCCAGGCCGTGTCCTTGCTC
ACGGAGTACGCGGCCAGCCTGGGCATCTTCTGCTCTTCCGGGAGGACCAGCCACCAGGT
CCCTGCTTCCCCTTCTCGGTGAGCGCGAACTGGATGGGGTGGTCTGCCCTGCGGGCACT
GCGAATAGCAAGACGGAGGCCAAACAGCAGGCAGCGCTCTTGCCTCTGCTACATCCGG
AGTCAGCTGGAGAACCAGAGTCCCCCAGACCTCCAGCCGGCCTCCACTGGCCCCCTG
AGCGTAGAGAATCCTGACCCATGAGCAGCGCTGCGCAGCGTTGGTGAAGCGCCGCTTT
GACCTCCTGTTGACGAGCGCTCGCCATACTGGGCCTGTAAAGGGACTGTGGCTGGAGTC
ATCCTGGAGAGGGAGATCCCGCTGCCAGGGGCCACGTGAAGGAGATCTACAAGTGGTG
GCTCTGGGCACCGCAGCAGCTGCTGTGCTGGCTGGCTGGAGTTCTCGGGCCAGCAGCTC
CAGACTGCCATGGCCTGGTCATCGCCCGAGGGCCCTGCTGAGGTTCTTGTCCGGCAG
CTCCTGCTGGCCACACAGGGGGGCCAAAGGGCAAGGAGCAGTCCGTGCTGGCCCCCAG
CCAGGGGCCCGACCCCCATTACCCCTCAAGCCCCGCGTCTTCTGCACCTCTACATCAGC
AACACCCCAAGGGCGCGGCCCGTGACATCTACCTGCCCCCACCTCGGAAGGTGGCCTC
CCGCACAGCCACCCATGCGCCTGCAGGCCATGTGCTCGGGCAGCTGAAGCCTGTGTGC
TACGTGGCGCCCTCGCTCTGTGACACCCACGTGGGCTGCCTGTCAGCCAGTGACAAGCTG
GCACGCTGGGCCGTGCTGGGGCTGGGTGGTGGCCTGCTGGCCACCTGGTGTCCCACTC
TACAGCACCCAGCCTCATCCTGGCTGACTCATGCCACGACCTCCGACTCTGAGCAGGGCC
ATCCACACCCGGCCCTGCCTGGACAGTGTCTGGGGCCATGCCTGCCACCTCCCTACGTC
CGGACCGCCTGCACCTGTTTGAGGGCCCCGGTGGCCCTTCCGAACCCACCCCTGAC
ACCTGCCGTGGCCTGAGCCTCAACTGGAGCCTGGGGGACCCTGGCATCGAGGTTGTGGAT
GTGGCCACCGGGCGTGTGAAGGCAATGCCGCCCTGGGGCCTCCCTCCCGTCTCTGCAAG
GCCTCCTTTCTCCGGCCTTTACCAGGGCGCCAGGGCTGTGGGGAAGCCCTACCTCCTG
GCCTTGAAGACCTACGAGGCTGCCAAGGCTGGGCCCTACCAGGAGCTCGCAGGCAGCTG
TCTCTCCTCTGGACCAGCAGGGCCTGGGGCTTGGCCCTCGAAGCCACTGGTGGGCAA
TTCAGAAAC

```

- Restriction Sites:** Please inquire
- ACCN:** NM\_001145400
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001145400.1](#), [NP\\_001138872.1](#)

**RefSeq Size:** 2076 bp

**RefSeq ORF:** 1752 bp

**Locus ID:** 161931

**UniProt ID:** [Q8NCV1](#)

**Cytogenetics:** 16q24.1