

Product datasheet for **SC326561**

ZDHC7 (NM_001145548) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHC7 (NM_001145548) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZDHC7
Synonyms:	DHC7; SERZ-B; SERZ1; ZNF370
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Restriction Sites:	SgfI-MluI
ACCN:	NM_001145548
Insert Size:	1038 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145548.1</u>


[View online »](#)

RefSeq Size: 3279 bp

RefSeq ORF: 1038 bp

Locus ID: 55625

UniProt ID: [Q9NXF8](#)

Cytogenetics: 16q24.1

Protein Families: Transmembrane

MW: 38.9 kDa

Gene Summary: Palmitoyltransferase with broad specificity (PubMed:22031296, PubMed:28196865). Palmitoylates JAM3 (PubMed:28196865). Palmitoylates SNAP25 and DLG4/PSD95 (By similarity). Palmitoylates sex steroid hormone receptors, including ESR1, PGR and AR, thereby regulating their targeting to the plasma membrane and their function in rapid intracellular signaling upon binding of sex hormones (PubMed:22031296). May play a role in follicle stimulation hormone (FSH) activation of testicular Sertoli cells (By similarity).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.