

Product datasheet for **SC326579**

RFX8 (NM_001145664) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RFX8 (NM_001145664) Human Untagged Clone
Tag:	Tag Free
Symbol:	RFX8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326579 representing NM_001145664. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTATGAGATTTACGTGGAGACCTGTGGGCAAAACACTGAGAACCAAGTCAACCCGGCCACCTTTGGG
AAGTGTGAAGATCATTACCCGATGAAGACAGACCCAGTTGGATCCCTTTGTCTGAATTCAGGAGATGT
CCATTTCTGGAGCAAGAAGTGGCAAGAAATACTCCTGTAATATGATGGCCTTCTTGCTGACGAATAC
TGCAACTATTGTCGAGACATTTACGAAATGTGGAGGACTTGCTACTTCTTCTGGAAGTCTCTGCAG
CAAGACACAGTCATGCTGATGTCATTGCCTGACGTGTGCCAGCTTTAAATGCTACGACGTCCAGCTG
TACAAGGAATTGAGGATGTTCTCCTTCCATGACTTCTTGAAGATGTTTCTATTACCTGAACTCT
GTGCAGTTATTTAGTAAGAAATTTAAGCTGTGGCTCCTTAATGCTTTGGAAGGTGTTCCAGCCCTTTG
CAGATCTCCAAACTCAAAGAAGTTACTCTATTTGTCAAAAGACTAAGAAGAAAGACGTACCTTTCCAAC
ATGGCTAAGACTATGCGAATGGTATTGAAAAGTAAGAGGCGTGTGAGCGTTTTGAAGTCAGATCTACAG
GCCATCATCAATCAAGGCACTTTGGCTACTTCTAAGAAAGCCCTGGCAAGTGACCGGAGTGGCGCAGAT
GAACTGGAGAACAACCCAGAGATGAAATGTTTAAGAACTTAATTTCTTGTGGAAACATCAACAGAT
CTCAGGGTATTCCTCAGCTGTCTGTCTTACATCTCCAAGCATTGTGTTCCAGACAAGCAGAAGCAAA
GAAGAGTTTACCAAATTGGCCGACGCTCCAGCTGAGATGGAATCTTCTCTCACTGCTGAAGCAAA
GCCATGACCTCTGCCACAGAGATAGTTTTGGCTCCTGGCATCTGTTTCACTTGTGCTTTTGAATAT
ATGATTATATACTTCAGTCATGCCTAGAGGAGGAAGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGG
ATGCTACCAGATGACCCGACTCTCGGCCAGCCAGACCAGGCACTTTTCCATTCTCTGAATTCCTCACTG
TCGAGGCGTGTGCCAGCCCAGCATGGAGCCACTGGGGGTGATGCCACACACATGGGCCAGGGCCGA
TATCCCGTGGGTGTGAGCAACATGGTCCTCAGGATCCTGGGCTTCTGCTGGACTGCCATGGGCAAT
AAGCTCATCCAGGTGCTGTTGGAAGATGAAACCACTGAAAGCGCAGTTAAACTCAGCCTTCTATGGGA
CAAGAAGCCCTCATAACCCTAAAAGATGGACAACAATTTGTGATTTCAGATATCAGATGTACCCCAAAGC
TCTGAAGATATTTATTTTCAGAGAAAACAATGCTAATGTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001145664

Insert Size: 1422 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145664.1](#)

RefSeq Size: 1729 bp

RefSeq ORF: 1422 bp

Locus ID: 731220

UniProt ID: [Q6ZV50](#)

Cytogenetics: 2q11.2
MW: 53.7 kDa
Gene Summary: May be a transcription factor.[UniProtKB/Swiss-Prot Function]