

Product datasheet for **SC326584**

ZBTB8B (NM_001145720) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZBTB8B (NM_001145720) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZBTB8B
Synonyms:	ZNF916B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326584 representing NM_001145720.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGATGCAATCCTATTATGCCAAGCTTTTGGGGGAGCTGAATGAACAGAGAAAGAGGGACTTTTTC
TGTGACTGCAGCATCATTGTGGAAGGGCGGATCTTCAAGGCCACAGGAACATTTGTTTGCTAACAGC
GGCTACTTCCGAGCCCTGCTCATTCACTATATCCAGGACAGCGGGCGGCATAGCACCGCCTCCTTGGAC
ATTGTACCTCTGATGCCTTCTCCATCATCTTAGATTTCTCTATTCTGGGAAGTTGGATTGTGTGGG
GAGAATGTGATTGAAGTGATGTCGGCTGCCAGCTACCTGCAGATGAATGACGTGGTGAACCTCTGCAAG
ACATACATTAGGTCATCCCTCGACATTTGCCGAAAGATGGAGAAGGAGGCTGCTGTGGCTGCAGCAGTG
GCGGCGGCAGCGCGCGGCTGCAGCGGCGGCAGCAGCGGCGGCTCATCAGGTTGACAGTGAAAGCCCC
AGTTCAGGCCGGAGGGACCTCTGTGGTACCAAGAGCTTGGTCTCTCTCCAGCCGAGGGAGAAAAG
AGCGTGGAGTGCCTGAGAGAGTCCCCTTGCAGTACTGCGGAGACTGCCACCCCTTGAAGTGGTGGTG
AGAGACAGCCTTGGCGGTGGCTCGGCTGACAGCAACCTCTCTACTCCACCAAACGGATAGAGCCCAAG
GTGGAATTTGATGCTGATGAAGTGAGGTGGACGTTGGTGAACAGCTGCAGCAGTATGCTGCCCCGCTG
AACCTGGCCACGTGGAGGAGGCTTGGCAAGCGGCCAGGCGGTTGACTTGGCTTACAGCAACTACCAC
GTGAAGCAGTTCTCGAGGCGCTCTTGCACAACAGCGCTGCCCGAGCAAGGATGATGCAGACCATCAC
TTTTCTAGGAGTTTGAAGGAAGACCAGAAGGTGCAGGAGTAGCCATGAGTTCCATGATGGATGTCCAG
GCTGACTGGTATGGAGAGGACTCAGGTGATGTGCTGGTGGTCCCATCAAGCTCCACAAGTGTCTTTTC
TGCCCTTACACTGCCAAACAGAAGGGCATCCTCAAGCGGCACATCCGTTACACACAGGCGAGCGGCC
TACCCCTGTGAGACCTGCGGCAAGAGGTTCACTCGACAAGAGCACCTGCGGAGCCACGCACTGAGTGTC
CACAGATCCAACCGTCCAATCATCTGCAAGGGCTGCAGGAGAACATTACGAGTCACTGTCCCAGGGG
CTGCGGGCGCTTCGGGCTGTGTGACAGCTGCACCTGCGTTACAGACACACCCGATGATGATGATTTG
ATGCCCATCAACCTTAGCTTGGTGGAGGCTTCTGAAAGCCAAGAAAAGAGCGACACAGACAATGAC
TGGCCAATCTATGTGGAGTCGGGTGAGGAAAATGACCCTGCTGGAGATGATTCTGATGACAAACCACAA
ATTCAGCCTAACTTATCAGACCGAGAGACATTACGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001145720

Insert Size: 1488 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145720.1</u>
RefSeq Size:	4019 bp
RefSeq ORF:	1488 bp
Locus ID:	728116
UniProt ID:	<u>Q8NAP8</u>
Cytogenetics:	1p35.1
MW:	54.2 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]