

Product datasheet for **SC326592**

ZNF447 (ZSCAN18) (NM_001145542) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF447 (ZSCAN18) (NM_001145542) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN18
Synonyms:	ZNF447
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326592 representing NM_001145542.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCGGAGAGTCGGAGGGCGAGCTGACGAGGACGCTACGGCGGCCGGAAGTCTGTGGTTCTGGCACCA
CCTGAGCCCACTGGGCATCTGGTCATCCCTGGCACCTCTCCTTTGGAGCCACCTGTCCCTGGCTAGAC
AGTCACATTTTCCAGTGCCGTTTTGGAAAAGATGTTGCCTTTGGAGAAGGGGTTTGCTCCCCAGGAGC
TCCCCAGCCCCGCCGGATCTGCCACGCCGGGTGAGCAGCCGGAGTCCAGCAGGAAGAACCCGAGACC
ATCCCTGAGAGGACCCCTGCTGACCTGGAGTTCTCCCGCCTGCGTTTCCGGGAATTTGTCTACCAGGAG
GCTGCCGGGCCCCACCAGACCCTGGCCCGGCTGCATGAGCTGTGCCGCCAGTGGCTGATGCCTGAGGCG
CGCTCAAAGGAGCAGATGCTGGAGCTGCTGGTGTGGAGCAGTTCTGGGCATCTGCCTGATAAGGTC
CGGCCCTGGTGGTGGCACAGTACCCTGAGAGCTGCAAGAAGGCAGCCTCCCTGGTGGAGGGCTCGCT
GATGTCTGGAAGAGCCAGGGATGCTGCTGGGCTCCCTGCGGGCTCATCTCAATCTTAGCGATGGA
GTGTACGAGAGGCACATGGACCCTGTGCTACCAGGCGAGCTCGCAGCCCCAGCCAGGCCCTTGA
GCTGGGGAGATCCCGGCACCTTCTGAGACACCCTGGCTTTCTCCGGACCCCTGTTTCTGGAACAGAGG
AGGGTCAGAGAAGCAAAGACCGAAGAGGACGGCCCTGCCAACCCGAGCAGAAGCTGAAGTCTTTCCA
GAGGACCTCAGCACCTGGGGGAGTGGGGCCACCTGGACCCTGCCGAGGAGAACCTGAAGAGTACCGG
AAGCTGTCTGTGGGGTATCAGCTTTCCAGCCTGACGCTGCCTCCAGGCTGGACACTGAGGAACTC
CGGTTGGTGGAAAGAGATCCACAAGGAAGCAGCCTCCAGAAGGCGGGAGGCGGCAGGAGAGCGTGGG
TGCGCCTGCGAGGAGGCCGCCCGCGGGGTGCTGCCTGAGCTGCCTACGGAGGCGCCCCTGGGGAC
GCCCTTGCCGATCCCCGTCGGGCACCACTGAGGAGGAGGAAGAGCAGCCTGGGAAGGCCCGGACCCG
CAGGACCCCAAGGACGCGGAGTCCGACTCTGCCACCGGATCGCAGAGGCAGTCCGTATCCAGCAGCCT
GCCCGGACAGGGGCACGGCGAAACTGGGAACCAAGAGGCCGCACCCGAGGATGGGGACGGGCAGAGC
CTCGAGGGCGTCTCTAGCTCCGGCGACAGCGCAGGGCTGGAGGCCGGGCAGGGCCCTGGGCTGACGAG
CCGGGCTTGTCGCCGGGAAGCCCTATGCCTGCGGCGAGTGGGGGAGGCCTTCGCGTGGCTCTCGCAC
CTGATGGAGACCACAGCAGCCATGGCGGCCGGAAGCGCTACGCCTGTCAGGGCTGCTGGAAGACCTTC
CACTTCAGCCTGGCCCTAGCCGAGCACCAGAAGACCCACGAGAAGGAGAAAAGCTACGCGCTGGGGGC
GCCCGGGCCCCCAACCGTCCACCCGGAAGCCAGGCGGGGGCTAGGGCGGGCGGTCCCCCAGAGAGC
GTGGAGGGCGAGGCTCCCCCGCACCCCAAGGGCGCAGAGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001145542

Insert Size: 1701 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145542.1](#)

RefSeq Size: 2605 bp

RefSeq ORF: 1701 bp

Locus ID: 65982

UniProt ID: [Q8TBC5](#)

Cytogenetics: 19q13.43

Protein Families: Transcription Factors

MW: 60.8 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longest protein (isoform 1).