

## Product datasheet for **SC326714**

### CRCP (NM\_001142414) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CRCP (NM_001142414) Human Untagged Clone
Tag:	Tag Free
Symbol:	CRCP
Synonyms:	C17; CGRP-RCP; CGRPRCP; POLR3I; POLR3J; RCP; RCP9; RPC9
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001142414, the custom clone sequence may differ by one or more nucleotides ATGCAGCGCGGATCCCGCGGAGCACCTTGGCGCGGGAGCTGGCACCTTGGCGTGT GGTGGCGCGGAGACAGCTGTGAAGTGAAGGATGCCAATTCTGCGCTTCTCAGTAACTAC GAGACGTTAAAATACATATCAAAAACACCATGCAGGCACAGAGTCTGAAATTGCAGA GAATTTCTCACAGCATTGAAAAGCCACAAGTTGACCAAAGCTGAGAAGCTCCAGTGCTG AACCAACCGGCTGTGACTGCTGTGGAGATCCAGCTGATGGTGAAGAGAGTGAAGAGCGG CTCACGGAGGAGCAGATTGAAGCTCTTCTCCACACCGTCACCAGCATTCTGCCTGCAGAG CCAGAGGCTGAGCAGAAGAAGAATACAAACAGCAATGTGGCAATGGACGAAGAGGACCCA GCA
Restriction Sites:	Please inquire
ACCN:	NM_001142414
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001142414.1](#), [NP\\_001135886.1](#)

**RefSeq Size:** 2648 bp

**RefSeq ORF:** 426 bp

**Locus ID:** 27297

**UniProt ID:** [O75575](#)

**Cytogenetics:** 7q11.21

**Protein Families:** Druggable Genome

**Gene Summary:** This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (4) lacks an internal segment in the 5' region, which results in an upstream AUG start codon, and also lacks an exon in the middle region, as compared to variant 1. The reading frame is not affected, and the resulting isoform (d) is shorter and has a distinct N-terminus, as compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.