

## Product datasheet for **SC326765**

### **SFTPA1 (NM\_001164646) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SFTPA1 (NM_001164646) Human Untagged Clone
Tag:	Tag Free
Symbol:	SFTPA1
Synonyms:	COLEC4; PSAP; PSP-A; PSPA; SFTP1; SFTPA1B; SP-A; SP-A1; SP-A1 beta; SP-A1 delta; SP-A1 epsilon; SP-A1 gamma; SPA; SPA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326765 representing NM_001164646. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**  
 ATGTGGCTGTGCCCTCTGGCCCTCAACCTCATCTTGATGGCAGCCTCTGGTGCTGTGTGCCAAGTGAAG  
 GACGTTTGTGTTGGAACCCCTGGTATCCCTGGAGAGTGTGGAGAGAAGGGGGAGCCTGGCGAGAGGGGC  
 CCTCCAGGGCTTCCAGCTCATCTAGATGAGGAGCTCCAAGCCACACTCCACGACTTTAGACATCAAATC  
 CTGCAGACAAGGGGAGCCCTCAGTCTGCAGGGCTCCATAATGACAGTAGGAGAGAAGGTCTTCTCCAGC  
 AATGGGCAGTCCATCACTTTTGATGCCATTGAGGAGCATGTGCCAGAGCAGGCGGCCGATTGCTGTC  
 CCAAGGAATCCAGAGGAAATGAGGCCATTGCAAGCTTCGTGAAGAAGTACAACACATATGCCTATGTA  
 GGCCTGACTGAGGGTCCCAGCCCTGGAGACTTCCGCTACTCAGACGGGACCCCTGTAACTACACCAAC  
 TGGTACCAGGGGAGCCCGAGGTCGGGAAAAGAGCAGTGTGTGGAGATGTACACAGATGGGCAGTGG  
 AATGACAGGAAGTGCCTGTACTCCCGACTGACCATCTGTGAGTTC**TGA**  
**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001164646
Insert Size:	600 bp


[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001164646.1</a></u>
<b>RefSeq Size:</b>	2083 bp
<b>RefSeq ORF:</b>	600 bp
<b>Locus ID:</b>	653509
<b>Cytogenetics:</b>	10q22.3
<b>MW:</b>	21.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a lung surfactant protein that is a member of a subfamily of C-type lectins called collectins. The encoded protein binds specific carbohydrate moieties found on lipids and on the surface of microorganisms. This protein plays an essential role in surfactant homeostasis and in the defense against respiratory pathogens. Mutations in this gene are associated with idiopathic pulmonary fibrosis. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (6) lacks an in-frame segment in the central coding region, compared to variant 1, resulting in an isoform (4) that is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>