

Product datasheet for **SC326853**

VSTM2B (NM_001146339) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VSTM2B (NM_001146339) Human Untagged Clone
Tag:	Tag Free
Symbol:	VSTM2B
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001146339, the custom clone sequence may differ by one or more nucleotides

ATGGAACAGCGGAACCGCTCGGTGCCCTCGGATACCTGCCGCCTCTGCTGCTGCATGCCCTGCTGCTCT
 TCGTGGCCGACGCTGCATTACAGAAAGTCCCCAAAGATGTGACAGTACGGGAGGGAGACGACATCGAAAT
 GCCCTGCGCGTTCCGGGCCAGCGGAGCCACCTCGTATTCGCTGGAGATTAGTGGTGGTACCTCAAGGAG
 CCACCCCGGGAGCTGCTGCACGAGCTGGCGCTCAGCGTGCCGGGCGCCCGGAGCAAGGTAACAAATAAGG
 ATGCAACTAAAATCAGCACCGTACGCGTCCAGGGCAATGACATCTCACACCGGCTTCGGCTGTCTGCCGT
 GCGGCTGCAGGACGAGGGCGTGTACGAGTGCCGCGTGTGCGACTACAGCGACGACGACACGACGAGGACAC
 AAGGCCAGGCGATGCTGCGCGTGTCTCGCGCTTCGCGCCGCCAACATGCAGGCCGCCGAGGCCGTGT
 CCCACATCCAGAGCAGCGGCCCGCTCGCCACGGCCAGCCAGCGCCGCCAACGCCAACACGCGGGCGC
 CGCGAGCCGTACCACCTCCGAGCCCGGCCGCGGCGACAAGAGCCCGCCGCCCGGGAGCCCTCCCGCGGCC
 ATCGATCCCGCAGTCCCGAGGCCGCGGCAGCCTCGGCGGCCACACGCCACCACCACAGTCGCGGCAG
 CTGCTGCTGCCTCGTCAGCGTCGCCGCATCGGGACAGGCGGTCTGCTGCGCCAGAGGCACGGCTCGGG
 CACCGGCCGTAGCTACACCACAGACCCACTCTTGTCCCTGCTCCTGTAGCTCTGCATAAGTTCTGCGC
 CTGCTCTTGGGACATTGA

Restriction Sites:	Please inquire
ACCN:	NM_001146339
Insert Size:	1700 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001146339.1, NP_001139811.1</u>
RefSeq Size:	1243 bp
RefSeq ORF:	858 bp
Locus ID:	342865
UniProt ID:	<u>A6NLU5</u>
Cytogenetics:	19q12
Protein Families:	Transmembrane