

Product datasheet for SC326864

Sex Hormone Binding Globulin (SHBG) (NM_001146280) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Sex Hormone Binding Globulin (SHBG) (NM_001146280) Human Untagged Clone
Tag:	Tag Free
Symbol:	SHBG
Synonyms:	ABP; SBP; TEBG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326864 representing NM_001146280. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGAGCAGAGGCCACTGGCTACCTCGCGCTGCTGCTGTTGCTGCTGTTGCTACTACTGCGTCAC
ACCCGCCAGGGATGGGCCCTGAGACCTGTTCTCCACCAGAGTGCCACAGACCCTCCGGCTGTCCAC
CTCAGCAATGGCCAGGACAAGACCTATCGCTGTCATGACCTTTGACCTACCAAGATCACAAAAACC
TCCTCCTCCTTTGAGGTTGCAACCTGGGACCCAGAGGGAGTGATTTTTATGGGGATACCAACCCTAAG
GATGACTGGTTTATGCTGGGACTTCGAGACGGCAGGCCTGAGATCCAACGCACAATCACTGGGCCAG
CTTACGGTGGGTGCTGGACCACGGCTGGATGATGGGAGATGGCACCAGGTGGAAGTCAAGATGGAGGGG
GACTCTGTGCTGCTGGAGGTGGATGGGGAGGAGGTGCTGCGCCTGAGACAGGTCTCTGGGCCCTGACC
AGCAAACGCCATCCCATCATGAGGATTGCGCTTGGGGGGCTGCTTCCCGCTTCCAACCTTCGGTTG
CCGCTGGTTTCTGCCCTGGATGGCTGCCTGCGCCGGGATTCCTGGCTGGACAAACAGGCCGAGATCTCA
GCATCTGCCCCACTAGCCTCAGAAGCTGTGATGTAGAATCAAATCCCGGGATATTTCTCCCTCCAGGG
ACTCAGGCAGAATTCAATCTCCGAGACATTCCCAGCCTCATGCAGAGCCCTGGGCCCTCTCTTTGGAC
CTGGGACTCAAGCAGGCAGCAGGCTCAGGCCACCTCCTTGCTCTGGGACACCAGAGAACCATCTTGG
CTCAGTCTCCACCTCCAAGATCAAGAGAAGACTTCCACCTCTTTTGCCTGA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001146280

Insert Size: 882 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146280.2](#)

RefSeq Size: 1125 bp

RefSeq ORF: 882 bp

Locus ID:	6462
UniProt ID:	P04278
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Secreted Protein
MW:	32.4 kDa
Gene Summary:	<p>This gene encodes a steroid binding protein that was first described as a plasma protein secreted by the liver but is now thought to participate in the regulation of steroid responses. The encoded protein transports androgens and estrogens in the blood, binding each steroid molecule as a dimer formed from identical or nearly identical monomers. Polymorphisms in this gene have been associated with polycystic ovary syndrome and type 2 diabetes mellitus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (3) lacks an alternate exon in the 3' coding region which results in a frameshift compared to variant 1. Isoform 3 is shorter and has a distinct C-terminus compared to isoform 1.</p>