

Product datasheet for **SC326868**

ZNF705G (NM_001164457) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF705G (NM_001164457) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF705G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326868 representing NM_001164457. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGCATTCCTAAAGAACTGACTTTTGAAGATGTAGCTATTGACTTCACCCAGGAAGAGTGGGCCATG
 ATGGACACATCCAAGAGAAAGCTGTACAGAGATGTGATGCTGGAAAATATCAGTCACCTGGTGTCCCTC
 GGGTACCAGATAAGCAAATCCTATATAATTTTGCAGCTGGAGCAAGGAAAGAGCTGTGGAGGGAAGGA
 AGAGTATTTCTTCAAGACCAGAAATCCAAACAGGGAAAGTGCCCTTAAGAAAACACACATGATATCCATG
 CATCCTATCACCAGAAAAGACGCATCCACCAAGTATGACAATGGAGAACTCTCTATTCTGGAGGATCCT
 TTTGAATGTAATGATTCGGGAGAAGATTGCACTCGCAGTTCCACAATAACTCAGTGTGTTGTTAACTCAT
 AGTGGAAGAAACCCCTATGTGCAAGAACAGTGTGGAATCCCTTCGTAATCTTTGTCCACTGAACCA
 CATAAACAAATTCATACTAAAGGTAAATCATATCAGTGAATCTATGTGAAAAGGCCTATACTAATTGC
 TTTCACCTTAGACGGCACAAGATGACTCACACTGGAGAGAGGCCATATGCATGTCATCTATGTAGAAAA
 GCCTTCACTCAGTGTCTCACCTTAGAAGACACGAGAAAACCTCACACGGGACAGAGACCATATAAGTGT
 CATCAATATGGGAAAGTCTTTATTCAATCCTTTAACCTTCAAAGACATGAGAGAACTCACCTTGGAAAA
 AAGTGTATGAATGTGATAAAAGTGGGAAAGCCTTTAGTCAAAGCTCTGGCTTTAGAGGAAACAAAATA
 ATTCACACTGGAGAGAAACCACATGCTTGTCTTCTATGTGGAAGGCCTTCAGTCTGTCTCCAACCTT
 AGATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001164457
Insert Size:	903 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164457.1</u>
RefSeq Size:	925 bp
RefSeq ORF:	903 bp
Locus ID:	100131980
UniProt ID:	<u>A8MUZ8</u>
Cytogenetics:	8p23.1
MW:	34.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]