

Product datasheet for **SC326892**

ADAL (NM_001159280) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ADAL (NM_001159280) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAL
Synonyms:	HsMAPDA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326892 representing NM_001159280. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATAGAGGCAGAAGAGCAACAGCCTTGAAGACAGACTTCTATTCTGAATTGCCAAAAGTGGAACTT
CATGCCCACTTGAATGGATCCATTAGTTCTCATACCATGAAGAAATTAATAGCCAGAAGCCAGATCTT
AAAATCCACGATCAGATGACTGTGATTGACAAGGGAAAGAAAAGAAGCTTTGAAGAATGTTCCAGATG
TTTCAAATTCATCAGCTTACTAGTAGCCCTGAAGATATTCTAATGGTCACAAAAGATGTCATAAAA
GAATTTGCAGATGACGGCGTCAAGTACCTGGAATAAGGAGCACACCCAGAAGAGAAAATGCTACTGGA
ATGACTAAAAGACTTATGTGGAATCTATACTTGAAGGTATAAACAGTCCAACAAGAAAATTGGAC
ATTGATGTTAGGTATTTGATAGCAGTTGACAGAAGAGGTGGCCCTTTAGTAGCCAAGGAGACTGTAAAA
CTTGCCGAGGAGTTCTTCTTTCTACTGAGGGTACAGTCTTGGCCTTGACCTCAGTGGAGACCCTACT
ATTCAAAACAAAAAAGAAAACACAAATACTCCTGGATCTGCTTCTGACAGAATCGGGCATGGAACA
TTTCTCAACTCCGGTGAGGGAGGATCCCTGGATCTGGTGGACTTTGTGAGGCAACATCGGATACCACTG
GAACTCTGTTTGACCTCAAACGTCAAAGTCAAGCAGTCCATCTTATGACCAGCACCATTTCCGATTC
TGGTACAGCATTGCCATCCTTCTGTGATCTGTACTGATGATAAGGGTGTGTTTGAACACACCTTTCT
CAAGAGTACCAGCTGGCAGCTGAAACATTTAATTTGACCCAGTCTCAGGTGTGGGATCTGTCTTATGAA
TCCATCAACTACATCTTTGCTTCTGACAGCACCAGATCTGAACTGAGGAAGAAAATGGAATCACCTGAAG
CCAGAGTGTTACATATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001159280

Insert Size: 987 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001159280.2](https://www.ncbi.nlm.nih.gov/nuccore/NM_001159280.2)

RefSeq Size: 4489 bp

RefSeq ORF: 987 bp

Locus ID: 161823

UniProt ID:	<u>Q6DHV7</u>
Cytogenetics:	15q15.3
Protein Families:	Druggable Genome
MW:	37.4 kDa
Gene Summary:	Putative nucleoside deaminase. May catalyze the hydrolytic deamination of adenosine or some similar substrate and play a role in purine metabolism (By similarity).[UniProtKB/Swiss-Prot Function]