

## Product datasheet for **SC326906**

### MANSC4 (NM\_001146221) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MANSC4 (NM_001146221) Human Untagged Clone
Tag:	Tag Free
Symbol:	MANSC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326906 representing NM_001146221. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCATGTGGCAGAGGTAGCAGTGAACGTGATATTGCTCCTAAGCATGGGGTGGACATCAGACTCTCTC
TGCTCACCCACAATTTTTACAGAGACTGCTGGATCCGTCGCTTCCCAGGTCTTCTAATCAATCTGGAG
GAGTCTCAGAAGCTGGGAGCCAGTTCTTGAAGTATTATCTGAAAGCACTGGCCAGAAGTGCAGTAGG
AGCTGTGTCTTCGGAAGGATGTTTCTGTAACCTGGCTGTCTTACCACAGTCCTATTCATGACAAT
ATCAACTGCCTCCATGTTCACTGCCAACACTGGAGAGCTGCATATTAGAGCCTGGAACCAAGTCCATT
TTGTACAACATAACAGACGGTATAGATCCGGATTTGCTGGTTTTGAACAATCTCCACATATCTAAAT
ACTCGTTTTCATCCAATAGATGGGACAGACTAAGGATTCTAAAAGCTATGAATTTAGATAAAACAACC
ACCACGATAAATGGTATGCTGCCATCCACAGAGGCTCCATCCTCAACCACGCATCAAGATTTGGTTGTA
AACACAAACAGTACCAGTTATTCTAAGGAATTAACCACAGATTTTTGGGCAAGATTTACTTCCCTGAAT
GAGTCCATTACCACAAAGATAAATAAGGTGTCACCAAGTACTGATTTTCATCAGCAATCCAGATAAAG
ACTATTTCTCCTTTCTTTGAACCCATAGACACAAAACCTTCTCATATGCCTGTTCCACCTGGACTCAAC
AGTAGCAAACAATTAATAACAAAACAAAGGATACAATAGCAGAAACCACACATCTGCAAATGAAGAT
GAGGTATCTGTGACTTCAAAGACTTGGCTGGTTTCTGTGGCCCTTGCACCTCTGTATCTTTCTCGGC
TGTGTATAGTCATCCTGGCATCTGGATGCTGTGAAAGCAGCAGGGCCAGTATAAACCAGGACAGAGA
AAATCAGGATCCTTGCAATTAACAAACCGTAACCATATGAAGGAGAACTCTTCATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001146221
Insert Size:	1023 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001146221.1</a></u>
<b>RefSeq Size:</b>	1095 bp
<b>RefSeq ORF:</b>	1023 bp
<b>Locus ID:</b>	100287284
<b>UniProt ID:</b>	<u><a href="#">A6NHS7</a></u>
<b>Cytogenetics:</b>	12p11.22
<b>MW:</b>	37.9 kDa