

Product datasheet for SC326911

Aminomethyltransferase (AMT) (NM_001164711) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aminomethyltransferase (AMT) (NM_001164711) Human Untagged Clone
Tag:	Tag Free
Symbol:	Aminomethyltransferase
Synonyms:	GCE; GCST; GCVT; NKH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326911 representing NM_001164711. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGAGGGCTGTAAGTGTGGTGGCCCGTCTGGGCTTTCGCCTGCAGGCATTCCCCCGGCCTTGTGT
CGTCCACTTAGTTGCGCACAGACCAAGATACTTGGTAGTGACCGGGTGAAGCTGATGGAGAGTCTAGTG
GTTGGAGACATTGCAGAGCTAAGACCAACCAGGGGACTGTCGCTGTTACCAACGAGGCTGGAGGC
ATCTTAGATGACTTGATTGTAACCAATACTTCTGAGGGCCACCTGTATGTGGTGTCCAACGCTGGCTGC
TGGGAGAAAGATTTGCCCTCATGCAGGACAAGGTCAGGGAGCTTCAGAACCAGGGCAGAGATGTGGGC
CTGGAGGTGTTGGATAATGCCCTGCTAGCTCTGCAAGGCCCACTGCAGCCAGGTAACAGGCCGGC
GTGGCAGATGACCTGAGGAACTGCCCTTCATGACCAGTGTGTGATGGAGGTGTTTGGCGTGTCTGGC
TGCCCGTGCACCCGCTGTGGCTACACAGGAGAGGATGGTGTGGAGATCTCGGTGCCGGTAGCGGGGCA
GTTACCTGGCAACAGCTATTCTGAAAAACCCAGAGGTGAAGCTGGCAGGGCTGGCAGCCAGGGACAGC
CTGCGCCTGGAGGCAGGCCTCTGCCTGTATGGGAATGACATTGATGAACACACTACACCTGTGGAGGGC
AGCCTCAGTTGGACACTGGGGAAGCGCCGCGAGCTGCTATGGACTTCCCTGGAGCCAAGGTCATTGTT
CCCCAGCTGAAGGGCAGGGTGCAGCGGAGGCGTGTGGGTTGATGTGTGAGGGGGCCCCATGCGGGCA
CACAGTCCCATCCTGAACATGGAGGGTACCAAGATTGGTACTGTGACTAGTGGCTGCCCTCCCCCTCT
CTGAAGAAGAATGTGGCGATGGGTTATGTGCCCTGCGAGTACAGTCGTCAGGGACAATGCTGCTGGTA
GAGGTGCGGGCGGAAGCAGCAGATGGCTGTAGTCAGCAAGATGCCCTTTGTGCCCAAACTACTATACC
CTCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001164711



[View online »](#)

Insert Size:	1044 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164711.1
RefSeq Size:	2047 bp
RefSeq ORF:	1044 bp
Locus ID:	275
UniProt ID:	P48728
Cytogenetics:	3p21.31
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Nitrogen metabolism, One carbon pool by folate
MW:	37.4 kDa
Gene Summary:	This gene encodes one of four critical components of the glycine cleavage system. Mutations in this gene have been associated with glycine encephalopathy. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011] Transcript Variant: This variant (3) lacks an in-frame exon compared to variant 1. The resulting protein (isoform 3) is shorter compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.