

Product datasheet for SC326920

DISC1 (NM_001164555) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	DISC1 (NM_001164555) Human Untagged Clone
Tag:	Tag Free
Symbol:	DISC1
Synonyms:	C1orf136; SCZD9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326920 representing NM_001164555. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCAGGCGGGGTCTCAGGGCGCCACGCCGCCGGCGGGCGGCGTGTAGCCACCGCGCAGGC
AGCCGGGATTGCTTACCACCTGCAGCGTCTTTTCGGAGGCGGGCTGGCACGGAGGCCGGGCTACATG
AGAAGCTCGACAGGGCTGGGATCGGGTCTTTCCAGCAGTGGGCACACTGTTCCGGTTCACAGGA
GGGGTGTCTGGCAGGAGTCCCACACTCGGAGTCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGC
CTCTTGGTCCGGAGCCCTGTTTCCAAGAGTGCAGCAGCCCTACTGTGACCTCTGTGAGAGGAACCTCG
GCGCACTTTGGGATTCAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCCT
GGGAGTGTGGTGGCAGCAAGAGTTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGGGAG
GCAGCCTGCAGCGATGGAGCAAGCGTGTCCGGGCGAGCAGGCTCTCTGCCATCAGCAGAGTTGAGTAGC
AACAGTGCAGCCCTGGCTGTGGCCCTGAGGTCCCCCAACCCCTCTGGCTCTCACAGTGCCTTTACC
TCAAGCTTTAGCTTTATTCGGCTCTCGCTTGGCTCTGCCGGGAACTGGAGAAGCAGAAGGCTGCCCA
CCATCCAGAGAGGCTGAGTCCCATTGCCAGAGCCCCAGGAGATGGGAGCCAAAGCTGCCAGCTTGGAC
GGCCCTCAGGAGACCCGCGATGTCTCTCGGCCCTTCACTCTCTTGGCTACACGGGTCTCTGCAGAC
TTGGCCAGGCCCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTTTACCAGACATGGACCCT
GGCTCCTCCAGTCTCTGGATCCCTCACTGGCTGGCTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT
GCCCACTTTGGGACACCCTGCTCAGGAAATGGGAGCCAGTGTCTCGGGACTGCCTGCTGAGAAACCGG
AGGCAGATGGAGCTGAGACGTTACAACAAGATTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001164555

Insert Size: 1071 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164555.1](#)

RefSeq Size: 1794 bp

RefSeq ORF: 1071 bp

Locus ID: 27185

UniProt ID: [Q9NRI5](#)

Cytogenetics: 1q42.2

MW: 37.2 kDa

Gene Summary: This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (r) lacks an internal exon and multiple 3' exons but has an alternate 3' segment, as compared to variant L. The resulting isoform (r, also known as isoform 42) is much shorter and has a distinct C-terminus, as compared to isoform L.