

Product datasheet for SC326952

NPTN (NM_001161363) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	NPTN (NM_001161363) Human Untagged Clone
Tag:	Tag Free
Symbol:	NPTN
Synonyms:	GP55; GP65; np55; np65; SDFR1; SDR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326952 representing NM_001161363. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCGGTTCGTGCTGCCAGCGCCCTGGCCCTCTCGCTGTTGCTGGTCTCTGGCTCCCTCCTCCCA
GGGCCAGGCGCCGCTCAGAACGCTGGGTTTGTCAAGTCGCCCATGTCAGAACTAAGCTCACGGGGAC
GCCTTTGAGCTGTACTGTGACGTGGTGGGAGCCCCACGCCAGAGATCCAGTGGTGGTACGCAGAAGTC
AACCGGGCAGAGTCTTTAGACAGCTGTGGGACGGTGTCTCGAAGCGCCGTGTACCCTAAACCCGCC
TACGGGTCAAACGGCGTGAGTGTGCTGAGAATAACCCGGCTCACCTTGGAGGACTCTGGGACTTACGAG
TGCAGGGCCAGCAACGACCCCAAGAGGAATGACTTGAGGCAAAACCCCTCCATAACATGGATTTCGAGCC
CAGGCCACCATAAGCGTCTTTCAGAAAGCAAGGATTGTCACCAGTGAAGAGGTCATTATTCGAGACAGC
CCTGTTCTCCCTGTACCTGCAGTGTAACTCACCTCCAGCTCTCACACCTTACATACAGCTACTGG
ACAAAGAATGGGGTGGAACTGAGTGCCACTCGTAAGAATGCCAGCAACATGGAGTACAGGATCAATAAG
CCGAGAGCTGAGGATTCAGGCGAATACCACTGCGTATATCACTTTGTCAGCGCTCCTAAAGCAAACGCC
ACCATTTGAAGTAAAGCCGCTCCTGACATCACTGGCCATAAACGGAGTGAGAACAAGAATGAAGGGCAG
GATGCCACTATGTATTGCAAGTCAGTTGGCTACCCCAACCCAGACTGGATATGGCGCAAGAAGGAGAAC
GGGATGCCCATGGACATTGTCAATACCTCTGGCCGTTCTTTCATCATCAACAAGGAAAATTACACTGAG
TTGAACATTGTGAACCTGCAGATCACGGAAGACCCTGGCGAGTATGAATGTAATGCCACCAACGCCATT
GGCTCCGCTCTGTTGCTACTGTCTCAGGGTGGGAGCCACTGGCCCACTCTGGCCTTTCTTGGGA
ATTCTGGCTGAAATTATCATCTTGTGGTGTATGTTGTGTATGAGAAGAGGAAGAGGCCAGATGAG
GTTCTGACGCTGGACCAATGAAAACCACTCTACCAACAATCACAAAGATAAAAACTTGCGCCAGAGA
AACACAAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001161363

Insert Size: 1185 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161363.1](#)

RefSeq Size: 2442 bp

RefSeq ORF: 1185 bp

Locus ID: 27020

UniProt ID: [Q9Y639](#)

Cytogenetics: 15q24.1

Protein Families: Druggable Genome, Transmembrane

MW: 43.9 kDa

Gene Summary: This gene encodes a type I transmembrane protein belonging to the Ig superfamily. The protein is believed to be involved in cell-cell interactions or cell-substrate interactions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2009]
Transcript Variant: This variant (c) uses an alternate in-frame splice site in the 3' coding region compared to variant b. This results in a shorter protein (isoform c), compared to isoform b.