

Product datasheet for **SC326957**

SPDYE6 (NM_001146210) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPDYE6 (NM_001146210) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPDYE6
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001146210, the custom clone sequence may differ by one or more nucleotides ATGGACAGAACGGAGACTAGGTTCCGTAAGAGGGGACAGATTACGGGAAAGATCACGACC AGCCGTCAACCGCACCCCAGAAATGAGCAGAGTCCCCAGCGGAGCACCTCGGGGTACCCC CTCCAGGAGGTGGTGGATGATGAAATGTTGGGACCATCAGCCCCCTGGGGTAGATCCCAGC CCCCATGTAGTCCCTTGGCTGGAAAAGGAAGAGGGAGTGGTCAGATGAATCTGAGGAG GAGCCGGAGAAGGAGCTCGCCCCGAGCCTGAGGAGACCTGGGTAGTGGAGATGCTGTGT GGGCTCAAGATGAAGCTGAAGCAACAGCGAGTGTATCCATCCTCCCTGAGCACCACAAG GACTTCAACAGTCAGCTTGCCCCCTGGGGTAGATCCCAGCCCCCGCATAGGTCCTTTTGC TGGAAAAGGAAGATGGAGTGGTGGGACGAATCTGAGGAGTCGTTGGAGGAGGAGCCACGG AAGGTGCTCGCCCCTGAGCCTGAGGAGATCTGGGTGGCGGAGATGCTGTGTGGCCTCAAG ATGAAGCTGAAGCGACGGCGAGTGTGCTCGTGTCTCCCTGAGCACCACGAGGCCTTCAAC AGGCTGCTTGAGGATCCTGTCTATTAAGATTCTTGGCCTGGGACAAAGATCTGAGGGTG TCGGACAAGTATCTCCTGGCTATGGTCATAGCGTATTTTACGCCGGGCTGGCTTCCCTCC TGGCAATACCAACGCATTCAATTTCTTCTGGCTCTCTACCTGGCCAATGACATGGAGGAG GACGACGAGGACTCCAAACAAAACATCTTCCACTTCTGTATAGGAAGAACCCTCTCGC ATACCCTTGCTCCGTAAGCGTTGGTTCAGTTAGGCCATTCCATGAACCCGAGGGCCAGG AAGAACCCTCTCGCATACCTTGCTCCGTAAGCGTCGGTTCAGTTATACCGTTCCACG AACCCGAGGGCCAGGAAGAACCCTCTCGCATACCTTGCTCCGTAAGCGTCGGTTCAG TTATACCGTTCCATGAATCGAGGGCCAGGAAGAACCCTCTCAGATAGTCTGTTCCAG AAACGACGGTTCCACTTCTTCTGTTCCATGAGCTGCAGGGCTTGGGTTTCCCCAGAGGAG TTGGAGGAGATCCAGGCTTATGACCCAGAGCACTGGGTGTGGGCGCGAGATCGCGCTCAC CTTTCC
Restriction Sites:	Please inquire
ACCN:	NM_001146210



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001146210.1</u> , <u>NP_001139682.1</u>
RefSeq Size:	3145 bp
RefSeq ORF:	1209 bp
Locus ID:	729597
UniProt ID:	<u>P0CI01</u>
Cytogenetics:	7q22.1
Gene Summary:	This gene encodes a cell cycle regulatory protein which plays an important role in cell cycle progression by binding and activating cyclin-dependent kinases. The encoded protein belongs to a family of cyclin-dependent kinase regulators that contain a speedy box domain which is required for cyclin-dependent kinase activation. [provided by RefSeq, Jul 2017]