

Product datasheet for **SC326959**

TICAM2 (TMED7-TICAM2) (NM_001164468) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TICAM2 (TMED7-TICAM2) (NM_001164468) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMED7-TICAM2
Synonyms:	MyD88-4; TICAM-2; TICAM2; TIRAP3; TIRP; TRAM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326959 representing NM_001164468. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGCGCCGGGGTCCGCGCAGCGCTGGGCGGCCGTGCGGGCCGTTGGGGTGCAGGCTGCTCGCA
CTGCTGCTACTGGTGCCTGGACCCGGCGCGCCTCTGAGATCACCTTCGAGCTTCTGACAACGCCAAG
CAGTGCTTCTACGAGGACATCGCTCAGGGCACCAAGTGCACCCTGGAGTTCAGGTGATTACTGGTGGT
CACTATGATGTAGATTGTCGATTAGAAGATCCTGATGGTAAAGTGTATACAAAGAGATGAAGAAACAG
TATGATAGTTTTACCTTCACAGCCTCAAAAATGGGACATACAAATTTTGCTTCAGCAATGAATTTTCT
ACTTTACACATAAAAAGTATATTTTGTATTTCAAGTTGGAGAAGACCCACCTTTGTTTCTAGTGAG
AACCGAGTCAGTGCTCTTACCCAGATGGAATCTGCCTGTGTTTCAATTCACGAAGCTCTGAAGTCTGTC
ATCGATTATCAGACTCATTCCGTTTAAGAGAAGCTCAAGGCCGAAGCCGAGCAGAGGATCTAAATACA
AGAGTGGCCTATTGGCACAGTGTGGATACAAGTCCAGGATATCATGAGTCAGATTCGAAGAAGTCTGAA
GATCTATCCTTGTGTAATGTTGCTGAGCACAGCAATACAACAGAGGGGCCAACAGGAAAGCAGGAGGGA
GCTCAGAGCGTGGAGAGATGTTGAAGAAGAAGCTGAAGAAGAGGTGTTCTCAAAATTTGTGATATTG
CATGCAGAAGATGACACAGATGAAGCCCTCAGAGTCCAGAATCTGCTACAAGATGACTTTGGTATCAAA
CCCGGAATAATCTTTGCTGAGATGCCATGTGGCAGACAGCATTTACAGAATTTAGATGATGCTGTAAT
GGGTCTGCATGGACAATCTTATTACTGACTGAAAACCTTTTTAAGAGATACTTGGTGAATTTCCAGTTC
TATACGTCCCTAATGAACCTCCGTTAACAGGCAGCATAAAACAACCTCTGTTATACCCATCGGCCCCCTG
AACAATCCCTTCCCGAGAAAGGACTCCCTTTGCCCTCAAACCATCAATGCCTTAGAGGAAGAAAGT
CGTGGATTTCTACACAAGTAGAAAGAATTTTTCAGGAGTCTGTGTATAAGACACAACAACTATATGG
AAAGAGACAAGAAATATGGTACAAAGACAATTTATTGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

Plasmid Map:


ACCN: NM_001164468

Insert Size: 1215 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164468.3](#)

RefSeq Size: 3655 bp

RefSeq ORF: 1215 bp

Locus ID: 100302736

UniProt ID: [Q86XR7](#)

Cytogenetics: 5q22.3

MW: 46.2 kDa

Gene Summary: This locus represents naturally occurring read-through transcription between the neighboring transmembrane emp24 protein transport domain containing 7 (TMED7) and toll-like receptor adaptor molecule 2 (TICAM2) genes. Alternative splicing results in multiple transcript variants, one of which encodes a fusion protein that shares sequence identity with the products of each individual gene. This fusion product functions to negatively regulate the adaptor MyD88-independent toll-like receptor 4 pathway. [provided by RefSeq, Nov 2010]
Transcript Variant: This variant (1) represents the shorter transcript and encodes the longer isoform (1), also known as TAG. This isoform shares sequence identity with the products of both the TMED7 and TICAM2 genes. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.