

Product datasheet for **SC326972**

PIP4K2C (NM_001146258) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIP4K2C (NM_001146258) Human Untagged Clone
Tag:	Tag Free
Symbol:	PIP4K2C
Synonyms:	PIP5K2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326972 representing NM_001146258. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGTCCTCCTCGGTCCCACCAGCCACGGTATCGGCGGCGACAGCAGGCCCGGCCAGGTTTCGGC
TTCGCCTCCAAGACCAAGAAGAAGCATTTCGTGCAGCAGAAGGTGAAGGTGTTCCGGGCGGCCGACCCG
CTGGTGGGTGTGTTCTGTGGGGCTAGCCCACTCGATCAATGAGCTCAGCCAGGTGCCTCCCCGGTG
ATGCTGCTGCCAGATGACTTTAAGGCCAGCTCCAAGATCAAGGTCAACAATCACCTTTTCCACAGGGAA
AATCTGCCAGTCATTTCAAGTTCAAGGAGTATTGTCCCAGGTCTTCAGGAACCTCCGTGATCGATTT
GGCATTGATGACCAAGATTACTTGGTGTCCCTTACCCGAAACCCCCAGCGAAAGTGAAGGCAGTGAT
GGTCGCTTCTTATCTCCTACGATCGGACTCTGGTCATCAAGAAGTATCCAGTGAGGACATTGCTGAC
ATGCATAGCAACCTCTCCAATATCACCAGTACATTGTGAAGTGCCATGGCAACACGCTTCTGCCCCAG
TTCTGGGGATGTACCGAGTCAGTGTGGACAACGAAGACAGCTACATGCTTGTGATGCGCAATATGTTT
AGCCACCGTCTTCTGTGCACAGGAAGTATGACCTCAAGGGTTCCTAGTGTCCCGGGAAGCCAGCGAT
AAGGAAAAGGTTAAAGAATTGCCACCCTTAAGGATATGGACTTTCTCAACAAGAACCAGAAAGTATAT
ATTGGTGAAGAGGAGAAGAAAATATTTCTGGAGAAGCTGAAGAGAGATGTGGAGTTTCTAGTGCAGCTG
AAGATCATGGACTACAGCCTTCTGCTAGGCATCCACGACATCATTCGGGGCTCTGAACCAGAGGAGGAA
GCGCCCGTGGGAGGATGAGTCAGAGGTGGATGGGACTGCAGCCTGACTGGACCTCCTGCTCTGGTG
GGCTCCTATGGCACCTCCCAGAGGGTATCGGAGGCTACATCCATTCCTCCATCGGCCCTGGGCCAGGA
GAGTTTGAAGTCTTCAATTGATGTCTATGCCATCCGGAGTGTGAAGGAGCCCCCAGAAGGAGGTCTAC
TTCATGGGCTCATTGATATCCTTACACAGTATGATGCCAAGAAGAAAGCAGCTCATGCAGCCAAAAT
GTCAAGCATGGGGCTGGGCAGAGATCTCTACTGTCCATCCGGAGCAGTATGCTAAGCGATTCTGGAT
TTTATTACCAACATCTTTGCC TAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001146258

Insert Size: 1266 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146258.1](#)

RefSeq Size: 2940 bp

RefSeq ORF: 1266 bp

Locus ID: 79837

UniProt ID: [Q8TBX8](#)

Cytogenetics: 12q13.3

Protein Families: Druggable Genome

Protein Pathways: Inositol phosphate metabolism, Phosphatidylinositol signaling system, Regulation of actin cytoskeleton

MW: 47.3 kDa

Gene Summary: May play an important role in the production of Phosphatidylinositol bisphosphate (PIP₂), in the endoplasmic reticulum.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 3' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (a).