

Product datasheet for **SC326973**

ODR4 (C1orf27) (NM_001164246) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ODR4 (C1orf27) (NM_001164246) Human Untagged Clone
Tag:	Tag Free
Symbol:	ODR4
Synonyms:	C1orf27; odr-4; TTG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326973 representing NM_001164246. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAAGAACCTACATTGTAGAAGAGACTGTTGGCCAGTATCTTTCAAACATAAATCTCCAAGGAAAG
GCTTTTGTCTCTGGCCTTTAATAGGACAGTGTTCGTCACAAAAGGATTATGTGATTCTTGCCACTAGA
ACGCCACCCAAAGAGGAGCAAAGTGAGAACCTCAAACATCCCAAAGCTAAGTTGGATAAATTGGATGAA
GAATGGGCCACAGAACATGCCTGCCAGCTAATGTTTGTGTGGAAAAGTCTATAAATAGAAAGAGATTG
TGAATTTACAGAGGAGGAAGTCTCAGAACGAGTGACACTTCACATTTGTGCTTCTACAAAAAATA
TTTTGTGCAACTATGATATCCATGATCCAAAGAGTTCAGCAAGACCAGCAGATTGGAAGTATCAAAGT
GGATTATCATCTCATGGCTTTCTTTAGAGTGTACAGTTCACATTAATATTCACATCCCACCTTTCTGCT
ACTTCTGTGAGCTATACTCTGGAGAAAAATACAAAGAATGGACTTACACGCTGGGCCAAGGAAATAGAA
AATGGTGTATTTGATTAATGGACAAGTTAAAGATGAAGATTGTGACCTATTAGAAGGACAGAAAAAA
TCTTCTAGAGGAAATACTCAAGCAACTAGTCATTTCTTTGATGTGAGAGTGCTAACGCAGTTGCTCCTG
AATTCAGACCACAGATCCACAGCCACAGTCCAGATATGTAGCGGTTCTGTAAACCTTAAGGGTGTGTG
AAATGCAGAGCTTATATCCACAGCAGTAAACCCAAAGTAAAGATGCTGTGCAGGCAGTAAAGAGGGAT
ATATTGAACACAGTTGCTGATCGTTGTGAAATGCTATTTGAGGATCTGCTTTTGAATGAAATCCAGAA
AAAAAGATTCTGAAAAAGAGTTCACAGTCCCTCCCTTATCGAGTCTTTGTTCCCTTCTTGATCCACT
GTAATGTTGTGATTATAAATTTGACGATGAGTCAGCTGAAGAAATCAGGGACCATTTTATGGAGATG
TTGGATCACACAATTCAAATAGAAGATTTGAAAATTCAGAGGAAACAAACACAGCTTGTATGAGTTCT
TCTATGAATAGTCAAGCTTCATTGGACAACACAGATGATGAACAACAAAAACAACCAATAAAACTACA
ATGTTATTGAAAATTCAGCAAAACATAGGTGTGATTGCAGCATTTACAGTTGCAGTCCTTGTGCGGGT
ATCTCCTTTCATTACTTCAGTGATTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001164246

Insert Size: 1269 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164246.1](#)

RefSeq Size: 3830 bp

RefSeq ORF: 1269 bp

Locus ID: 54953

UniProt ID: [Q5SWX8](#)

Cytogenetics: 1q31.1

Protein Families: Transmembrane

MW: 47.6 kDa

Gene Summary: May play a role in the trafficking of a subset of G-protein coupled receptors.

[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 5' coding region, compared to variant 1. The resulting isoform (3) lacks an internal segment, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.