

Product datasheet for **SC326987**

GRAMD2B (NM_001146321) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GRAMD2B (NM_001146321) Human Untagged Clone
Tag:	Tag Free
Symbol:	GRAMD2B
Synonyms:	GRAMD3; NS3TP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326987 representing NM_001146321. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGAAGACACAAAGCCTGCGAAAGTCTCGGGAAGAGGGAGAGCAAACCTGGCTCAGCCACTCAG
AGGCTGAGAAATGGTGTGGAGGAGAAAAAGAAAGCCTGCAGGTCGCCAACAGCCCAATCCCCTACCCCAT
CTGTGGAGGCGGACTCCCAGACCAGAAGAAAATCATTAGCCTATGCCCATATCTTCTCTCCAAATCTA
TGCTTCATTGAGGACCAAGAGTCAAATCCAGTTTGTGGTGCCTCTTAGCAAGTGATAAGAACGAC
TGTAAGCAAGCAAAATGACCCTAAGACTGAAAGAAAAAGTCTTCATCTCCAGCCAGTACAAG
GCCAATATGCACCTTCCAAAGTTGTTTCTAGTGTCCCAACGGAGGAACCACTGAAGCAAAGCTTTACC
TGTGCTCTACAGAAAGAAATACTATACCAAGGAAAGCTCTTGTATCAGAAAACCTGGATTGTTTTTCAT
TCCAAAGTCTTTGGAAAAGACACAAAGATCTCTATTCCAGCTTTCTCGGTAACCCCTAATAAGAAAACC
AAAAGTCTCTTAGTGCCAAACGCCCTGATCATAGCAACAGTCACAGACAGGTACATATTTGTCTCC
TACTCTCCAGAGATTCAACTTACAAACTACTAAAATCTGTGTGTGGACACTTAGAAAAACAAGTGT
GGTAACAGTCCCAATCCATCTTCTGTGAAAACAGTTTCCGAGCAGACCGCCCTTCATCTCTGCCTCTG
GATTTCAATGATGAATTCTCAGATCTGGATGGAGTGGTTCAACAAAGAAGGAAGACATGGAAGGATAT
AGCAGTTCTGGTCTCAAACCTGAATCTGAGAACTCTGAGATTTCCATGCGACAGAATCCCAAAACA
GTTCTGAATGTCTCCAAGGGAGAAGCAAAGCAACTCGGGCAGATGCCCATGTGAACAGAGTACCTGAA
GGAAAAGCCAAGAGTCTCCCTGTACAGGGTCTGTACAGAACTGTTGGAATCTTACATAAAGTCAAGTCT
CAGAAAATGTCGATGCTTACCATATTCTATATTCTATGCAATTGTTGTCTGTGCACTAATCATCTCG
ACCTTCTACATGAGATACAGAATTAATACTCTGGAGGAGCAGCTGGGGTACTAACCTCCATTGTGGAC
ACCCATAATACTGAACAGGCAGCACCATCTGGCCTGAGGTACAAGTACAATTAATGTGGAGGTTCTC
TGTCAAGAGCTTACAGCTAACATAGTGAATAGAAAAGATACAAAATAACTTACAGAAGTTGCTTGG
AATGGTGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001146321

Insert Size: 1323 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146321.1](#)

RefSeq Size: 2977 bp

RefSeq ORF: 1323 bp

Locus ID: 65983

UniProt ID: [Q96HH9](#)

Cytogenetics: 5q23.2

Protein Families: Transmembrane

MW: 49.6 kDa