

Product datasheet for **SC326988**

AGXT2L1 (ETNPPL) (NM_001146627) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGXT2L1 (ETNPPL) (NM_001146627) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETNPPL
Synonyms:	AGXT2L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326988 representing NM_001146627. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGACACTGTCACCCAGGAGTGGTCAAAGCTGCCCTGAAACAGATGGAAGTCTAAATACAAATTC
CGATTCTCCACGACAACATTGTTGAGTATGCCAAACGCCTTCAGCAACTCTGCCGGAGAACTCTCT
GTTTGTATTTTACAAATTCAGGATCCGAAGCCAACGACTTAGCCTTACGCCTGGCTCGGCAGTTTACA
GGCCACCAGGATGTGATCACTCTTGACCATGCTTACCATGGTACCTATCATCCTAATTGAGATTAGC
CCATATAAGTTTCAGAAAGGAAAAGATGTCAAAAAAGAATTTGTACATGTGGCACCAACTCCAGATACT
TACAGAGGAAAAATAGAGAAGACCATGCAGACTCAGCCAGTGTATGCAGATGAAGTGAAGAAAAATC
ATTGAAGATGCTCATAACAGTGGAAAGGAAGATTGCTGCCTTTATTGCTGAATCCATGCAGAGTTGTGGC
GGACAAATAATTCTCCAGCAGGCTACTTCCAGAAAGTGGCAGAATATGTACACGGTGCAGGGGGTGTG
TTTATAGCTGATGAAGTTCAAGTGGGCTTTGGCAGAGTTGGGAAACATTTCTGGAGCTTCCAGATGTAT
GGTGAAGACTTTGTTCCAGACATCGTCACAATGGGAAAACCGATGGGCAACGGCCACCCGGTGGCATGT
GTGGTAACAACCAAGAAATTCAGAAAGCCTTCAGCAGCTCTGGGATGGAATATTTAATACGTATGGA
GGAAATCCAGTATCTTGTGCTGTTGGTTGGCTGTCTGGATATAATTGAAAATGAAGACCTTCAAGGA
AATGCCAAGAGAGTAGGGAATATCTCACTGAGTACTGAAAAACAGAAGGCTAAACACACTTTGATA
GGAGATATTAGGGCATTGGCCTTTTTATTGGAATTTAGTGAAGGACCATGAAAAGGCCCTT
GCCACAGCTGAAGCTCAGCACATCATCTACAAGATGAAAGAAAAACGAGTGTCTCAGTGCCGATGGA
CCTCATAGAAATGTACTTAAAATAAAACCCTATGTGCTTCACTGAAGAAGATGCAAAGTTCATGGTG
GACCAACTTGATAGGATTTCTAACAGTTTTAGAAGAAGCTATGGGAACAAAACGAAAGTGTGACCTCT
GAGAATACTCCATGCAAAACAAAGATGCTGAAAGAAGCCACATAGAAGTCTTAGGGACAGCACCCT
GACTCAAAGAAAATCCCAGCAGAAAGAGAAATGGAATGTGCACGGATACACATTCAGTCTCAGTAAG
AGGCTCAAGACATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



View online »

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001146627

Insert Size: 1326 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146627.1](#)

RefSeq Size: 2036 bp

RefSeq ORF: 1326 bp

Locus ID: 64850

UniProt ID: [Q8TBG4](#)

Cytogenetics: 4q25

Protein Families: Druggable Genome

MW: 48.9 kDa

Gene Summary: Catalyzes the pyridoxal-phosphate-dependent breakdown of phosphoethanolamine, converting it to ammonia, inorganic phosphate and acetaldehyde.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (3) is shorter than isoform 1.