

Product datasheet for **SC326995**

AADACL1 (NCEH1) (NM_001146276) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	AADACL1 (NCEH1) (NM_001146276) Human Untagged Clone
Tag:	Tag Free
Symbol:	NCEH1
Synonyms:	AADACL1; NCEH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC326995 representing NM_001146276. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCAGCTGCCGCGGGCAGAAAGTTGCCGGAGGTCTCCGGTGGTATCGCCCTTCTCTTTGCCAG
CCCCTGGCGAGCCGAGCCAGGGCAAGATGAGGTCTGCTGCTCACCAGCCCTGGTGGCGCTG
GCCGCTATTACGTCTACATCCCGCTGCCTGGCTCCGTGTCGACCCCTGGAAGCTGATGCTGCTGGAC
GCCACTTCCGGGTGCACAGCAAGTGAAGTAACTGATCCACTACCTGGGACTGAGCCATCACCTGCTG
GCACTGAATTTTATCATTGTTTCTTTGGCAAAAAAGCGCGTGGTCTTCTGCCAAGTGAAGGTGACC
GACACAGACTTTGATGGTGTGGAAGTCAAGTGTGGAAGGCCCTCCGAAGCCCGAAGAGCCACTGAAA
CGCAGCGTCGTTTATATCCACGGAGGAGGCTGGGCCTTGCAAGTGAAGTGCCTGCTGACCTTCA
GATGAAATCAGGTATTATGATGAGCTGTGTACAGCAATGGCTGAGGAATTGAATGCTGTCATTGTTTCC
ATTGAATACAGGCTAGTCCAAAGGTTTATTTTCTGAGCAAAATTCATGATGTTGTACGGGCCACAAAG
TATTTCTGAAGCCAGAAGTCTTACAGAAGTATATGGTTGATCCAGGCAGAAATTTGCATTCTGGTGAC
AGTGTGGTGGAAATCTGGCTGCTGCCCTTGGACAACAGTTTACTCAAGATGCCAGCCTAAAAATAAG
CTCAAACATAAGCTTTAATTTATCCAGTCTTCAAGCTTTAGATTTTAAACACACCATCTTATCAGCAA
AATGTGAACACCCCAATCCTGCCCGCTATGTCATGGTGAAGTATTGGGTGGACTACTTCAAAGGCAAC
TATGACTTTGTGCAGGCAATGATCGTTAAACAATCACACTTCACTTGAATGTTGGAAGAGGCTGCTGCTGTC
AGGGCCCGTCTAAACTGGACATCCCTCTTGCCTGCATCCTTACAAAGAAGTACAAGCCTGTTGTACAG
ACCACAGGCAATGCCAGGATTGTCCAGGAGCTTCTCAGTTGCTGGATGCCCGCTCCGCCCACTCATT
GCAGACCAGGCAGTGCTGCAGCTCCTCCAAAGACCTACATTCTGACGTGTGAGCATGATGCTCCTCAGA
GACGATGGCATCATGTATGCCAAGCGTTTGGAGAGTGCCGGTGTGGAGGTGACCTGGATCACTTTGAG
GATGGCTTTCACGGATGTATGATTTTCACTAGCTGGCCCACTTCTCAGTGGGAATCCGGACTAGG
AATAGTTACATCAAGTGGCTAGATCAAAACCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001146276

Insert Size: 1347 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146276.1](#)

RefSeq Size: 4318 bp

RefSeq ORF: 1347 bp

Locus ID: 57552

UniProt ID: [Q6PIU2](#)

Cytogenetics: 3q26.31

Protein Families: Transmembrane

MW: 49.9 kDa

Gene Summary: Hydrolyzes 2-acetyl monoalkylglycerol ether, the penultimate precursor of the pathway for de novo synthesis of platelet-activating factor. May be responsible for cholesterol ester hydrolysis in macrophages, thereby contributing to the development of atherosclerosis. Also involved in organ detoxification by hydrolyzing exogenous organophosphorus compounds. May contribute to cancer pathogenesis by promoting tumor cell migration.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. CCDS Note: The coding region has been updated to shorten the N-terminus to one that is more supported by available conservation data.