

Product datasheet for SC327004

ER81 (ETV1) (NM_001163147) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ER81 (ETV1) (NM_001163147) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETV1
Synonyms:	ER81
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327004 representing NM_001163147. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGGATTTTATGACCAGCAAGTGCCTTACATGGTCACCAATAGTCAGCGTGGGAGAAATTGTAAAC
GAGAAACCAACAAATGTCAGGAAAAGAAAATTCATTAACAGAGATCTGGCTCATGATTGAGAAGAACTC
TTTCAAGATCTAAGTCAATTACAGGAAACATGGCTTGCAAGCTCAGGTACCTGACAATGATGAGCAG
TTTGTACCAGACTATCAGGCTGAAAGTTTGGCTTTTCATGGCCTGCCACTGAAAATCAAGAAAGAACCC
CACAGTCCATGTTACAGAAATCAGCTCTGCCTGCAGTCAAGAACAGCCCTTAAATTCAGCTATGGAGAA
AAGTGCCTGTACAATGTCAGTGCCTATGATCAGAAAGCCACAAGTGGGAATGAGGCCCTCCAACCCCC
ACACCATCCAGCAGCCAGTGTCCCACTGCATCATGCATCTCCAACTCAACTCATACCCGAAACCT
GACCGGGCCTTCCAGCTCACCTCCCTCCATCGCAGTCCATACCAGATAGCAGCTACCCATGGACCAC
AGATTTGCGCCGAGCTTTCTGAACCTGTAACCTTTCTCCTTTGCGGACGATGCCAAGGGAAGGA
CGTCCTATGTACCAACGCCAGATGTCTGAGCCAAACATCCCTTCCCACCACAAGGCTTTAAGCAGGAG
TACCACGACCCAGTGTATGAACACAACACCATGGTTGGCAGTGGCGCCAGCCAAAGCTTTCCCTCCT
CTGATGATTAACAGGAACCCAGAGATTTGCATATGACTCAGGCTGTATGTTTAAAAGGGCCCAAG
CAGTTTTATGATGACACCTGTGTGTCCAGAAAAATTCGATGGAGACATCAAACAAGAGCCAGGAATG
TATCGGGAAGGACCCACATACCAACGGCGAGGATCACTTCAGCTCTGGCAGTTTTGGTAGCTCTTCTG
GATGACCTTCAAATTCATTTTATTGCCTGGACTGGTGGAGGATGGAATTTAACTGATTGAGCCT
GAAGAGGTGGCCGACGTTGGGGCATTAGAAAAACAGGCCAGCTATGAACTATGATAAACTTAGCCGT
TCACTCCGCTATTACTATGAGAAAGGAATTATGCAAAAGGTGGCTGGAGAGAGATATGTCTACAAGTTT
GTGTGTGATCCAGAAGCCCTTTCTCCATGGCCTTTCCAGATAATCAGCGTCCACTGCTGAAGACAGAC
ATGGAACGTCACATCAACGAGGAGGACACAGTGCCTTTTCTCACTTTGATGAGAGCATGGCCTACATG
CCGGAAGGGGGTCTGCAACCCCAACCCCTACAACGAAGGCTACGTGTATTAACGCGT
ACGCGTACGCGGGCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001163147

Insert Size: 1365 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163147.1](#)

RefSeq Size: 6558 bp

RefSeq ORF: 1365 bp

Locus ID: 2115

UniProt ID: [P50549](#)

Cytogenetics:	7p21.2
Protein Families:	ES Cell Differentiation/IPS, Transcription Factors
MW:	52.5 kDa
Gene Summary:	<p>This gene encodes a member of the ETS (E twenty-six) family of transcription factors. The ETS proteins regulate many target genes that modulate biological processes like cell growth, angiogenesis, migration, proliferation and differentiation. All ETS proteins contain an ETS DNA-binding domain that binds to DNA sequences containing the consensus 5'-CGGA[AT]-3'. The protein encoded by this gene contains a conserved short acidic transactivation domain (TAD) in the N-terminal region, in addition to the ETS DNA-binding domain in the C-terminal region. This gene is involved in chromosomal translocations, which result in multiple fusion proteins including EWS-ETV1 in Ewing sarcoma and at least 10 ETV1 partners (see PMID: 19657377, Table 1) in prostate cancer. In addition to chromosomal rearrangement, this gene is overexpressed in prostate cancer, melanoma and gastrointestinal stromal tumor. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2016]</p> <p>Transcript Variant: This variant (2) lacks an alternate exon in the 5' UTR and an alternate in-frame exon in the 3' coding region, compared to variant 1. The resulting isoform (b) lacks an internal segment, compared to isoform a. This isoform contains the N-terminal TAD.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>