

Product datasheet for SC327015

FAM90A19 (NM_001164449) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM90A19 (NM_001164449) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM90A19
Synonyms:	FAM90A18
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001164449, the custom clone sequence may differ by one or more nucleotides <pre> ATGATGGCACGTGCGGACCCCAAATCTTGGGCCAAGAGACTGGTGAGAGCCCAGACCTC CAGAAGCAGCGGAGGCCCCAGTTGGGCCAAGGGCTCCCCGCCCGATGAAGAAGATCCC AGGCTCAAGTGCAAAACTGCGGGCCCTTTGGCCACACGCCAGAAAGTACCAGGTGCCCC ATGAAGTGCTGGAAGGCAGCCCTGGTTCAGCGACCTTGGGAAAAAGGAAGGAAGGAA AACCTGAAACCATGGAAGCCCCGGTGAAGCCAACCCGGGGCCCTTGAACAAGGATAAG GGAGAGAAGGAAGAGAGACCAAGGCAACAAGACCCGCAGAGGAAGGCTCTCCTCCACATG TTTTCCGGGAAACCTCCAGAGAAGCCGCTGCCGAATGGAAGGATCCACGGAGTCTTCT GATCATCTGAGGGTTGCAAGCGGGCCAATGCCGGTCCACACAACCAGTAAGAGGCCGCGC GTGGACCCTGTCTCGCTGATCGCTCAGCTGCCGAAATGTCTGGCAGGGGCTCCGTCTTG GCTTCACTGTCTCCCTCAGAAAAGCCAGCCTGAGCTCCTCCTCAAGTCTTGACCAAAG GAAAGACAGACAGGGGCTGCGGCCGACATGCCTCAGCCTGCAGTCAGGCACCAGGGCCGC GAGCCTCTCCTCGTGGTGAAGCCGACACACAGCAGCCCCGAGGGTGGTGCCGAGAAGTT CCCCAGGCTGCCTCCAAAACCCAGGCCTGCTCCAGGCCGCCAGACCCAGGCACAAGAC AAACGTCCTGCGGTGACCTCGCAGCCCTGCCGCCAGCCGCCACACAGCTTGGGCCTA GGCTCCAATCTCAGCTTCGGGCCAGGAGCAAGAGACCTGCCAGGCTCCGATTGAGGCT TGCCTGAAATTCCCCAAGAAACCGAGACTGGGTCCCTTCCAGATCCCCGAAAGCGCCATC CAGGGAGGTGAGCTGGGGGCCCCGGGAATCTCAACCTCCGCCAGCCGCAACCGAACTT GGACCAAGTACGTGCGCCCGAGATGGGCAGGAGGACACCGGCCAGGTGCCAGCGTCGAC TGGCAGCCTCCGCACAGCACACCTTGCCTGCCTACTGCCAGGCCTGCACCATGTCCCAT CACTCAGCGGCCAGCCATGATGGGGCCAGCCTCTCAGAGTGCTCTTCCGGAGACTGGAA AACGGACGCTGGAGCTCCAGCCTCCTGGCGGCCCTCATTTCACTCTCCTGAGAAGCCG GGAGCCTTCTCGCTCAGAGCCCTCATGTGTCAGAGAAGTCTGAGGCTCCCTGTGTTCTG GTCCACCGAGCGTCCTCTATGAGGACCTTCAGGTTTCTCCTCCTCAGAGGACAGCGAT TCTGACCTGGAG </pre>


[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001164449
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164449.1, NP_001157921.1</u>
RefSeq Size:	1395 bp
RefSeq ORF:	1395 bp
Locus ID:	728753
Cytogenetics:	8p23.1
Gene Summary:	FAM90A19 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, Oct 2009]