

## Product datasheet for **SC327030**

### AMHR2 (NM\_001164691) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	AMHR2 (NM_001164691) Human Untagged Clone
Tag:	Tag Free
Symbol:	AMHR2
Synonyms:	AMHR; MISR2; MISRII; MRII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC327030 representing NM\_001164691.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTAGGGTCTTTGGGGCTTTGGGCATTACTTCCCACAGCTGTGGAAGCACCCCAACAGGCCAACC
TGTGTGTTCTTTGAGCCCTGGAGTGCGGGAAGCACAAAGACACTGGGAGAGCTGCTAGATACAGGC
ACAGAGCTCCCCAGAGCTATCCGCTGCCTCTACAGCCGCTGCTGCTTTGGGATCTGGAACCTGACCCAA
GACCGGGCACAGGTGAAATGCAAGGATGCCGAGACAGTGATGAGCCAGGCTGTGAGTCCCTCCACTGT
GACCCAAGTCCCCGAGCCACCCAGCCCTGGCTCCACTCTTTACCTGCTCCTGTGGCACTGACTTC
TGCAATGCCAATTACAGCCATCTGCCTCTCCAGGGAGCCCTGGGACTCCTGGCTCCCAGGGTCCCAG
GCTGCCCCAGGTGAGTCCATCTGGATGGCACTGGTGTCTGCTGGGGCTGTTCTCTCTCTCTGCTGTG
CTGGGCAGCATCATCTTGCCCTGCTACAGCGAAAGAATAACAGAGTGCAGAGGTGAGCCAGTGCCAGAG
CCAAGGCCAGACTCAGGCAGGGACTGGAGTGTGGAGCTGCAGGAGCTGCCTGAGCTGTGTTCTCCCAG
GTAATCCGGGAAGGAGGTGATGCAGTGGTTTGGGCCGGGCAGCTGCAAGGAAAACGTTGCCATCAAG
GCCTTCCACCGAGGTCTGTGGCTCAGTCCAAGCTGAGAGAGCATTGTACGAACCTCCAGGCCTACAG
CACGACCATTGTCCGATTTATCACTGCCAGCCGGGGGGTCTGGCCGCTGCTCTGGGCCCTG
CTGGTACTGGAATGCATCCCAAGGGCTCCCTGTGCCACTACTTGACCCAGTACACCAGTACTGGGGA
AGTTCCCTGCGGATGGCACTGTCCCTGGCCAGGGCTGGCATTCTCCATGAGGAGCGCTGGCAGAAT
GGCCAATATAAACAGGTATTGCCACCCAGAGATCTGAGCAGCCAGAATGTGCTCATTGCGGAAGATGGA
TCGTGTGCCATTGGAGACCTGGGCCTTGCTTGGTGTCCCTGGCCTCACTCAGCCCCCTGCCTGGACC
CCTACTCAACCACAAGGCCAGCTGCCATCATGGAAGACCCTGATGGGCTGAGGGAGCTCCTAGAAGAC
TGTGGGATGCAGACCCAGAAGCACGGCTGACAGCTGAGTGTGTACAGCAGCGCCTGGCTGCCTTGCC
CATCTCAAGAGAGCCACCCCTTTCCAGAGAGCTGTCCAGTGGCTGCCACCTCTCTGCCAGAAGAC
TGTACTTCAATTCTGCCCTACCATCCTCCCCTGTAGGCCTCAGCGGAGTGCCTGCCACTTCAGCGTT
CAGCAAGGCCCTTGTCCAGGAATCCTCAGCCTGCCTGTACCCTTCTCTGTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001164691

**Insert Size:** 1437 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001164691.1</a></u>
<b>RefSeq Size:</b>	1578 bp
<b>RefSeq ORF:</b>	1437 bp
<b>Locus ID:</b>	269
<b>UniProt ID:</b>	<u><a href="#">Q16671</a></u>
<b>Cytogenetics:</b>	12q13.13
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>Protein Pathways:</b>	Cytokine-cytokine receptor interaction, TGF-beta signaling pathway
<b>MW:</b>	51.8 kDa
<b>Gene Summary:</b>	<p>This gene encodes the receptor for the anti-Mullerian hormone (AMH) which, in addition to testosterone, results in male sex differentiation. AMH and testosterone are produced in the testes by different cells and have different effects. Testosterone promotes the development of male genitalia while the binding of AMH to the encoded receptor prevents the development of the mullerian ducts into uterus and Fallopian tubes. Mutations in this gene are associated with persistent Mullerian duct syndrome type II. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Sep 2009]</p> <p>Transcript Variant: This variant (3) lacks two exons in the 3' coding region compared to variant 1. The resulting protein (isoform 3) is shorter compared to isoform 1. Variants 2 and 3 encode isoforms that are the same length, but have distinct protein sequences.</p>