

Product datasheet for **SC327031**

GOLGA6L10 (NM_001164465) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA6L10 (NM_001164465) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA6L10
Synonyms:	GOLGA6L18
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001164465, the custom clone sequence may differ by one or more nucleotides

```

ATGTGGCCCCAACCCCGCTCCCTCCCCACCCGCGATGTCAGAAAAACACAACAGGGG
AAATTGGCCGCAGCCAAGAAAAAGTTAAAAGCATATTGGCAGAGGAAGAGCCCTGGCATT
CCAGCAGGAGCTAACAGGAAAAAGAAAGTCAATGGCAGTAGCCCTGACACAGCCACTTCT
GGTGGTTACCACTCACCTGGGGATTGAGCAACAGGTATCTACGGGGAGGGCCGTGCATCC
TCTACTACCTGCAGGATCTGGAGAGCCAGTACCAAGAACTAGCAGTGGCCCTGGATTCA
AGCTCCGCAATAATCAGTCAACTCACTGAAAACATCAATTCACTGGTTCGCACATCTAAG
GAGGAGAAGAAGCATGAGATACATCTGGTACAGAAGCTTGGGAGGAGCTTGTTCAAACCTC
AAAAACCAGACGGCTGAACCCCTGGCCCCAGAGCCCCAGCAGGGCCATCTAAGGTGGAG
CAGCTACAAGATGAGACCAACCACCTAAGGAAGGAGCTAGAGAGTGTGGGAAGACAGCTC
CAGGCTGAGGTGGAAAACAATCAGATGTTGAGTCTCCTGAACAGGAGACAGGAGGAGAGG
CTACGTGAACAGGAGGAGAGGCTACATGAACAGGAGGAGAGGCTACATGAGCAGGAGGAG
AGGCTGTGTGAACAGGAGGAGAGGCTGTGTGAACAGGAGGAGAGGCTACGTGAACAGGAG
GAGAGGCTACGTGAACAGGAGGAGAGGCTACATGAACAGGAGGAGAGGCTGTGTGAACAG
GAGGAGAGGCTATGTGAACAGGAGGAGAGGCTGTGTGAACAGGAGAAGCTGCCAGGGCAG
GAGAGGCTGCTGGAAGAGGTGGAGAAGCTGTTAGAACAGGAGAGGCGGCAGGAGGAGCAG
GAGAGGCTGCTGGAGAGGAGAGGCTGCTGGACGAGGTGGAGGAGCTCCTGGAGCAGGAG
AGGCTTCGGCAACAGGATGAGAGGCTGTGGCAGCAGGAGACTCTCGGGAGCTGGAGAGG
CTGCGGGAGCTGGAGAGGCTGCGGGAGCTGGAGAGGCTGCGGGAGCTGGAGAGGATGCTG
GAGCTGGGGTGGGAAGCCCTGTACGAGCAGCGGGCCGAGCCACGAGCGGCTTCGAGGAG
CTGAACAACGAGAACAAGAGCACACTGCAGTTGGAGCAGCAAGTAAAGGAGCTGGAGAAG
TCGGGTGGAGCTGAAGAGCCAAGAGGCTCCGAGTCTGCAGCAGCAGCCAGACCAGTACCT
GGAGCCCCAGTCCACAAGGAGCTTGATGTGCGGACAAGCAGGGTGGACTCCCCAGGAG
CACCCAGGCTTGAGTGGAGAAGCTGTTGGTACAGGAGAGGCGGCAGGAGGAGCAGAAGAG
GCTGCATGCCATTCTTTTCGGGCTGCGGAGAACAGGGAGCTAAACATCACCATCATC

```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001164465
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164465.1, NP_001157937.1</u>
RefSeq Size:	3548 bp
RefSeq ORF:	1440 bp
Locus ID:	647042
UniProt ID:	<u>A6NI86</u>
Cytogenetics:	15q25.2