

Product datasheet for SC327039

TCF7L2 (NM_001146283) Human Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: TCF7L2 (NM_001146283) Human Untagged Clone

Tag: Tag Free Symbol: TCF7L2

Synonyms: TCF-4; TCF4

Mammalian Cell Neomycin

Selection:

Vector:

pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn





Fully Sequenced ORF:

>SC327039 representing NM_001146283.

Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC


ATGCCGCAGCTGAACGGCGGTGGAGGGGATGACCTAGGCGCCAACGACGAACTGATTTCCTTCAAAGAC GAGGGCGAACAGGAGGAGAAGAGCTCCGAAAACTCCTCGGCAGAGAGGGGATTTAGCTGATGTCAAATCG TCTCTAGTCAATGAATCAGAACGAATCAAAACAGCTCCTCCGATTCCGAGGCGGAAAGACGGCCTCCG CCTCGCTCCGAAAGTTTCCGAGACAAATCCCGGGAAAGTTTGGAAGAAGCGGCCAAGAGGCAAGATGGA GGGCTCTTTAAGGGGCCACCGTATCCCGGCTACCCCTTCATCATGATCCCCGACCTGACGAGCCCCTAC CTCCCCAACGGATCGCTCTCGCCCACCGCCCGAACCTATCTCCAGATGAAATGGCCACTGCTTGATGTC CAGGCAGGGAGCCTCCAGAGTAGACAAGCCCTCAAGGATGCCCGGTCCCCATCACCGGCACACATTGTC AGCCCCTCCCTTGCTGCACTCAGGGACATGACTGTCAGCACTTCTACCCCCCCTCAGACTTCACTGTC AGCACTCAAGTCTTCAGGGACATGAAAAGGAGCCACTCCTTACAAAAAGTTGGGGAGCCCTGGTGTATT GAGTCTAACAAAGTGCCAGTGGTGCAGCACCCTCACCATGTCCACCCCTCACGCCTCTTATCACGTAC AGCAATGAACACTTCACGCCGGGAAACCCACCTCCACACTTACCAGCCGACGTAGACCCCAAAACAGGA ATCCCACGGCCTCCGCACCCTCCAGATATATCCCCGTATTACCCACTATCGCCTGGCACCGTAGGACAA ATCCCCCATCCGCTAGGATGGTTAGTACCACAGCAAGGTCAACCAGTGTACCCAATCACGACAGGAGGA CCACCACATCATACGCTACACACGACGGGCATTCCGCATCCGGCCATAGTCACACCAACAGTCAAACAG GAATCGTCCCAGAGTGATGTCGGCTCACTCCATAGTTCAAAGCATCAGGACTCCAAAAAAGGAAGAAGAA AAGAAGAAGCCCCACATAAAGAAACCTCTTAATGCATTCATGTTGTATATGAAGGAAATGAGAGCAAAG GTCGTAGCTGAGTGCACGTTGAAAGAAAGCGCGGCCATCAACCAGATCCTTGGGCGGAGGTGGCATGCA CTGTCCAGAGAAGAGCAAGCGAAATACTACGAGCTGGCCCGGAAGGAGCGACAGCTTCATATGCAACTG TACCCCGGCTGGTCCGCGCGGGATAACTATGGAAAGAAGAAGAAGAGGAAAAAGGGACAAGCAGCCGGGA CTTCAGATGGAAGCTTACTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:

Sgfl-Mlul

Plasmid Map:



ACCN: NM 001146283



TCF7L2 (NM_001146283) Human Untagged Clone - SC327039

Insert Size: 1470 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a

> point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative

RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning

> into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube Components:

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001146283.1

RefSeq Size: 4021 bp RefSeq ORF: 1470 bp Locus ID: 6934

UniProt ID:

Q9NQB0

Cytogenetics: 10q25.2-q25.3

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Acute myeloid leukemia, Adherens junction, Arrhythmogenic right ventricular

cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer,

Melanogenesis, Pathways in cancer, Prostate cancer, Thyroid cancer, Wnt signaling pathway

54.6 kDa MW:



Gene Summary:

This gene encodes a high mobility group (HMG) box-containing transcription factor that plays a key role in the Wnt signaling pathway. The protein has been implicated in blood glucose homeostasis. Genetic variants of this gene are associated with increased risk of type 2 diabetes. Several transcript variants encoding multiple different isoforms have been found for this gene.[provided by RefSeq, Oct 2010]

Transcript Variant: This variant (3) has multiple differences in the coding region, compared to variant 1, one of which results in a translational frameshift. The resulting protein (isoform 3) has a distinct C-terminus and is shorter than isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.