

Product datasheet for **SC327043**

AGXT2L1 (ETNPPL) (NM_001146590) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGXT2L1 (ETNPPL) (NM_001146590) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETNPPL
Synonyms:	AGXT2L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327043 representing NM_001146590.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGCGAGCTGTACAGTAAGCGGGACACTCTGGGGCTGAGGAAGAAGCACATCGGGCCCTCATGCAAA
GTTTTCTTGCATCGGATCCCATCAAAATAGTGAGAGCCAGAGGCAGTACATGTTTGATGAGAACGGT
GAACAGTACTTGGACTGCATCAACAATGTTGCCATGGAGTGGTCAAAGCTGCCCTGAAACAGATGGAA
CTGCTAAATACAAATTCTCGATTCTCCACGACAACTTGTGATGATGCCAAACGCCTTTCAGCAACT
CTGCCGGAGAAACTCTCTGTTTGTATTTTACAAATTCAGGATCCGAAGCCAACGACTTAGCCTTACGC
CTGGCTCGGCAGTTCAGAGGCCACCAGGATGTGATCACTCTTGACCATGCTTACCATGGTCACCTATCA
TCCTTAATTGAGATTAGCCCATATAAGTTTCAGAAAGGAAAAGATGTCAAAAAGAATTTGTACATGTG
GCACCACTCCAGATACTTACAGAGGAAAATATAGAGAAGACCATGCAGACTCAGCCAGTCTTATGCA
GATGAAGTGAAGAAAATCATTGAAGATGCTCATAACAGTGAAGGAAGATTGCTGCCTTTATTGCTGAA
TCCATGCAGAGTTGTGGCGGACAATAATTCCTCCAGCAGGCTACTTCCAGAAAGTGGCAGAATATGTA
CACGGTGCAGGGGTGTGTTTATAGCTGATGAAGTTCAAGTGGGCTTTGGCAGAGTTGGGAAACATTC
TGGAGCTTCCAGATGTATGGTGAAGACTTTGTTCCAGACATCGTCACAATGGGAAAACCGATGGGCAAC
GGCCACCCGGTGGCATGTGTGTTAAACAACCAAGAAAATTGCGAAGCCTTCAGCAGCTCTGGGATGGAA
TATTTTAAACGTATGGAGGAAATCCAGTATCTTGTGCTGTTGGTTTGGCTGTCTGGATATAATTGAA
AATGAAGACCTTCAAGGAAATGCCAAGAGAGTAGGGAATTATCTCACTGAGTTACTGAAAAACAGAAAG
GCTAAACACACTTTGATAGGAGATATTAGGGGCATTGGCCTTTTTATTGGAATTGATTTAGTGAAGGAC
CATCTGAAAAGGACCCCTGCCACAGCTGAAGCTCAGCACATCATCTACAAGATGAAAGAAAAACGAGTG
CTTCTCAGTGCCGATGGACCTCATAGAAATGTACTTAAAAATAAAACCACCTATGTGCTTCACTGAAGAA
GATGCAAAGTTCATGGTGGACCAACTTGATAGGATTCTAACAGTTTTAGAGAAGCTATGGGAACCAAAA
ACCGAAAGTGTGACCTCTGAGAATACTCCATGCAAAAACAAGATGCTGAAAGAAGCCACATAGAATCG
CTTAGGGACAGCACCCTGACTCCAAGAAAATCCCAGCAGAAAGAGAAATGGAATGTGCACGGATACA
CATTCACTGCTCAGTAAGAGGCTCAAGACATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001146590

Insert Size:	1482 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001146590.1
RefSeq Size:	2115 bp
RefSeq ORF:	1482 bp
Locus ID:	64850
UniProt ID:	Q8TBG4
Cytogenetics:	4q25
Protein Families:	Druggable Genome
MW:	55 kDa
Gene Summary:	<p>Catalyzes the pyridoxal-phosphate-dependent breakdown of phosphoethanolamine, converting it to ammonia, inorganic phosphate and acetaldehyde.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 5' coding region compared to variant 1. This results in a shorter protein (isoform 2) compared to isoform 1.</p>