

Product datasheet for **SC327048**

ZNF727 (NM_001159522) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF727 (NM_001159522) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF727
Synonyms:	ZNF727P
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001159522, the custom clone sequence may differ by one or more nucleotides

```

ATGCGAGTGCTAACATTCAGGGATGTGGCTGTAGAATTCTCCCCAGAAGAGTGGAATGC
CTGGACTCTGCTCAGCAGCGTTTGTATAGGGATGTGATGTTAGAGAACTACGGAAACCTG
TTCTCCTTGGGTCTTGCTATCTTTAAGCCAGACTTGATTACCTATCTGGAGCAAAGAAAA
GAGCCTTGAATGCGAGGAGACAGAAGACAGTAGCCAAACACCCAGCTGGCTCTTTGCAT
TTTACTGCAGAGATATTGCTGGAGCAGCAGATAAACGATTCATTTCAAAAAGTGATCCTG
AGAAAACTGGAAGCTGTGACCTTAATACTTTACGTTTAAAGAAAGACTACCAACGTGTG
GGTAATTGCAAGGGGCAGAAAAGCAGTTATAATGGCATTCAATGTTTGTGAGCTACC
CGTAGCAAACCTGTCAATATAATAAATGTGGCAAAGCTTTTGGGTTGTGCTCAATCTTC
ACTGAACATAAGAAAATTTTTCAGAGAGAGAAATGCTACAAATGTGAAGAATGTGGCAA
GACTGTAGGTTGTGAGATTTTACCATACAGAAGAGAAATTCATACTGCAGATAGAAGTTAC
AAATGTGAAGAATGTGGCAAAGCCTGTAAAAAGTTCTCAAACCTTACTGAACATAATAGA
GTTCACTACTGGAAAGAAACCTACAAATGTGAAGAATGTGGCAAACATTTACCTGTTCC
TCAGCCCTTACTAAACACAAGAGAAATCATACTGGAGACAGACCTACAAATGCGAAGAA
TGTCACAAAGCCTTTAGGTGTTGCTCAGACCTTACTAAACATAAGAGAATTCATACTGGA
GAGAAACCTACAAATGTAAAGAATGTCACAAAGCCTTTAGGTGTTGCTCAGACCTTACT
AAACATAAGAGAATTCATACTGGAGAGAAACCTACAAATGTAATGAATGTGAAAAGCT
TTTATGTGGATCTCGGCCCTTAGTCAACATAACAGAATTCATACTGGAGAGAAACCTAC
ATTTGTGAAGAATGTGGCAAAGCCTTTACCTACTCCTCAACCTTATTAGCCACAAGAGA
ATTCATATGGAATTGAGACCTTACAAATGTGAAGAATGTGGCAAACCTTTAAGTGGTTC
TCAGACCTGACTAATCATAAGAGAATTCACACTGGAGAGAAACCTACAAATGTGAAGAA
TGTGGCAAAGCTTTACCTGCTCCTCAAACCTTATTAACACAAGAGAATTCATATGGAA
GTGAGACCTTACAAATGTGAAGAATGTGGCAAACCTTTAAGTGGTTCAGACCTGACT
AATCATAAGAGAATTCACACTGGAGAGAAAGCCCTACAAATGTGAAGAATGTGGCAAAC
TTTACCTGCTCCTCAAGCCTTATTAACACAAGAGAAGTCATACTGGAGACAGACCTACA
AGTGCAAAGAATGTGGCAAAGCCTTTAGGTGTTCTCAGACCTTACTAAACATAAGG

```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001159522
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159522.1</u> , <u>NP_001152994.1</u>
RefSeq Size:	1679 bp
RefSeq ORF:	1500 bp
Locus ID:	442319
UniProt ID:	<u>A8MUV8</u>
Cytogenetics:	7q11.21
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]