

Product datasheet for **SC327050**

TULP3 (NM_001160408) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TULP3 (NM_001160408) Human Untagged Clone
Tag:	Tag Free
Symbol:	TULP3
Synonyms:	TUBL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327050 representing NM_001160408.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGCTTCGCGTCCCGGCTCAGTCCAGCGGGACAGTGTCTTCCATGAAGAAATGATGAAGATG
CGACAGGCTAAGCTGGATTATCAGAGGCTACTACTTGAAGAGGCAAAGGAAAAAGCGCCTTGAGCCA
TTTATGGTGCAGCCCAATCCAGAAGCCAGGCTAGTTCGGCAAAGCCAAGGGCCAGTGATGAGCAGACT
CCCTTGGTGAAGTGCATACTCCCCACAGCAATGTCATCTTACATGGTATTGATGGTCCAGCTGCTGTC
CTGAAACCAGACGAAGTTCATGCTCCATCAGTAAGCTCCTCTGTTGTGGAAGAAGATGCTGAAAACACC
GTGGATACTGCTTCCAAGCCAGGACTTCAGGAGCGTCTCCAAAAGCATGATATCTCTGAAAGTGTGAAC
TTCGATGAGGAGACTGATGGAATATCCCAGTCAGCATGTTTAGAAAGACCCAATTCTGCATCAAGCCAG
AATTCAACCGATACAGGCACTCCGGTTCTGCTACTGCCGCCAACCAGCTGATAACCTCCTGGGAGAC
ATAGACGACCTGGAGACTTTGTGTATAGTCTGCCCTCAAGGTGTACAGTAAGATGTCGGATAATC
CGGGATAAAAGGGGAATGGATCGGGGTCTTTCCCCACCTACTATATGTACTTGGAAAAAGAAGAAAT
CAGAAGATATTTCTTCTTGCAGCTAGAAAAGCGAAAAAGAGCAAAACAGCCAACCTACCTTATCTCCATT
GATCCAGTTGATTTATCTCGTGAAGGAGAAAGTTATGTCGGCAAGCTTAGATCCAACCTCATGGGGACC
AAGTTTACAGTTTATGACCGTGGCATCTGCCCCATGAAGGGCCGGGGTTTGGTAGGAGCGGCCACACC
CGGCAGGAGCTGGCTGCCATCTCCTATGAAACAACGTAAGTTGAAAGTCTAGGAAAATGTCT
GTGATCATTCTTGAATGACACTGAATCATAAGCAGATCCCCTATCAGCCACAAAACAACCATGACAGT
TTGCTCTCAAGGTGGCAGAACAGAACTATGAAAATCTGGTTGAGCTGCACAACAAGGCCCGTCTGG
AACAGTGACACTCAGTCTATGTCCTCACTCCGTGGCCGGTCACTCAGGCGTCTGTGAAGAATTC
CAGATAGTCCACAAAAATGACCCTGATTATATAGTCATGCAGTTTGGACGTGTGGCAGATGACGTGTC
ACACTGGATTACAACCTACCACCTTTGTGCAGTACAGGCCTTTGGCATCGGTCTTTCTAGCTTTGACAAA
TGTATCCAAACCTTGAAGTGCAGGAGCTCTGTGAGCTCCACCCTCAGCACCATTAGCTGCATCCCTT
GTGCACAGAACTGTTTGCAGCGTTGGGTGGGACACCCGTGGCGGCTGCTCCCTCAGACCTCCCTCTG
TGGACTGACCTCTCACCTCCGCCTGTTGTTCTGCACCACATCAGATAAGCATGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001160408

Insert Size:	1506 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001160408.1
RefSeq Size:	2001 bp
RefSeq ORF:	1506 bp
Locus ID:	7289
UniProt ID:	O75386
Cytogenetics:	12p13.33
Protein Families:	Druggable Genome, Transcription Factors
MW:	56.5 kDa
Gene Summary:	<p>This gene encodes a member of the tubby gene family of bipartite transcription factors. Members of this family have been identified in plants, vertebrates, and invertebrates, and they share a conserved N-terminal transcription activation region and a conserved C-terminal DNA and phosphatidylinositol-phosphate binding region. The encoded protein binds to phosphoinositides in the plasma membrane via its C-terminal region and probably functions as a membrane-bound transcription regulator that translocates to the nucleus in response to phosphoinositide hydrolysis, for instance, induced by G-protein-coupled-receptor signaling. It plays an important role in neuronal development and function. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, May 2009]</p> <p>Transcript Variant: This variant (2) differs in the 3' coding region and 3' UTR, compared to variant 1. The encoded isoform (2) has a longer and distinct C-terminus compared to isoform 1.</p>