

Product datasheet for **SC327055**

PODNL1 (NM_001146254) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	PODNL1
Synonyms:	SLRR5B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



Fully Sequenced ORF: >SC327055 representing NM_001146254.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGAGGCCCACCCAGTGGCCGAGCCTGCTGCTGCTCCTGCTGTTGCCGGGGCCCCGCCCGTCGCC
GGCTTGGAAGACGCTGCCTTCCCCACCTGGGGAGAGCTTGCAGCCCCGCCCGGGCCTGTCCCTG
CGCTGCTCCTGCCCGGAGTCGACACTGTGGACTGTGATGGCTTGGACCTTCGAGTGTTCGCGACAAC
ATCACCAGAGCCGCTCAGCACCTCTCCCTGCAGAACAAACAGCTCCAGGAATCCCCACAATGAGCTG
TCCCGCCTCAGTGGCCTGCGAACCTCAACCTCCACAACAACCTCATCTCTCCGAAGGCCTGCCTGAC
GAGGCCTTCGAGTCCCTCACCAGCTGCAGCACCTCTGCGTGGCTCACAACAAGCTCTCAGTGGCCCT
CAGTTTCTGCCCGGTCCTCCGTGTCGCGGATCTGGCTGCCAACCAAGTATGGAGATCTTCCCCCTC
ACCTTTGGGAGAGAAGCCGCACTCAGGTCCGTGTACCTCCACAACAACAGCTGAGCAACGCTGGCCTG
CCCCCGACGCTTCCGCGGCTCCGAGGCCATCGCCACCTCAGCCTCTCCAACAACAGCTCAGCTAC
CTGCCCGCAGCCTGCCCGCTCACTCGAGCGGCTCCACCTGCAGAACATCTCATCTCAAGGTGCC
CGAGGAGCCCTGAGCCCGCAGACTCAACTCCGTGAGCTCTACCTCCAGCACAACAGCTGACAGACAGT
GGCCTGGATGCCACCACCTTCAGCAAGCTGCATAGCCTTGAATACCTGGATCTCTCCACAACAGCTG
ACCACAGTGCCCGCGGCTGCCCGGACCTGGCTATCTGCACCTGGGCCGAACCGCATCCGGCAG
GTGGAGGCGGCTCGGCTGCACGGGGCGGTGGTCTGCGCTATTTGTTGCTGCAGCACAACAGCTGGG
AGCTCAGGGCTGCCCGCGGGGCTCTGCGGCCGTGCGGGGCTGCACAGCTGCACCTCTATGGCAAT
GGGCTGGACCGCTGCCTCCAGCCCTGCCCGCGCCTGCGTGGCTGCTGCCCAACAACAGCTG
GCCGCGCTGGGTGCCGTGACCTGGTGCACACCGGGCTGACGGAGCTTAACCTGGCTATAACCGC
CTGGCAGCGCCGTGTGCACACCGGGCTTCCGCGGTTGCGTGGCTGCGCAGCCTCGACCTGGCA
GGGAATCAGCTAACCGGCTGCCATGGGCTGCCACTGGCTGCGCACCTGCAGCTGCACGCAAC
CAGCTGCGGATGCTCGAGCCCGAGCCTCTGGCCGGCTGGACCAACTGCGGGAGCTCAGCCTGGCGCAC
AACCGGCTCCGGGTGCGGACATCGGGCCAGGCACCTGGCATGAGCTCCAAGCCCTCCAGGTGAGGCAC
AGGCTGGTTAGCCACACTGTCCCCAGGGCCCTCCATCCCCCTGCTGCCCTGCCACGTCCCAACATT
CTAGTTAGCTGGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001146254

Insert Size: 1533 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001146254.1</u>
RefSeq Size:	3308 bp
RefSeq ORF:	1533 bp
Locus ID:	79883
UniProt ID:	<u>Q6PEZ8</u>
Cytogenetics:	19p13.12
MW:	56.4 kDa