

Product datasheet for **SC327060**

SARS2 (NM_001145901) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SARS2 (NM_001145901) Human Untagged Clone
Tag:	Tag Free
Symbol:	SARS2
Synonyms:	mtSerRS; SARS; SARSM; SerRS; SerRSmt; SERS; SYS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327060 representing NM_001145901.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGCGTCCATGGCGGGCGCTTGTGGCCTTTGCTGACTCGTGGGGTTCCGGCCCCGGGGAGGC
TGCATCTCCAACGATAGTCCAAGGAGAAGTTTCACTACAGAGAAACGAAACCGGAACCTCCTGTACGAG
TATGCGCGGAGGGCTACAGCGCACTCCCTCAGCTGGACATAGAGCGTTCTGCGCATGCCAGAAGAG
GCCGCACACGCCCTGGAGCTCCGCAAGGGGAGCTGCGCTCGGGCGACCTGCCCGCATCATCTCGACA
TGGCAGGAGCTGAGGCAGCTGCAGGAGCAGATCCGGAGCCTGGAGGAAGAGAAGGCAGCTGTGACTGAG
GCAGTGGGGCCCTGCTGGCAAACAGGACAGTGGTGAAGTGCAGCAGGTGCGGCTGGATCCAGGTGCT
GGCTCCATATTTGGTCTACGTTCTCCATTCCAGGCCAGTTTCTCTCTTGTGGAGGCCAGCTT
GAGGAGCAGTTTACCTGCAGGCGCTGAAGCTGCCAACCAGACCCAGACGTGCCCGTGGGGAT
GAGAGCCAGGCTCGAGTGTCCACATGGTCGGAGACAAGCCAGTTTCTCCTTCCAACCTCGGGGCCAC
CTGGAAATTTGGCGAGAACTCGACATCATCCGTCAGAAGCGCCTGTCCCAGTGTCTGGCCACCGGTCC
TATTACCTGCGCGGGGCTGGAGCCCTCCTGCAGCACGGCCTGGTCAACTTCAATTCAACAAGCTTCTC
CGCCGGGGCTTACCCCCATGACGGTGCCAGACCTTCTCCGCGGAGCAGTGTGTTGAAGGCTGTGGGATG
ACACCAAATGCCAACCCATCCAAATTTACAACATCGACCCTGCCCGCTTCAAAGATCTCAACCTGGCT
GGAACAGCGGAGGTGGGGCTTGCAGGCTACTTCATGGACCACACCGTGGCCTTCAGGGACCTGCCAGTC
AGGATGGTTTGTCCAGCACCTGCTACCGGGCAGAGACAAACACGGGACAGGAACCCCGGGGCTGTAT
CGAGTACACCACTTACCAAGGTGGAGATGTTTGGGGTGCAGGCCCTGGGCTGGAGCAGAGCTCACAG
CTGCTGGAGGAGTTCCTGTCCCTCAGATGGAGATCTTGACAGAGCTGGGCTTGCACCTCCGGGCTCTG
GATATGCCACCAAGAAGTGGGCTCCCGCCTACCGCAAGTTTGACATTGAGGCTGGATGCCAGGC
CGAGGCCGCTTTGGAGAGGTACCAGTGTCTCCAAGTGCACAGACTTCCAGAGCCCGCCTCCACATC
ATGTTCCAGACCGAGGCTGGGGAGCTGCAGTTGCCACACGGTGAACGCCACCGCCTGTGCTGTCCCC
CGCCTTCTCATCGCGCTCCTGGAGAGTAACCAGCAGAAGGACGGCTCAGTGTGTCGCCCTGCCCTC
CAGTCTACCTCGGCACTGATCGGATCACAGCCCTACCCAGTGCCTCTCCAGTACATCGGCCCAAC
CAGCCCCGAAGCCTGGGCTGCCTGGCCAGCCTGCTGTAAGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001145901

Insert Size:	1563 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145901.1
RefSeq Size:	2083 bp
RefSeq ORF:	1563 bp
Locus ID:	54938
UniProt ID:	Q9NP81
Cytogenetics:	19q13.2
Protein Pathways:	Aminoacyl-tRNA biosynthesis
MW:	58 kDa
Gene Summary:	<p>This gene encodes the mitochondrial seryl-tRNA synthetase precursor, a member of the class II tRNA synthetase family. The mature enzyme catalyzes the ligation of Serine to tRNA(Ser) and participates in the biosynthesis of selenocysteinyl-tRNA(sec) in mitochondria. The enzyme contains an N-terminal tRNA binding domain and a core catalytic domain. It functions in a homodimeric form, which is stabilized by tRNA binding. This gene is regulated by a bidirectional promoter that also controls the expression of mitochondrial ribosomal protein S12. Both genes are within the critical interval for the autosomal dominant deafness locus DFNA4 and might be linked to this disease. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Mar 2009]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>