

Product datasheet for **SC327062**

ROD1 (PTBP3) (NM_001163788) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ROD1 (PTBP3) (NM_001163788) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTBP3
Synonyms:	ROD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327062 representing NM_001163788.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATAGTTCTACTCCTTCTACAGGTGTGTATGCTAATGGGAATGACAGCAAGAAATTTAAACGAGAT
AGACCCCTGTTCCGCTTCCCGTGTTCATCTCGAAAAATCCATGTGATGTCACCGAAGCAGAG
ATCATATCATTAGGTCTACCATTTGGCAAAGTAACATAATCTTTTATGTTGAAAGGAAAAAGCCAGGCT
TTCTTAGAAATGGCTTCTGAGGAAGCTGCCGTTACTATGGTGAATTATTACACTCCTATTACTCCTCAC
CTTCGAAGCCAGCCTGTTTATATTCAGTATTCCAATCACAGAGAAGCTTAAGACTGACAATCTACCTAAT
CAAGCTCGAGCCCAAGCTGCACTGCAGGCTGTCAGTGCCGTCCTCAATCAGGAAGCCTGGCCCTTCTGGA
GGTCTTCCAATGAAGGCACAGTCTACCTGGGCAGAGCCCTGTGCTTGAATAATTATTGAAAACCTC
TTTTACCTGTTACCCTGGAAGTCTTTCATCAGATATTTCTAAATTTGGCACAGTCTTGAAGATTATC
ACCTTTACAAAGAATAATCAGTTTCAAGCCTTGCTTCAGTATGCTGACCCAGTAAATGCACATTATGCC
AAAATGGCTCTGGATGGCCAGAATATCTATAATGCATGCTGCACTCTGCCATTGACTTCTCCAAGCTC
ACCAGCCTAATGTGAAATATAATAATGACAAAAGCAGAGACTTCACTCGCTTAGACCTTCTACTGGT
GATGGCCAGCCATCCCTTGAACCCCTATGGCTGCTGCTTTTGGTGCACCGGGTATAATTTCTTACCA
TATGCAGGGGCTGCTGGATTTGCCCCAGCCATTGGATTTCTCAAGCTACAGGTCTATCAGTTCCAGCT
GTTCTGGAGCTCTTGGTCTCTCAAAATCACCTCTTCTGCTGCTCACTGGAAGGATGGCCATTCTGGG
GCTAGTGGTATACCAGGAAATCTGTTCTACTCGTCACAAATCTCAATCTGATCTTATCACACCACAT
GGGCTTTTATCCTATTTGGAGTCTATGGTGTATACATCGAGTGAAGATTATGTTTAATAAGAAAGAA
AATGCCTTGGTTCAGATGGCGGATGCAAATCAAGCTCAGCTAGCAATGAACCATTAAGTGGTCAGAGA
CTTTATGGGAAAGTCTTCGTGCTACACTGTCCAAACATCAAGCAGTACAGCTTCCCTCGAGAGGGACAA
GAAGACCAAGGTCTGACTAAGGATTTTCAAGCAATAGTCTTTGCATCGCTTTAAAAAGCCGGGCTCTAAA
AACTTCCAGAATATCTTTCCACCATCAGCCACTCTGCATCTTTCCAACATTCCCCCTTCTGTTACAGTG
GATGATCTGAAGAACCTTTTATAGAAGCTGGATGTTCAAGTGAAGGCTTTTAAATCTTTTCAAGAAAT
CGCAAAATGGCGCTATTCAATTGGGATCTGTGGAAGAAGCAATTCAGGCCCTCATTGAGCTTCATAAC
CATGACCTTGAGAAAATCACACCTCAGAGTTTCCTTCTCAAAATCTACAATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001163788

Insert Size:	1575 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001163788.2
RefSeq Size:	7970 bp
RefSeq ORF:	1575 bp
Locus ID:	9991
UniProt ID:	O95758
Cytogenetics:	9q32
Protein Families:	Druggable Genome
MW:	56.8 kDa
Gene Summary:	<p>The protein encoded by this gene binds RNA and is a regulator of cell differentiation. The encoded protein preferentially binds to poly(G) and poly(U) sequences in vitro. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence and uses an alternate in-frame splice junction at the 3' end of an exon compared to variant 6. The resulting isoform (2) is shorter at the N-terminus and contains an alternate short internal segment compared to isoform 6. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>