

Product datasheet for **SC327076**

DISC1 (NM_001164549) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DISC1 (NM_001164549) Human Untagged Clone
Tag:	Tag Free
Symbol:	DISC1
Synonyms:	C1orf136; SCZD9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327076 representing NM_001164549.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCAGGCGGGGTCTCAGGGCGCCAGCCGCCGGCGGGCGCGGTGAGCCACCGCCAGGC
AGCCGGGATTGCTTACCACCTGCAGCGTGCTTTTCGGAGGGCGGCTGGCACGGAGGCCGGGCTACATG
AGAAGCTCGACAGGGCCCTGGGATCGGGTTCCTTTCCCCAGCAGTGGGCACACTGTTCCGGTTCACAGGA
GGGGTGTCTGGCAGGAGTCCCACCCTCGGAGTCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGC
CTTTGGTCCGGAGCCCTGTTTCCAAGAGTGCAGCAGCCCTACTGTGACCTCTGTGAGAGGAACCTCG
GCGCACTTTGGGATTCAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCCT
GGGAGTGTGGTGGCAGCAAGAGTTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGGGAG
GCAGCTGCAGCGATGGAGCAAGGCGTGTCCGGGCAGCAGGCTCTCTGCCATCAGCAGAGTTGAGTAGC
AACAGCTGCAGCCCTGGCTGTGGCCCTGAGGTCCCCCAACCCCTCTGGCTCTCACAGTGCCTTTACC
TCAAGCTTTAGCTTTATTCGGCTCTCGTTGGCTCTGCCGGGAACGTGGAGAAGCAGAAGGCTGCCCA
CCATCCAGAGAGGCTGAGTCCCATTGCCAGAGCCCCAGGAGATGGGAGCCAAAGCTGCCAGCTTGAC
GGGCTCACGAGGACCCGCGATGTCTCTCGGCCCTTCACTCTTTGGCTACACGGGTCTCTGCAGAC
TTGGCCAGGCCCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTCTTTACCAGACATGGACCCT
GGCTCCTCCAGTTCTCTGGATCCCTCACTGGCTGGCTGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT
GCCCACTTTGGGACACCCTGCTCAGGAAATGGGAGCCAGTGTGCGGGACTGCCTGTGAGAAACCGG
AGGCAGATGGAGGTAATATCCTTAAGATTAACACTTCAGAACTTCAGAAAGTGCAGTTGAGAATGAT
GATTATGATAAAGCTGAGACGTTACAACAAGATTAGAAGACCTGGAACAAGAGAAAAATCAGCCTGCAC
TTTCAACTTCTTCAAGGCAGCCAGCTCTTAGCAGTTCCTGGGTACCTGGCAGCACAAGTCCAGGCT
GCCTTGGCCCGTGGGCCACTCAGCAGGCCAGCGGAGATGACACCCACACCCCACTGAGAATGGAGCCG
AGGCTGTTGGAACCCACTGCTCAGGACAGCTTGCACGTGTCCATCACGAGACGAGACTGGCTTCTCAG
GAAAAGCAGCAGCTACAGAAAGAAATCGAAGCTCTCAAGCAAGGATGTTTGTGCTGGAAGCCAAAGAT
CAACAGCTGAGAAGGGAAATAGAGGAGCAAGAGCAGCAACTCCAGTGGCAGGGCTGCGACCTGACCCCA
CTGGTGGCCAGCTGTCCCTGGGTGAGTGCAGGAGGTCAGCAAGGCTTGCAGGACACCCTGGCCCTCA
GCCGTCAGATTCCTTCCATGCAGAGCCACCGAAACCATAAGGAGTACTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001164549
Insert Size:	1644 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164549.1</u>
RefSeq Size:	2279 bp
RefSeq ORF:	1644 bp
Locus ID:	27185
Cytogenetics:	1q42.2
MW:	58.9 kDa
Gene Summary:	<p>This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11)(q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (I) lacks multiple 3' exons but has an alternate 3' segment, as compared to variant L. The resulting isoform (I, also known as isoform 37) is much shorter and has a distinct C-terminus, as compared to isoform L.</p>